

QUICK START GUIDE

HI 902 Color

AUTOMATIC POTENTIOMETRIC TITRATOR

Revision 2.00



QUICK START GUIDE

Dear customer,

Congratulations on choosing a Hanna Instruments product.

This guide has been written for HI 902 titrators with color display, USB interface and software version 2.00 and higher.

Please read this Quick Start Guide carefully before using the instrument. This guide will provide you with the necessary information for the correct use of the instrument.

The purpose of this guide is to present a quick overview of setting up and using the instrument.

For detailed information illustrating the extensive capabilities of your titrator, please refer to the Instruction Manual.

© 2010 Hanna Instruments

All rights are reserved. Reproduction in whole or in part is prohibited without the written consent of the copyright owner, Hanna Instruments Inc., 584 Park East Drive, Woonsocket, Rhode Island 02895, USA.

Contents

INTRODUCTION4

SAFETY MEASURES.....5

TITRATOR CONNECTIONS6

USER INTERFACE7

HOW TO SELECT YOUR LANGUAGE8

HOW TO USE THE CONTEXTUAL HELP8

METHODS8

HOW TO CALIBRATE A pH ELECTRODE9

HOW TO PERFORM A TITRATION10

QUICK START GUIDE

INTRODUCTION

The HI 902 automatic titrator is designed to perform a wide variety of potentiometric titrations with high accuracy, flexibility and reproducibility, allowing the user to obtain both accurate results and high-speed analysis.

The titrators can perform fixed endpoint or equivalence point titrations by measuring the pH/mV and temperature of the sample.

Titration reports and methods can be transferred to a PC via a USB interface, saved to a USB storage device or printed directly from the titrator. An external monitor and keyboard can also be attached for added convenience.

How can I find certain information?

- The **Quick Start Guide** will help the user learn how to operate the titrator within a short period of time.
- The **Instruction Manual** provides a complete description of the operating principles (user interface, general options, methods, titration mode, pH, mV and ISE mode, maintenance, etc.).
- The **Titration Theory** outlines the basic concepts of titration.
- The contextual **Help** screens contain detailed explanations of every screen.

SAFETY MEASURES


The following safety measures must be followed:

1. Always ensure that the power supply cable is connected to a grounded main power plug.
2. Never connect or disconnect the pump assembly with the titrator turned on.
3. Verify that the burette and the attached tubing are assembled correctly.
4. Always check that the titrant bottle and the titration beaker are placed on a flat, stable surface.
5. Always wipe up spills and splashes immediately.
6. Avoid the following environmental working conditions:
 - Severe vibrations
 - Direct sunlight
 - Atmospheric relative humidity above 95% non-condensing
 - Environment temperatures below 10°C and above 40°C
 - Explosion hazards
7. Have the titrator serviced by qualified service personnel only.


QUICK START GUIDE

TITRATOR CONNECTIONS

Front View




Rear View



USER INTERFACE


Keypad

The titrator keypad has 29 keys grouped in five categories, as follows:



Display

The titrator has a 7.5" graphical backlit display. The main screen is shown below with short explanations.



The user interface contains several screens. In each screen, many information fields are present at the same time. The information is displayed in an easy-to-read manner.


Virtual option keys describe the function performed when the corresponding option key is pressed.

QUICK START GUIDE

HOW TO SELECT YOUR LANGUAGE

To change the language, press **General Options** from the main screen. Highlight the *Language* option and then press **Select**. Using the \triangle and ∇ keys, select the language from the options listed in the **Set Language** screen and press **Select**.

Restart the titrator in order to apply the new language setting.



HOW TO USE THE CONTEXTUAL HELP

Information about the titrator can be easily accessed by pressing **?**. The contextual help can be accessed at any time and it provides useful information about the current screen.

METHODS

The HI 902 titrator can store up to 100 methods.

Standard Methods

Each titrator is supplied with a package of standard methods. Standard method packs are developed at Hanna Instruments laboratories to meet analysis requirements of specific industries (e.g., water treatment, wine, dairy, etc.).

User Defined Methods

User defined methods allow the user to create and save their own methods. Each new method is based on an existing method which is altered to suit a specific application.

HOW TO CALIBRATE A pH ELECTRODE

To enter pH calibration mode, press then to enter pH mode. Once in pH mode, press to enter pH calibration mode.

PREPARATION

Pour small quantities pH 4.01, pH 7.01 and pH 10.01 buffer solutions into clean beakers. If possible, use plastic beakers to minimize any EMC interferences.

For accurate calibration and to minimize cross-contamination, use two beakers for each buffer solution: one for rinsing the electrode and one for calibration.

CALIBRATION PROCEDURE

Three buffer entry types are available: Automatic, Semi-automatic and Manual Selection.

The default option is Manual Selection.

- If the instrument has been previously calibrated and calibration was not cleared, the old calibration can be cleared by pressing .

Note: *It is very important to clear calibration history when a new electrode is used because most errors and warning messages that appear during calibration depend on calibration history.*

- Use the or to select pH 4.01 buffer solution.
- Use the second beaker of pH 4.01 buffer solution to rinse the pH electrode, temperature probe and propeller stirrer.
- Immerse the pH electrode, temperature probe and propeller stirrer in the pH 4.01 buffer solution. The pH electrode's bulb must be completely immersed in the buffer solution and the reference junction needs to be 5-6 mm below the surface. Add additional buffer if necessary.
- Press to turn on the propeller stirrer.
- Once the reading has stabilized, press to update the calibration.
- Repeat this procedure for pH 7.01 and 10.01 buffer solutions.
- Press to accept and exit pH calibration mode.

QUICK START GUIDE

HOW TO PERFORM A TITRATION

Required Solutions

- Titrant - 500 mL of 0.1 M (mol/L) Sodium Hydroxide (NaOH) in a titrant bottle.
- Sample - 0.1 mol/L Hydrochloric Acid (HCl).
- Distilled or deionized water.

Note: Analytical grade reagents and water should be used for accurate results.

Priming the Burette

- Insert the aspiration tube in the titrant bottle and the dispensing tube in a waste beaker.
- From the main screen press .
- Highlight the *Prime Burette* option and then press .
- Enter the number of burette rinses. At least 3 rinses are recommended.
- Press to start.
- The message "Executing..." will be displayed.

Note: Make sure you have continuous liquid flow inside the burette. For accurate results, the aspiration tube, the dispensing tube and the syringe must be free of air bubbles.

Method Selection

For this analysis we will use the **HI1009EN Neutralization w/ NaOH**.

To select this method:

- Press . Use the \triangle and ∇ keys to highlight **HI1009EN Neutralization w/ NaOH**.
- Press .

Setting Method Parameters

To display the method parameters, press . The **View/Modify Method** screen will be displayed.

Only certain parameters can be changed.

For this titration, the NaOH titrant concentration and the size of the HCl sample need to be entered.

To accomplish this:

- Highlight *Titrant Conc.* option, then press . The **Titrant Concentration** screen will be displayed.
- Enter the correct value, then press .
- Highlight *Analyte Size* option, then press .
- Enter the volume of the sample (ex: 5 mL), then press .
- Press , highlight *Save Method* option and then press .

Titrant Concentration

Enter the titrant concentration.

0.1000 M (mol/L)

ACCEPT	Escape	Delete Digit	
--------	--------	--------------	--

Sample Volume

Enter the initial sample volume in milliliters.

5 mL

This volume will be used when Fixed sample size is selected.

ACCEPT	Escape	Delete Digit	
--------	--------	--------------	--

Setup Titration Report

Users can select the information that is stored for each titration.

To setup the titration report, follow the procedure below:

- From the main screen, press . The **Data Parameters** screen will be displayed.
- Highlight *Setup Titration Report* and press .
- Mark the fields to be included in the titration report with the "*" symbol. Use the  and  keys to highlight a field and  /  to toggle the field.
- Press  to save the customized report.

QUICK START GUIDE

Preparing the Sample

- Add 50 to 65 mL of distilled / deionized water to the titration beaker.
- Use a pipette or burette to add 5.0 mL of the sample (0.1M Hydrochloric Acid (HCl)) into the same beaker.
- Slide the stirrer assembly up.
- Place the beaker under the stirrer assembly.
- Lower the stirrer assembly until it rests on the positioning collar.
- Adjust the height of the stirrer assembly so it is as close as possible to the bottom of the beaker.
- Adjust the level of the sample solution with distilled / deionized water so that the pH electrode bulb is completely immersed in the sample solution and the reference junction of the electrode is 5-6 mm below the surface.

Note: Make sure that the pH electrode, temperature probe and propeller do not touch each other or the beaker.

Performing a Titration

- From the main screen, press . You will be prompted to enter the analyte size. Enter 5 mL and press . The titrator will start the analysis.
- At the end of the titration, the message "Titration Completed" will appear on the display with the final concentration of the analyte in the sample and the equivalence endpoint volume.

Understanding the Displayed Information

During a titration the following screen is displayed:



11:37:43 AM May 16, 2010			
Analog 1			
Neutralization w/ NaOH			
0.40 meq/L			
Stirrer 1 1400/1400 RPM			
In Progress			
Pump 1 Running Burette: 25 mL			
ATC	Volume Delivered	pH	
20.8 °C	2.425 mL	5.775	
	View Curve	Suspend	Stop

Viewing Graph During Titration

After a few doses are dispensed View Curve will become active. Press View Curve to display the real-time titration graph.

The curves displayed are plots of the pH and the 1st derivative versus Titrant Volume (for details, see the Instruction Manual).



The two graphs are scaled to fit in the same screen window. Press Select Trace to change the y-axis scale to either the pH values or the 1st derivative values.



Titration Termination

The titration is normally terminated when the first equivalence endpoint is detected according to the selected algorithm. To ensure the correct detection and interpolation of the equivalence endpoint, the titrator will dispense a few additional doses after the endpoint was reached.

The titration result can be displayed either in the main screen or in the **Graph of Titration Data** screen:



When the titration has ended, the titrator will display the equivalence endpoint volume and the final concentration of the analyte together with the **Titration Completed** message.

To view the titration graph and/or results, press results.

When the titration ends, an "x" will mark the endpoint on the pH versus titrant volume curve in the **Graph of Titration Data** screen. The value of the endpoint volume is also displayed next to the endpoint.

QUICK START GUIDE

Results

The results obtained from a titration are stored in a report file that can be viewed, transferred to a USB Storage Device or PC, or printed.

Viewing the last titration data

- From the main screen, press . The **Data Parameters** screen will be displayed.
- From the **Data Parameters** screen highlight the *Review Last Titration Report* option and press . The **Review Result** screen will be displayed.
- Use the  and  keys to display information related to the last titration performed. See *titration report* on page 15.

Printing the titration report

Connect a DOS / Windows compatible printer directly to the DB 25-pin connector located on the back of the titrator.

Note: When connecting the printer, please turn off the titrator and the printer.

Printing out the report:


- From the **Review Report** screen, press .
- During the information transfer to the printer, the message "Printing" will be displayed on the screen.
- Press  to return to the **Data Parameters** screen.
- Press  again to return to the main screen.

Saving data to USB Storage Device

This feature allows saving the results of titrations or pH / mV / ISE logging sessions on a USB storage device.

- From the main screen, press , the **General Options** screen will be displayed.
- Highlight the *Save Files to USB Storage Device* option using the  and  keys.
- Insert the USB storage device into the USB socket.
- Press , the **List of Files on Titrator** screen will be displayed.
- Use the  or  keys to select the file type: "report files".
- Press  to transfer all available reports to the USB storage device, or highlight the name of the report file to be transferred and press





- Transferring a report file will automatically transfer the corresponding log file and titration graph (.BMP file, if applicable).
- Press Escape to return to the **General Options** screen.
- Press Escape again to return to the main screen.

Titration report

While scrolling with the Page Up and Page Down keys, the fields below can be seen on the titrator display or printed. The same information is available on the saved report file (Ti_00007.rpt in this example).

```
Titration Report

Method Name:                Neutralization w/NaOH
Time & Date:                10:43 Apr 16, 2010
Titration ID:              Ti_00007

Standardization Data
Buffer Potential Efficiency Temperature
Time and Date
4.006pH    169.9mV    100.7%    A  22.0°C
           10:20 Apr 16, 2010
7.020pH    -7.8mV     96.5%    A  22.0°C
           10:23 Apr 16, 2010
10.040pH  -178.6mV    96.5%    A  21.9°C
           10:25 Apr 16, 2010

GLP & Instrumentation Data
Sample Name:                Sample HCl-1
Company Name:              Hanna Instruments
Operator Name:
Electrode Name:           HI 1131 NO -2
Field 1:                  Any text
Field 2:                  Any text
Field 3:                  Any text
Titrator Software Version    v2.00
Base Board Software Version: v2.00
Pump 1 Software Version:    v1.4
Base Board Serial Number:   01040409
Analog Board Serial Number: 30040409
Pump 1 Serial Number:      70040207
Factory Calibration Date:   Jan 28, 2010

Method Parameters
Name:                      Neutralization w/NaOH
Method Revision:           1.0
Pump Configuration:
  Titrant Pump :           Pump 1
Analog Board:             Analog 1
```

QUICK START GUIDE

```

Stirrer Configuration:           Stirrer 1
Dosing Type:                     Dynamic
  min Vol:                       0.050 mL
  max Vol:                       0.500 mL
  delta E:                       20.000 mV
End Point Mode:                  pH 1EQ point, 1st Der
Recognition Options:
  Threshold:                      50 mV/mL
  Range:                          NO
  Filtered Derivatives:          NO
Pre-Titration Volume:           0.000 mL
Pre-Titration Stir Time:       15 Sec
Measurement Mode:               Signal Stability
  delta E:                       1.0 mV
  delta t:                       2 Sec
  t min wait:                    2 Sec
  t max wait:                    15 Sec
Electrode Type:                 pH
Calculations:                   Sample Calc. by Volume
Titrant Name:                   NaOH
Titrant Conc.:                  0.1000 M (mol/L)
Analyte Size:                   5.000 mL
Analyte Entry:                  Manual
Maximum Titrant Volume:        20.000 mL
Stirring Speed:                 1400 RPM
Potential Range:                -2000.0 to 2000.0 mV
Volume/Flow Rate:              25 mL / 50.0 mL/min
Signal Averaging:               1 Reading
Final Result Format:             X.XX
M (mol/L) -> M (mol/L)

```

```

V mol mol
-_*_*__

```

```

  L mol
  _____

```

```

mL L
-_*___

```

```

  1000mL

```

```

V = volume dispensed in liters.
0.100 mol/L -> titrant conc.
1.000 mol/mol -> (sample/titrant)
5.000 mL -> sample volume

```

Nr	Volume[ml]	mV	pH	Graphic	Temp[°C]	Time
0	0.000	235.2	2.857	0.0	A 19.1	00:00:00
1	0.050	234.6	2.866	-10.2	A 19.0	00:00:21
2	0.100	233.9	2.880	-15.8	A 19.1	00:00:27
3	0.200	232.2	2.908	-16.7	A 19.1	00:00:39
4	0.390	231.1	2.928	-6.0	A 19.1	00:00:45
5	0.590	228.6	2.970	-12.3	A 19.1	00:01:04
6	0.790	226.9	3.000	-8.7	A 19.1	00:01:20
7	0.990	225.5	3.024	-6.9	A 19.1	00:01:37
8	1.190	224.7	3.038	-4.0	A 19.1	00:01:43
9	1.390	223.9	3.051	-4.0	A 19.1	00:01:49

QUICK START GUIDE

10	1.590	223.0	3.066	-4.3	A	19.1	00:01:55
11	1.790	222.1	3.082	-4.6	A	19.1	00:02:01
12	1.990	221.2	3.098	-4.6	A	19.1	00:02:06
13	2.190	220.1	3.115	-5.1	A	19.1	00:02:11
14	2.390	219.0	3.134	-5.6	A	19.1	00:02:17
15	2.590	217.8	3.155	-6.0	A	19.1	00:02:23
16	2.790	216.5	3.177	-6.6	A	19.1	00:02:29
17	2.990	215.1	3.202	-7.3	A	19.1	00:02:34
18	3.190	213.4	3.231	-8.4	A	19.1	00:02:40
19	3.390	211.5	3.263	-9.3	A	19.1	00:02:46
20	3.590	209.2	3.302	-11.4	A	19.1	00:02:51
21	3.790	206.6	3.348	-13.4	A	19.1	00:02:57
22	3.990	203.2	3.406	-16.8	A	19.1	00:03:02
23	4.190	198.9	3.479	-21.4	A	19.1	00:03:08
24	4.390	193.1	3.578	-29.0	A	19.1	00:03:14
25	4.556	186.2	3.697	-41.7	A	19.1	00:03:20
26	4.670	179.6	3.810	-57.8	A	19.1	00:03:25
27	4.753	172.9	3.925	-81.2	A	19.1	00:03:31
28	4.812	166.4	4.036	-110.0	A	19.2	00:03:37
29	4.856	160.1	4.144	-143.5	A	19.2	00:03:43
30	4.889	153.7	4.253	-189.9	A	19.2	00:03:54
31	4.915	147.1	4.367	-259.9	A	19.2	00:04:00
32	4.934	141.0	4.471	-322.7	A	19.2	00:04:11
33	4.949	135.2	4.571	-388.0	A	19.2	00:04:17
34	4.964	127.5	4.702	-512.0	A	19.2	00:04:23
35	4.979	117.3	4.877	-680.0	A	19.2	00:04:29
36	4.994	104.2	5.102	-875.3	A	19.2	00:04:35
37	5.009	87.9	5.381	-1088.0	A	19.2	00:04:41
38	5.024	69.6	5.695	-1221.3	A	19.2	00:04:50
39	5.039	51.2	6.010	-1226.0	A	19.2	00:05:08
40	5.054	31.6	6.344	-1301.3	A	19.2	00:05:36
41	5.069	7.3	6.762	-1625.3	A	19.2	00:06:07
42	5.084	-37.9	7.557	-3010.0	A	19.2	00:06:38
43	5.099	-120.0	9.024	-5476.0	A	19.2	00:06:48
44	5.114	-144.7	9.464	-1642.7	A	19.2	00:06:54
45	5.129	-158.2	9.705	-900.7	A	19.2	00:07:01
46	5.144	-168.1	9.883	-664.0	A	19.2	00:07:08

Titration Results

Method Name: Neutralization w/NaOH
 Time & Date: 10:43 Apr 16, 2010
 Analyte Size: 5.000 mL
 End Point Volume: 5.090 mL
 pH Equivalence Point: 8.131
 Results: 0.10 meq/L
 Initial & Final pH: 2.857 to 9.884
 Titration Duration: 7:09 [mm:ss]
 Operator Name:

Analyst Signature: _____


QS 902C
10/10

INSTRUCTION MANUAL

HI 902 Color

AUTOMATIC POTENTIOMETRIC TITRATOR

Revision 2.00



Dear customer,

Congratulations on choosing a Hanna Instruments product.

Please read this instruction manual carefully before using the instrument. This manual will provide you with the necessary information for the correct use of the instrument.

This instruction manual has been written for the HI 902 titrators with color display, USB interface and software version 2.00 or later.

TABLE OF CONTENTS

- Chapter 1. INTRODUCTION**
 - Chapter 2. SETUP**
 - Chapter 3. USER INTERFACE**
 - Chapter 4. GENERAL OPTIONS**
 - Chapter 5. METHODS**
 - Chapter 6. TITRATION MODE**
 - Chapter 7. pH, mV and ISE MODE**
 - Chapter 8. AUXILIARY FUNCTIONS**
 - Chapter 9. MAINTENANCE, PERIPHERALS**
 - Appendix 1. TECHNICAL SPECIFICATIONS**
 - Appendix 2. BACK TITRATION**
 - Appendix 3. MULTIPLE EQUIVALENCE POINTS**
 - Appendix 4. LINKED METHOD**
 - Appendix 5. ACCESSORIES**
-

1 INTRODUCTION

HI 902 is an automatic potentiometric titrator with high accuracy, great flexibility and repeatability.

The titrator is designed to perform a variety of potentiometric titrations, allowing the user to obtain both good results and high-speed analysis.

The main attributes of this titrator are:

Flexibility	Support up to 100 titration methods (standard and user defined).
High accuracy	Precise dosing system. Precise mV and pH measurements (± 0.1 mV and ± 0.001 pH accuracy). Interpolated endpoint volume.
Repeatability	Powerful built-in algorithms for endpoint detection, including first derivative and second derivative detection algorithms, filtered derivatives, selectable range for equivalence point detection. Fixed mV or pH endpoint.
Quick results	Pre-defined titration methods. Pre-titration dosing feature. Dynamic / Linear dosing feature.
Complete report	The results are displayed directly on the LCD in the selected units. Titration graph can be displayed on the LCD and saved as a bitmap. Customizable titration reports can be printed, saved on a USB storage device or transferred to a PC via the USB interface.
Direct measurements	The titrator can also be used as a mV, pH, or ISE meter. Data logging is available for mV, pH, and ISE measurements.
GLP features	Up to 5 calibration points. Reminders for titrant age and standardization expiration. Fields for specific annotations.
Large graphical display	7.5" (320x240 pixels) graphical color display. Easy-to-view text and graphs. User-friendly interface.
Self diagnosis and integrated help	Integrated help screens are available. Self diagnosis features for peripheral devices including pump, valve, burette, stirrer. Error management with warning and error messages. Pre-defined titration methods for troubleshooting.

This manual provides information regarding installation and functionality of the titrator and refined operation suggestions.

Before using the titrator it is recommended you become familiar with its various features and functionality.

Chapter 2. Contents

2	SETUP	2 - 3
2.1	Unpacking	2 - 3
2.2	Safety Measures	2 - 4
2.3	Installation	2 - 5
2.3.1	Titrator Front View	2 - 5
2.3.2	Titrator Rear View	2 - 5
2.3.3	Titrator Left-side View	2 - 6
2.3.4	Titrator Assembly	2 - 6
2.3.4.1	Assembling Stirrer Stand and Support	2 - 6
2.3.4.2	Attaching Stirrer	2 - 7
2.3.4.3	Connecting the Pump	2 - 8
2.3.4.4	Attaching Burette Blank Support	2 - 8
2.3.4.5	Attaching the Burette	2 - 9
2.3.4.6	Electrical Connections	2 - 10

2 SETUP

2.1 Unpacking

The titrator and the accessories are shipped in a single box containing:

ITEM	QUANTITY
1 Titrator	1 pcs
2 Pump Assembly.....	1 pcs
3 Burette Assembly	1 pcs
• Burette (with 25-mL syringe)	
• Aspiration Tube with Fitting and Protection Tube	
• Dispensing Tube with Normal Dispensing Tip, Fitting, Protection Tube and Tube Guide	
• Tube Locks	
• Tool for Valve Fitting and Burette Cap Removal	
• Light Spectrum Protection Screen	
4 Stirrer Assembly.....	1 pcs
• Overhead Stirrer	
• Propeller (3 pcs)	
• Stirrer Stand	
• Stirrer support with positioning collar and positioning screw	
5 Burette Blank Support	1 pcs
6 Pump and Burette Locking Screws with Plastic Head	2 pcs
7 Temperature Sensor	1 pcs
8 Shorting Cap	1 pcs
9 Power Cable	1 pcs
10 USB Cable	1 pcs
11 Instruction Manual Binder	1 pcs
12 USB Memory Stick	1 pcs
13 HI 900 PC Application (Installation Kit on USB Stick)	1 pcs
14 Quality Certificate	1 pcs

See **Appendix 4** section **A 4.3 Titrator components** for pictures.

If any of the items are missing or damaged, please contact your sales representative.

Note: *Save all packing materials until you are sure that the instrument functions correctly. Any damaged or defective items must be returned in their original packing materials together with the supplied accessories.*

SETUP


2.2 Safety Measures

The following safety measures must be followed:


1. Always ensure that the power supply cable is connected to a grounded main power plug.
2. Never connect or disconnect the pump assembly with the titrator turned on.
3. Verify that the burette and the attached tubing are assembled correctly (see Section 9.1, *Burette Maintenance* for more details).
4. Always check that the titrant bottle and the titration beaker are on a flat surface.
5. Always wipe up spills and splashes immediately.
6. Avoid the following environmental working conditions:
 - Severe vibrations
 - Direct sunlight
 - Atmospheric relative humidity above 95% non-condensing
 - Environment temperatures below 10°C and above 40°C
 - Explosion hazards
7. Have the titrator serviced only by qualified service personnel.

2.3 Installation

2.3.1 Titrator Front View




2.3.2 Titrator Rear View



SETUP

2.3.3 Titrator Left-side View




2.3.4 Titrator Assembly

Note: Assembly operations must be completed before connecting the titrator to the power supply!

2.3.4.1 Assembling Stirrer Stand and Support

To assemble the stirrer stand and support:


- Remove the screws (B) from the titrator base (A).
- Attach the stirrer stand (C) to the titrator.
- Tighten the stirrer stand (C) using the previously removed screws (B).
- Screw the stirrer support (D) in the stirrer stand (C).




2.3.4.2 Attaching Stirrer

To attach the stirrer to the titrator, follow these steps:

- Attach the propeller (B) to the stirrer (A) by pressing it onto the stirrer shaft.
- Insert the electrode (C) and temperature sensor (D) into the dedicated holes on the stirrer. Push them in until they are in a stable position.



- Slide the stirrer on the stirrer support and set the height by tightening the screw located on the positioning collar (E).




SETUP

2.3.4.3 Connecting the Pump

To connect the pump, follow these steps:

- Retrieve the pump cable from inside the bay. The pump 1 connector is located in the left bay.
- Connect the cable (A) to the pump as shown below. The pump connector (B) is located in the lower part of the pump, near the motor.
- Lower the pump into the titrator, then slide it towards the front of the titrator chassis (C) until it is firmly latched.
- Secure the pump with the locking screw (D).

Figure 1



This procedure can be repeated to connect a second pump.

2.3.4.4 Attaching Burette Blank Support


To attach the burette blank support, follow these steps:

- Insert burette blank support in to the bay. Lower the burette blank support into the titrator, then slide it towards the front of the titrator chassis until it is firmly latched.
- Secure the burette blank support with the locking screw.

2.3.4.5 Attaching Burette

Make sure that the mark from the valve actuating cap and from the burette body are aligned as shown in Figure 2.

Figure 2



While ensuring the correct coupling between the syringe plunger (1) and the pump piston (2) (Figure 3), slide the burette into the support on the burette pump (Figure 4).




Figure 3





Figure 4

SETUP

2.3.4.6 Electrical Connections

- Connect the electrode to the BNC connector (K).
- Connect the temperature sensor to the RCA connector (C).
- Connect the stirrer to the MINI-DIN connector (D).
- Connect the power-supply cable to the power supply connector (B).



Nr	Function	Type of Connector
A	Power switch	
B	Power supply (115 or 230 VAC, 50-60 Hz) Check on titrator rear panel (M)	IEC Power Line Connector
C	Temperature sensor	RCA Socket
D	Stirrer	4-pin Mini DIN
E	RS232 interface (Balance Interface)	Standard DB 9 Pin Socket
F	Connector for expansion device	8-pin Mini DIN
G	External PC keyboard	6-pin Mini DIN (Standard PS2)
H	External display	Standard VGA Display 15-pin Socket
I	Printer	DB 25-pin Socket
J	USB interface	USB Standard B
K	Connection for pH, ORP and ISE half-cell or combination electrodes	BNC Socket
L	Reference electrode	Ø 4 mm Banana Socket


Chapter 3. Contents

3	USER INTERFACE	3-3
3.1	Start Up	3-3
3.2	Description	3-4
3.2.1	Keypad	3-4
3.2.1.1	Function Keys	3-4
3.2.1.2	Option Keys	3-4
3.2.1.3	Arrow Keys	3-5
3.2.1.4	Numeric Keys	3-5
3.2.1.5	Enter Key	3-5
3.2.2	Display	3-5
3.2.3	The Main Screen	3-6
3.3	Menu Navigation	3-7
3.3.1	Selecting an Option.....	3-7
3.3.2	Selecting a Menu Item	3-7
3.3.3	Entering Text	3-7
3.3.4	Saving Modifications	3-8

3 USER INTERFACE**3.1 Start Up**

Once the instrument is assembled and installed, follow the steps below to start the titrator:

- Connect the instrument to an outlet equipped with a ground wire. Make sure that the voltage of the main power matches that specified by the titrator.
- Turn on the titrator from the power switch located on the back of the instrument.
- Wait until the titrator finishes the initialization process.
- Press when prompted or wait a few seconds for titrator to start.



Note: All the performed initialization processes must be successfully completed. If one of them is terminated by a "Failed" message, restart the titrator from the power switch. If the problem persists, contact your dealer.


USER INTERFACE

3.2 Description

This chapter describes the basic principles of navigating through the user interface, selecting fields and entering values from the keypad.

3.2.1 Keypad

The titrator's keypad is grouped into five categories, as follows:



3.2.1.1 Function Keys

If one of these keys is pressed, the associated function is immediately performed. Some of the keys are active only in specific screens:

	Starts or stops a titration
	Turns the stirrer ON and OFF
	Reserved
	Access the results menu
	Displays contextual Help

3.2.1.2 Option Keys

These keys are assigned to the virtual keys on the display. Their functions are listed in the boxes above the buttons and vary depending on the displayed screen.

An underlined virtual key can also be activated by pressing .

3.2.1.3 Arrow Keys

These keys have the following functions:

- Move the on-screen cursor.
- Increase and decrease the stirrer speed and other settings.
- In the alphanumeric screen, to select a character.
- Navigate through menu options.

3.2.1.4 Numeric Keys

- Keys 0 to 9 Used for numeric entries.
- +/- Toggles between positive and negative values.
- . Decimal point.
- e Initiates entry of exponent for scientific notation.


3.2.1.5 Enter Key

Both , keys perform the same functions:

- Accept alphanumeric data entry.
- Executes the default (underlined) virtual option key.

3.2.2 Display

The titrator has a large color graphical display. The main screen is shown below with short explanations of the screen segments.



USER INTERFACE

The user interface contains several screens. For each titrator function, one or more screens are used.







3.2.3 The Main Screen

After start up and initialization, the first screen displayed is the main screen.

Main screen fields:

Method name:	Displays the name of the selected method.
Time and date:	Displays the current date and time.
Temperature reading:	Displays the measured temperature.
ATC:	Automatic temperature compensation
Manual:	Manual temperature compensation
Manual:	Temperature probe is not connected, manual temperature compensation
Stirrer information:	Actual / Set stirrer speed is displayed in RPM. When stirrer is off, the stirrer information is not displayed.
End point volume:	Displays the volume delivered to reach the titration end point. When no titration has been performed, the displayed volume is "0.000 mL".
Titration result:	Displays the titration result.
mV or pH reading:	Displays the current readings. The reading will be in mV or pH.
mV:	Indicates actual potential reading.
rel mV:	Indicates relative potential reading.
pH:	Indicates actual pH value.
Titration status:	Displays the status of the selected titration. No results is displayed when a titration has not been performed.
Reminders:	Indicates when a task needs to be performed and displays error or warning messages.
Pump 1 Selected:	Displays the active pump.
Analog 1 / Analog 2	When two analog boards are present, shows the active one.

3.3 Menu navigation



3.3.1 Selecting an Option

To select an option, simply press the option key below the virtual key. For example, to access the **Method Options** screen press the option key below it.

3.3.2 Selecting a Menu Item

To select an item from the menu screen use the arrow keys \triangle and ∇ to move the cursor.

When the menu is larger than the display, a scroll bar is active on the right side. The Page Up and Page Down keys can be used to scroll through the pages.

To activate the selected menu item, press enter or Select.

3.3.3 Entering Text

To enter text in an alphanumeric input box, first erase the previous text by using Delete Letter.

To enter a letter, highlight it using the arrow keys then press enter. Use the same procedure to enter the whole name.

For editing, use the Cursor Left and Cursor Right keys.


When editing is complete, press Accept.

The method name will be updated and displayed in the name field of the **View/Modify Method** screen.

When all the desired parameters have been set, press

Escape.

USER INTERFACE



3.3.4 Saving Modifications

The **Saving Method** screen allows the user to save the modifications. To exit from **Saving Method** screen without saving, press or highlight the *Exit Without Saving Method* option and then press . To save the modifications highlight the *Save Method* option and then press .


Note: To access the contextual help menu, press at any time. Help is related to the displayed screen. Press or press again to return to the previous screen.

Chapter 4. Contents

4	GENERAL OPTIONS	4 - 3
4.1	Temperature	4 - 3
4.1.1	Temperature Source	4 - 4
4.1.2	Manual Temperature Setting	4 - 4
4.1.3	Temperature Units	4 - 5
4.2	Date and Time Setting	4 - 5
4.3	Display Settings	4 - 6
4.4	Beeper	4 - 7
4.5	Stirrer	4 - 7
4.6	Language	4 - 8
4.7	Total Volume Alert	4 - 8
4.8	Titrant Age Reminder	4 - 9
4.9	Save Files to USB Storage Device	4 - 9
4.10	Restore Files from USB Storage Device	4-10
4.11	USB Link with PC	4-11
4.12	Setup Balance Interface	4-12
4.13	Printer Mode	4-13
4.14	Reset to Default Settings	4-13
4.15	Update Software	4-14


4 GENERAL OPTIONS

The **General Options** screen gives access to options that are not directly related to the titration process or pH / mV / ISE measurement. To access this screen, press  from the main screen. The available menus are described below:



4.1 Temperature

The **Temperature Menu** allows access to all of the settings related to temperature.



GENERAL OPTIONS

4.1.1 Temperature Source

Select the temperature source used for temperature compensation.

Temperature Source				
Select the temperature source.				
Automatic Temperature Manual Temperature				
Select	Escape			

When *Automatic Temperature Compensation* is selected "ATC" is displayed on the main screen and the temperature is read by the temperature probe.

When *Manual Temperature* is selected "Manual" is displayed on the main screen and a preset temperature value is used for temperature compensation.

Note: The selected temperature source will be indicated in the report files, A for Automatic and M for Manual.

4.1.2 Manual Temperature Setting

If the temperature probe is not connected, the user can manually set the temperature used by the titrator for compensation. This can be done when the *Manual Temperature* option is selected (see Section 4.1.1 *Temperature Source*).

Manual Temperature				
Enter the manual temperature to be used when the temperature probe is being overridden or no temperature probe.				
25.0 °C				
The temperature range is from -5.0 to 105.0°C.				
Accept	Escape	Delete Digit		

The temperature value can be set between -5 and 105 °C.

4.1.3 Temperature Units

The following temperature units can be selected.

Temperature Units				
Select the temperature units to be displayed.				
Celsius -5.0 to 105.0°C				
Fahrenheit 23.0 to 221.0°F				
Kelvin 268.2 to 378.2 K				
Select	Escape			

The temperature ranges are as displayed in the **Temperature Units** screen.

4.2 Date and Time Setting

This screen allows the user to set the date and time.

Date and Time Setting				
Enter the date.				
5	5	2010		
month	day	year		
Enter the time.				
13	5	33		
hour	minute	second		
Press Next to move to the next entry.				
ACCEPT	Escape	Delete Digit	Next	AM/PM

Use the  and  keys or the numeric keys to modify the date and time.


Press  to move the cursor to the next field.

Press  or  to change the time format.

GENERAL OPTIONS

4.3 Display Settings

This screen allows the user to customize the display settings.



Option Keys:



Increases the backlight saver time interval



Decreases the backlight saver time interval

The backlight intensity can be adjusted using \triangle and ∇ keys.

There are 8 levels of backlight intensity, ranging from 0 to 7.

A color palette is displayed in the center of the screen allowing an easy selection of the appropriate backlight intensity.

The backlight saver option protects the display during standby periods when no keys have been pressed for a set amount of time.

If the display backlight is off, any keystroke will activate the backlight without performing any action.

The range for the backlight saver timer is 1 to 60 minutes. To disable the backlight saver, increase the time to the maximum allowed. The "Off" indication will appear.

4.4 Beeper

This screen allows the user to be turn the Beeper On (Enable) or Off (Disable).

Beeper														
Select the option.														
<table border="1"><tr><td>Beeper Off</td><td colspan="4" style="background-color: black;"></td></tr><tr><td>Beeper On</td><td colspan="4"></td></tr></table>					Beeper Off					Beeper On				
Beeper Off														
Beeper On														
Select	Escape													

The beeper will sound after a titration is completed, when an invalid key is pressed or when a critical error occurs during titration.

4.5 Stirrer

This screen allows the stirrer to be enabled or disabled.

Stirrer														
Select the option.														
<table border="1"><tr><td>Disabled</td><td colspan="4" style="background-color: black;"></td></tr><tr><td>Enabled</td><td colspan="4"></td></tr></table>					Disabled					Enabled				
Disabled														
Enabled														
Select	Escape													

GENERAL OPTIONS

4.6 Language

Select an available language.

Set Language				
Select the language.				
English Portuguese Español				
Select	Escape			

4.7 Total Volume Alert

This screen allows a programmable reminder to appear when the titrant reservoir is below 100 mL. The titrant volume will decrease as the titrant is used.

Total Volume Alert				
Enter the amount of titrant available to the titration system from its reservoir. The mLs will decrease as the titrant is depleted.				
1000.0 mL				
A reminder will appear when less than 100 mLs of titrant volume is left.				
Accept	Escape	Delete Digit		Off

The "Low Titrant Volume" reminder message will appear when the available titrant volume is under 100 mL.

After the new titrant volume has been set on the titrator (in the **Total Volume Alert** screen), a warning message appears reminding the user to perform titrant re-standardization. The volume of titrant can be set from 0 to 10,000 mL.

4.8 Titrant Age Reminder

A programmable reminder will appear when it is time to verify the titrant concentration or to change the titrant.

<p>Titrant Age Reminder</p> <p>Enter the number of days to pass since the last Titr. Vol. updating or the last Start pressing, whereafter the reminder appears.</p> <p style="text-align: center;">30 days</p> <p>The range is from 0 to 31 days.</p>				
Start	Escape	Delete Digit		Off

The "Check Titrant Concentration" reminder will appear when the set number of days has passed since the total volume alert was set or since the timer was started. The reminder can be disabled by pressing . The range is from 0 to 31 days.


4.9 Save Files to USB Storage Device

This option allows the user to save files from the titrator to a USB storage device. On the titrator, the available file types are:

- | | |
|-----------------------|--|
| Standard Method Files | - HIXXXYY.MTD (e.g.: HI0001EN.MTD, HI1004EN.MTD) |
| User Method Files | - USERXXXX.MTD (e.g.: USER0001.MTD) |
| Report Files | - Ti_XXXXX.RPT, mV_XXXXX.RPT, pH_XXXXX.RPT, ISEXXXXX.RPT, mVrXXXXX.RPT (e.g.: Ti_00001.RPT, mV_00001.RPT, pH_00001.RPT, ISE00001.RPT, mVr00001.RPT) |

Use the ◀ and ▶ keys to select the file type. The number of files and each file name on the titrator will be displayed.

GENERAL OPTIONS



The option keys allow the following operations:

Delete File	Deletes the highlighted file
Delete All	Deletes all currently displayed files
Copy File	Copies the highlighted file from titrator to a USB storage device
Copy All	Copies all currently displayed files from titrator to a USB storage device
Escape	Returns to the General Options screen


The status of the transfer ("successful" / "unsuccessful") and the file name of the currently processed file are displayed during copying or deleting.

Note: The saved files will be stored on the USB key in the HI902 folder, as follows:

- Methods: **USB Drive:\HI902\Methods*.mtd**
- Reports: **USB Drive:\HI902\Reports*.rpt**

4.10 Restore Files from USB Storage Device

This screen allows the user to transfer files from the USB storage device to the titrator.



The file types that can be transferred are:

- | | |
|-----------------------|--|
| Standard Method Files | - HIXXXYY.MTD (e.g.: HI0001EN.MTD, HI1004EN.MTD) |
| User Method Files | - USERXXXX.MTD (e.g.: USER0001.MTD) |
| Report Files | - Ti_XXXXX.RPT, mV_XXXXX.RPT, pH_XXXXX.RPT, ISEXXXXX.RPT, mVrXXXXX.RPT (e.g.: Ti_00001.RPT, mV_00001.RPT, pH_00001.RPT, ISE00001.RPT, mVr00001.RPT) |

Use the ◀ and ▶ keys to select the file type.

The number of files and the name of each file found on the USB storage device is displayed on the screen.

The option keys allow the following operations:



- Deletes the highlighted file from the USB storage device.
- Deletes all currently displayed files from the USB storage device.
- Copies the highlighted file from diskette to the USB storage device.
- Copies all currently displayed files from diskette to the USB storage device.
- Returns to the **General Options** screen.

Note: In order to restore files from USB key, please ensure that the methods and / or reports you wish to transfer to the titrator are in the correct folder:

- Methods: **USB Drive: \HI902\Methods*.mtd**
- Reports: **USB Drive: \HI902\Reports*.rpt**

4.11 USB Link with PC

In order to use this feature, the USB cable needs to be connected from the titrator to the PC. Make sure that HI 900 PC application is running on the PC.


In the **USB Communication** screen:

“Active / Inactive”: shows the status of the USB link with the PC.

“Active” means that the titrator is using the USB communication with the PC and not with another device.

“Ready” shows that the titrator is able to communicate with the PC.

During transfer of any information between the PC and the titrator, “Transmit” and information about the percentage of current file already transferred are displayed.



GENERAL OPTIONS

4.12 Setup Balance Interface

This screen allows the users to connect an analytical balance for automatic acquisition of sample mass prior to titration or standardization.

```
Setup Balance Interface
Select the balance to be activated.
Lab Balance
Enable Balance  Escape  New Balance  Edit
```

The balance is connected to the titrator via RS 232 interface.

Press  to add a new balance to the list.

Press  to enable the balance interface feature.

Press  to disable the balance feature (automatic weight acquisition will be not available).

Press  to customize the serial communication parameters by accessing the **Balance Configuration** screen.

```
Balance Configuration
Select the option to be modified.
Balance Name Lab Balance
Baud Rate 9600
Data Bits 8 Bits
Parity No Parity
Stop Bit 1 bit
Edit Request Command B
Select  Escape  Test Balance
```

Be sure that the settings on the titrator *Balance Configuration* menu match the settings for your particular balance (baud rate, data bits, parity, stop bit number, request command syntax). It may be necessary to change settings on your balance. Users should consult their balance instruction manual.

Before leaving this screen be sure the connection with the balance is working properly by pressing the  key.

4.13 Printer Mode

This screen allows the users to select the printing mode: Ansi (default), Ascii and Text mode.

Printer Mode																			
Select the option.																			
<table border="1"><tr><td>Ansi</td><td></td><td></td><td></td><td></td></tr><tr><td>Ascii</td><td></td><td></td><td></td><td></td></tr><tr><td>Text</td><td></td><td></td><td></td><td></td></tr></table>					Ansi					Ascii					Text				
Ansi																			
Ascii																			
Text																			
Select	Escape																		

Ansi mode:

Use this mode when your printer is set as Ansi. In this case all the accented characters / symbols available in titrator will be printed on your printer.

Ascii mode:

Use this mode when your printer is set as Ascii. In this case only some of the accented characters / symbols available in titrator will be printed on your printer.

Text mode:

Use this mode when you don't need to print the accented characters.

4.14 Reset to Default Settings

This option restores the manufacturer settings.

Note: Please be careful !!! This will also delete all the user created methods and restore all manufacturer settings such as titrator configuration, standard method parameters, etc.

Confirmation of Reset				
Are you sure you want to reset the titrator to manufacturer settings?				
This will delete the standardization data, all the user methods and reports.				
Reset	Escape			

GENERAL OPTIONS

4.15 Update Software

This screen allows the user to update the titrator software from a USB storage device containing a software setup kit.

Update Software				
Current version:	HI902	v2.00		
New version:	HI902	v2.01		
Are you sure you want to update the current software with the new version?				
Accept	Escape	Refresh		

Chapter 5. Contents

5	METHODS	5 - 3
5.1	Selecting Methods	5 - 3
5.2	Standard Methods	5 - 4
5.2.1	Upgrading Standard Methods	5 - 4
5.2.2	Deleting Standard Methods	5 - 4
5.2.3	Restore the Standard Methods to the Manufacturer Settings	5 - 5
5.3	User Methods	5 - 5
5.3.1	Creating User Methods	5 - 5
5.3.2	Deleting User Methods	5 - 6
5.4	View / Modify Method	5 - 6
5.5	Method Options	5 - 7
5.5.1	Naming the User Method	5 - 7
5.5.2	Method Revision	5 - 8
5.5.3	Titration Type	5 - 8
5.5.4	Pump Configuration	5 - 9
5.5.5	Dosing Type	5 - 10
5.5.5.1	Linear Dosing	5 - 10
5.5.5.2	Dynamic Dosing	5 - 11
5.5.6	End Point Mode	5 - 13
5.5.6.1	Fixed End Point (pH or mV)	5 - 13
5.5.6.2	Equivalence End Point (pH or mV)	5 - 14
5.5.7	Recognition Options (Equivalence End Point (pH or mV) only)	5 - 17
5.5.7.1	Threshold	5 - 18
5.5.7.2	Range	5 - 19
5.5.7.3	Filtered Derivatives	5 - 20
5.5.8	Pre-Titration Volume	5 - 20
5.5.9	Pre-Titration Stir Time	5 - 21

METHODS


5.5.10 Measurement Mode	5 - 22
5.5.10.1 Signal Stability	5 - 22
5.5.10.2 Timed Increment	5 - 24
5.5.11 Electrode Type	5 - 25
5.5.12 Blank Option	5 - 25
5.5.13 Calculations	5 - 26
5.5.13.1 No Formula (mL only)	5 - 26
5.5.13.2 No Formula (L only)	5 - 26
5.5.13.3 Sample Calculation by Weight	5 - 27
5.5.13.4 Sample Calculation by Volume	5 - 28
5.5.13.5 Standardize Titrant by Weight	5 - 30
5.5.13.6 Standardize Titrant by Volume	5 - 30
5.5.13.7 Generic Formula	5 - 30
5.5.13.8 Back Titrations	5 - 33
5.5.13.8.1 Sample Calculation by Weight	5 - 33
5.5.13.8.2 Sample Calculation by Volume	5 - 35
5.5.13.8.3 Generic Formula	5 - 37
5.5.14 Dilution Option	5 - 37
5.5.15 Titrant Name	5 - 38
5.5.16 Titrant Concentration	5 - 38
5.5.17 Analyte Size	5 - 38
5.5.18 Analyte Entry	5 - 38
5.5.18.1 Fixed Weight or Volume	5 - 39
5.5.18.2 Manual Weight or Volume	5 - 39
5.5.18.3 Same as Previous (Linked Methods only)	5 - 39
5.5.19 Titrant 1 Entry (Back Titration only).....	5 - 39
5.5.19.1 Calculated by Formula	5 - 39
5.5.19.2 Fixed by User	5 - 39
5.5.20 Maximum Titrant Volume	5 - 40
5.5.21 Stirring Speed	5 - 40
5.5.22 Potential Range	5 - 41
5.5.23 Volume/Flow Rate	5 - 41
5.5.24 Signal Averaging	5 - 42
5.5.25 Final Result Format	5 - 42
5.5.26 Linked Method	5 - 43
6 Printing	5 - 43

5 METHODS

All of the parameters required to complete an analysis are grouped into a method. The titrator is supplied with a pack of standard methods. Standard and user methods can be upgraded, saved or deleted by connecting the titrator to a PC using the HI900 PC application or a USB storage device.


5.1 Selecting Methods

To select a method, press Select Method from the main screen. A list of available methods will be displayed.



In the **Titration Methods** screen, you can view the list of all available methods (standard and user methods).

To select a method, highlight the method and press Select. The name of the selected method will be displayed on the main screen.



METHODS

5.2 Standard Methods

The standard methods are developed for the most common types of analysis. Only specific method parameters can be modified by the user (see section 5.5 *Method Options*).

Also, standard methods can be used as models to create new user methods.

5.2.1 Upgrading Standard Methods

To upgrade the titrator with new standard methods, follow the steps below:

From USB Storage Device:

- Insert the USB storage device into the USB port, located on the left side of the titrator.
- Press  from the main screen.
- Using  and  keys, highlight the *Restore Files from USB Storage Device* option and choose .
- Using  and  keys, navigate through file types to find "standard method files". The list with available standard methods will be displayed.
- Press the  or  key to upgrade the titrator with the standard methods.
- Press  to return to **General Options** screen.

From PC:

You can upgrade the titrator with standard methods from a PC using the HI 900 PC application (see section 4.13 *USB Link with PC*).

5.2.2 Deleting Standard Methods

Unnecessary standard methods can be removed from the titrator by following the procedure below:

From General Options Screen:



- Using the  and  keys, highlight the *Save Files to USB Storage Device* option and press ;
- Using the  and  keys, navigate through the file types menu to find "standard method files". The available standard methods will be displayed.
- Press the  or  keys to remove unnecessary standard methods.
- Press  to return to the **General Options** screen.

From PC:

Unnecessary standard methods can be removed from the titrator using the HI 900 PC application (see section 4.11 *USB Link with PC*).

5.2.3 Restore the Standard Methods to the Manufacturer Settings

You can restore the standard methods to the manufacturer setting by highlighting a standard method and pressing .





5.3 User Methods

These methods are defined by the user (usually by modifying a standard method). The user methods can be developed in accordance with the requirements of the user. All method parameters can be modified by the user.

5.3.1 Creating User Methods

To create a new user method start from a standard or user method and follow these steps:

- Press  from the main screen.
- Using the  and  keys, highlight an existing method from the methods list.
- Press . A new user method will be generated.
- Press  to activate the new user method.




METHODS

Note: Only a limited number of user methods can be generated. The titrator can hold 100 methods (standard and user). When it is reached a warning message will be displayed.


5.3.2 Deleting User Methods

To remove a user method, press  from the main screen. Highlight the user method that you want to delete and press . A screen will appear in order to confirm the deletion. Press  again to confirm, or press  to cancel the operation.



5.4 View / Modify Method

To modify the method's parameters, press  from the main screen. A list of all the parameters for the selected method will be displayed. Using the  and  keys, highlight the option that you want to modify and choose .



If the *Back Titration* option is chosen, the following option must be set:

Break at Titrant Changing

Select the option.

NO
YES

"NO" - without break at titrant changing.
"YES" - with break at titrant changing.

Select	Escape			
--------	--------	--	--	--

Note: Selecting "YES" will stop titration temporarily between the first and the second phase of the back titration. A break in the titrant will allow you to perform a task related to the analysis (e.g.: boiling the sample to remove carbon dioxide).

5.5.4 Pump Configuration

This option allows you to choose the pump that will be used for the titration. For a Back Titration you can select the pump to be used for titrant 1 and titrant 2.

Titrant Pump Selection

Select the pump for titrant.


Pump 1
Pump 2

Select	Escape			
--------	--------	--	--	--

METHODS


5.5.5 Dosing Type

The Titrator has two dosing types: *Linear Dosing* and *Dynamic Dosing*.



5.5.5.1 Linear Dosing


Linear dosing dispenses a pre-defined volume of titrant with every addition (see graph).



The *Linear Dosing* option is recommended for titrations with a slower reaction rate, difficult nonaqueous titrations, and specific applications.

Note: For steep and normal titration curves, smaller volume increments are recommended, to obtain many measured points around the equivalence point.
For flat titration curves, larger volume increments are recommended for equivalence point detection.

To set the dosing volume, select *Linear Dosing* and enter the optimum dose.



Dosing volume ranges are:

5 mL burette	0.001	to	4.500 mL
10 mL burette	0.001	to	9.000 mL
25 mL burette	0.005	to	22.500 mL


5.5.5.2 Dynamic Dosing

The titrator determines the titrant dose by trying to maintain a certain potential change (ΔE) with each addition.

After a titrant dose, if the potential change is lower than the set ΔE , the next dose will be progressively increased until *max Vol* is attained. If the potential change is still lower than the set value, the titration will continue with the *max Vol* doses.

After a titrant dose, if the potential change is higher than the set ΔE , the next dose will be progressively decreased until *min Vol* is attained. If the potential change is still higher than the set value, the titration will continue with the *min Vol* doses.

The titrant is added in volumes that depend on the proximity of the end point as shown in the graph below.



METHODS

Dynamic dosing allows for larger doses far from the endpoint, reducing the total titration time. Closer to the end point, smaller doses are made, providing more data and improved accuracy.

Dynamic Dosing				
Enter min Vol, max Vol and delta E.				
0.015 mL - min Vol				
0.500 mL - max Vol				
4.500 mV - delta E				
Press Next to move to the next entry.				
Accept	Escape	Delete Digit	Next	

The following parameters must be set:

min Vol: The smallest dose to be dispensed during a titration.

The *min Vol* must be greater than or equal to:

0.001 mL for a 5 mL burette

0.001 mL for a 10 mL burette

0.005 mL for a 25 mL burette

max Vol: The largest dose to be dispensed during a titration.

The *max Vol* must be less than or equal to 4.000 mL.

delta E: Sets the fixed potential jump that has to be achieved after each titrant dose.

The allowed range is between 0.1 and 99.999 mV.

Recommendations for dosing parameters:

For steep and normal titration curves the recommended settings are:

delta E 3.5 to 9 mV

min Vol 0.010 to 0.025 mL (for a 25 mL burette)

max Vol 0.075 to 0.250 mL (for a 25 mL burette)

For flat titration curves the recommended settings are:


delta E 10 to 15 mV

min Vol 0.050 to 0.150 mL (for a 25 mL burette)

max Vol 0.400 to 0.600 mL (for a 25 mL burette)

To achieve the highest levels of accuracy and reproducibility, it is recommended that 20-80% of the nominal burette volume used for each titration is consumed. If lower volumes of titrant are required, a smaller burette can be used.


5.5.6 End Point Mode



5.5.6.1 Fixed End Point (pH or mv)

Fixed End Point (pH):

The titration is terminated when the preset pH value has been exceeded. The end point volume is a calculated value based on the dispensed volume when pH is under the preset value and the dispensed volume when pH exceeded the preset value.



The range is from - 2.000 to 20.000 pH.


METHODS

Fixed End Point (mV):

The end point detection algorithm is the same as for pH, but the threshold value is expressed in mV.

Preset mV End Point				
Enter the end point mV value.				
0.0 mV				
The range is from -2000.0 to 2000.0 mV.				
Accept	Escape	Delete Digit		

The range is from - 2000.0 to 2000.0 mV.



5.5.6.2 Equivalence End Point (pH or mV)

The titration is normally terminated when the equivalence point is detected (the point where the added quantity of titrant equals the quantity of analyte present in the sample).

Number of Equivalence Points

Up to 5 equivalence points can be detected.

Number of Equivalence Points				
Enter the number of equivalence points to be found.				
3 points				
The range is between 1 and 5 equivalence points.				
Accept	Escape	Delete Digit		

End Point Determination

The first and the second derivative of the S-shaped titration curve can be used to detect the equivalence point.

End Point Determination				
Select the end point determination.				
1st derivative				
2nd derivative				
Select	Escape			

The equivalence point detection algorithm requires three additional titrant doses to be dispensed after the equivalence point is reached.


The reported end point volume is a calculated value based on a number of points around the equivalence point.

The potentiometric S-shaped titration curve is the response in potential (mV) or pH between the indication of the electrode versus the volume of titrant added.

The inflection point of the S-shaped titration curve is assumed to be the equivalence point of the chemical reaction.


For non-symmetric titration curves, the theoretical error can be reduced by using the dynamic dosing.

METHODS



1st Derivative:


When first derivative is used to recognize the equivalence point, the titration curve inflection point (EQP) is the point where the first derivative reaches its maximum value.



The detection algorithm looks for the maximum value of the first derivative. The first derivative must be greater than the threshold value at the maximum point (see *Recognition Options* on page 5-17).

2nd Derivative:

When second derivative is used to recognize the equivalence point, the titration curve inflection point (EQP) is the point where the second derivative crosses zero.



The detection algorithm looks for the point where the second derivative changes sign. The checked point, or first derivative, must be greater than the threshold value (see *Recognition Options* on page 5-17).

5.5.7 Recognition Options (Equivalence End Point (pH or mV) only)

The **Recognition Options** screen is a set of parameters used to avoid false detection of the equivalence point due to the chemical system (titrant / sample species and concentrations) and / or electrode response.

The **Recognition Options** screen is available only when *Equivalence End Point (pH or mV)* option is selected.

Recognition Options

Select the options for equivalence point recognition.

Threshold	500 mV/mL
Range	NO
Filtered Derivatives	NO

Select
Escape

METHODS

5.5.7.1 Threshold

This parameter must be set by the user according to the analysis. The threshold represents the absolute value of the first derivative, expressed in mV/mL, below which the detection algorithm does not search for the equivalence point.

Threshold

Enter the threshold for equivalence point detection.


500 mV/mL

Recommended value is between:
 1 and 450 mV/mL for FLAT Curve,
 450 and 1800 mV/mL for NORMAL Curve,
 1800 and 9999 mV/mL for STEEP Curve.

Accept	Escape	Delete Digit	
--------	--------	--------------	--

Range is between 1 and 9999 mV/mL.

The recommended value is 40% of the absolute value of the first derivative.



Depending on the titration curve profile, the following guide can be used:

TITRATION CURVE PROFILE	THRESHOLD (mV/mL)
Flat	1 to 450
Normal	50 to 1800
Steep	1800 to 9999

5.5.7.2 Range

Range is an optional feature for equivalence point recognition. The titrator will only look for an equivalence point between the set values.


The *Range* option can be enabled by selecting *YES* in the **Range Options** screen.

<p style="text-align: center;">Range Options</p> <p>Select option for equivalence point range.</p> <div style="border: 1px solid black; padding: 2px;"><p>NO YES</p></div> <p>"NO" - without equivalence point range. "YES" - with equivalence point range.</p>					<p style="text-align: center;">Range Limits</p> <p>Enter Limit 1 and Limit 2 for range.</p> <p style="text-align: center;">-2.0 mV - Limit 1 20.0 mV - Limit 2</p>				
Select	Escape				Accept	Escape	Delete Digit	Next	

pH Range -2.000 to 20.000

mV Range -2000.0 to 2000.0

The Limit 2 value must not be equal to the Limit 1 value.



METHODS

5.5.7.3 Filtered Derivatives

This option adds a filtering procedure in the 1st and 2nd derivative computation algorithm that reduces the influence of pH or mV noise.

The *Filtered Derivatives* option can be enabled by selecting *YES* in the ***Filtered Derivatives Option*** screen.

Filtered Derivatives Option														
Select option for filtered derivatives.														
<table border="1"><tr><td>NO</td><td></td><td></td><td></td><td></td></tr><tr><td>YES</td><td></td><td></td><td></td><td></td></tr></table>					NO					YES				
NO														
YES														
"NO" - without filtered derivatives. "YES" - with filtered derivatives.														
Select	Escape													

Noise can be due to:

- Chemical system properties (sample, titrant, solvent), such as slow chemical reactions or unbuffered samples such as wastewater, tap water, wine
- Electrode response
- Incorrect method parameters settings such as *Signal Stability*, *Stirring Speed*, etc.
- Insufficient titrant additions

Note: A shift in the endpoint volume by 1 or 2 doses may be seen due to filtering.

5.5.8 Pre-Titration Volume

During a titration, the equivalence point is reached after many titrant doses. These doses take up extra time while having no relevance for equivalence point detection.

Pre-titration volume adds a large initial dose to jump directly to the proximity of the equivalence point.

This first dose occurs after the pre-titration stir time is completed.

The ranges for pre-titration volumes are shown below:

- 0.001 to 4.500 mL for a 5 mL burette
- 0.001 to 9.000 mL for a 10 mL burette
- 0.005 to 22.500 mL for a 25 mL burette

Pre-Titration Volume				
Enter the initial titrant volume to be dispensed.				
9.000 mL				
Press Help to view the valid ranges for the pre-titration volume.				
Accept	Escape	Delete Digit		

To disable a pre-titration volume, enter 0.000 mL.

Note: A pre-titration volume is highly recommended whenever possible. Fewer doses will considerably shorten the overall titration duration.

5.5.9 Pre-Titration Stir Time

When enabled, the sample is mixed for a set period of time before any titrant is added. This allows the sample to become homogeneous.

The range is from 0 to 180 seconds.


Pre-Titration Stir Time				
Enter the initial mixing time prior to the start of the titration.				
10 seconds				
The range is from 0 to 180 seconds.				
Accept	Escape	Delete Digit		

The *Pre-Titration Stir Time* option is disabled if 0 seconds is entered.

METHODS

5.5.10 Measurement Mode


During titration, the acquisition of the potential (mV) value of the solution can be done in two ways: by using either *Signal Stability* or *Timed Increment* option.



5.5.10.1 Signal Stability

When *Signal Stability* is selected, the titrator acquires the potential (mV) only when stable conditions are reached.

The principles of signal stability are plotted below:



The signal stability window (condition) represents the time interval (Δt) during which the potential measured in solution (mV) is confined inside the potential interval (ΔE). The new signal value is acquired if the stability condition is reached after the minimum (t_{\min}) wait time.

If the stability condition is not reached and the maximum (t_{max}) wait time has elapsed, the potential is acquired.

```
Signal Stability
Enter mV variation (delta E) in the time
interval (delta t) min and max wait time
period to the next sample measurement.

  0.3 mV      - delta E
  1.5 seconds - delta t
   5 seconds - t min wait
  30 seconds - t max wait

Accept  Escape  Delete  Next
       Digit
```

ΔE - maximum change in potential during Δt
The range is from 0.1 to 99.9 mV.

Δt - the time interval during which the potential is measured.
The range is from 0.5 to 10.0 seconds.

$t_{min\ wait}$ - the minimum elapsed time before a stability check. This is also the minimum elapsed time between two doses.
The range is from 2 seconds to $t_{max\ wait}$ time.


$t_{max\ wait}$ - the maximum elapsed time between two successive doses. If the $t_{max\ wait}$ has elapsed, a new dose is added even if the signal stability condition is not reached.
The range is from $t_{min\ wait}$ time to 180 seconds.

METHODS

5.5.10.2 Timed Increment

When *Timed Increment* is selected, the titrator acquires the potential (mV) at a fixed time interval (no signal stability check).

The time period between two acquisitions must be set according to the reaction and the response time of the electrode.



```
Timed Increment
Enter the period of time to wait until
the next dose.

      5 seconds

The range is from 2 to 180 seconds.

Accept  Escape  Delete
          Digit
```

The range is from 2 to 180 seconds.

5.5.11 Electrode Type

Enter the type of the electrode, up to 24 characters. The electrode type will appear in the titration report.

Electrode Type

Select the highlighted letter by using the arrow keys then press "Enter".
 Select the empty field for a space.
 Press Accept to save the electrode type.

B	C	D	E	F	G	H	I	J	K	L	M										
N	O	P	Q	R	S	T	U	V	W	X	Y	Z									
a	b	c	d	e	f	g	h	i	j	k	l	m									
n	o	p	q	r	s	t	u	v	w	x	y	z									
À	Á	Â	Ã	Ä	Å	Ç	È	É	Ê	Ë	Ï	Ñ	Ò	Ó							
Ô	Õ	Ö	Ù	Ú	Û	ä	å	à	á	â	ã	ç	è	é	ê	ë	ì	í	î	ï	\$
0	1	2	3	4	5	6	7	8	9	%	#	.									
,	?	!	()]	<	>	=	/	+	-									
*	\	_	`	^	'	:															

■ PH _____

Accept	Escape	Delete Letter	Cursor Left	Cursor Right
--------	--------	---------------	-------------	--------------

5.5.12 Blank Option

This feature allows the user to select the procedure for the blank calculations (where *V* is the volume of titrant dispensed during the titration and *Blank* is the volume of titrant consumed by the blank sample).

Blank Option

Select the option.

U - Blank

Blank - U

No Blank

Select	Escape		
--------	--------	--	--

If one of the options (*V-Blank* or *Blank-V*) is selected in the **View / Modify Method** screen, the *Blank Value* will be active on the View/Modify Method screen and the value of the blank can be set (in liters).

METHODS

Blank Value				
Enter the blank volume in liters.				
1.2530E-3 L				
ACCEPT	Escape	Delete Digit		

5.5.13 Calculations

The final result is computed using the end point volume (titrant volume at the equivalence point or at the fixed end point), and a formula selected by the user.

Calculations				
Select either the calculation to be performed or modify the variables.				
Edit Variable Values No Formula (mL only) No Formula (L only) Sample Calc. by Weight Sample Calc. by Volume Stdz. Titrant by Weight Stdz. Titrant by Volume Generic Formula				
Select	Escape			

Edit Variable Values

This option allows the user to edit the variables in a previously selected calculation. For each formula, selected variables can be changed.

5.5.13.1 No Formula (mL only)

If this option is selected, only the volume of titrant (mL) required to reach the end point is displayed.

5.5.13.2 No Formula (L only)

If this option is selected, only the volume of titrant (L) required to reach the end point is displayed.

5.5.13.3 Sample Calculations by Weight

This calculation is used when the concentration of an analyte is determined by the weight of the sample. The results are based on the initial sample weight (in grams). When you choose this formula, select the *Titrant Unit* first and then the *Final Result Unit*.

Titrant Units					Final Result Units				
Select the titrant unit.					Select the unit for your results.				
<div style="border: 1px solid black; padding: 2px;"> M (mol/L) N (eq/L) g/L mg/L </div>					<div style="border: 1px solid black; padding: 2px;"> ppt (g/kg) ppm (mg/kg) ppb (µg/kg) % = (g/100g) mg/g mg/kg mol/kg mmol/g eq/kg meq/kg </div>				
Select	Escape				Select	Escape			

The titrator will calculate the results based on the selected units.
 The titrator will provide the results based on the titrant and sample units selected.

Titrant Units:

- | | |
|-----------|-------------------|
| M (mol/L) | moles/liter |
| N (eq/L) | equivalents/liter |
| g/L | grams/liter |
| mg/L | milligrams/liter |

Final Result Units:

- | | |
|-------------|---|
| ppt (g/kg) | parts per thousand (grams/kilogram) |
| ppm (mg/kg) | parts per million (milligrams/kilogram) |
| ppb (µg/kg) | parts per billion (micrograms/kilogram) |
| % (g/100 g) | percentage in weight (grams/100 grams) |
| mg/g | milligrams/gram |
| mg/kg | milligrams/kilogram |
| mol/kg | moles/kilogram |
| mmol/g | millimoles/gram |
| eq/kg | equivalents/kilogram |
| meq/kg | milliequivalents/kilogram |

METHODS

A formula example is shown below using M (mol/L) as the titrant unit and g/kg as the final result unit:

Calculating Sample Concentration

M (mol/L) --> ppt (g/Kg)

The calculation is:

$$\frac{U \times \frac{\text{mol}}{\text{L}} \times \frac{\text{mol}}{\text{mol}} \times \frac{\text{g}}{\text{mol}}}{\frac{\text{g}}{\text{kg}} \times 1000\text{g}}$$

Select the variables to change value.
U = volume dispensed in liters.

1.000 mol/L -> titrant conc.

1.000 mol/mol -> (sample/titrant)

1.000 g/mol -> mw of sample

1.000 g -> sample weight

Select	Escape	Save / Exit	
--------	--------	----------------	--

Variables can be set according to the amount of sample and titrant used.

5.5.13.4 Sample Calculations by Volume

This calculation is used when the concentration of an analyte is determined in terms of the volume of sample. The results are based on the initial sample volume (in milliliters). When choosing the formula, select the *Titrant Unit* first and then the *Final Sample Unit*. The titrator will calculate the results based on the selected units.

Titrant Units

Select the titrant unit.

M (mol/L)

N (eq/L)

g/L

mg/L

Select	Escape		
--------	--------	--	--

Final Result Units

Select the unit for your results.

ppt (g/L)

ppm (mg/L)

ppb (µg/L)

M (mol/L)

N (eq/L)

g/L

mg/L

µg/L

mol/L

mmol/L

mg/mL

mg/100mL

g/100 mL

eq/L

Select	Escape	Page Up	Page Down
--------	--------	------------	--------------

Titrant Units:

M (mol/L)	moles/liter
N (eq/L)	equivalents/liter
g/L	grams/liter
mg/L	milligrams/liter

Final Result Units:

ppt (g/L)	parts per thousand (grams/liter)
ppm (mg/L)	parts per million (milligrams/liter)
ppb (µg/L)	parts per billion (micrograms/liter)
M (mol/L)	Molarity (moles/liter)
N (eq/L)	Normality (equivalents/liter)
mg/L	milligrams/liter
µg/L	micrograms/liter
mmol/L	millimoles/liter
mg/mL	milligrams/milliliter
mg/100 mL	milligrams/100 milliliters
g/100 mL	grams/100 milliliters
eq/L	equivalences/liter
meq/L	milliequivalences/liter

A formula example is shown below using N (eq/L) as the titrant units and ppt (g/L) as the final result units:

Calculating Sample Concentration

N (eq/L) --> g/L

The calculation is:

$$\frac{U \times \frac{\text{eq}}{\text{L}} \times \frac{\text{mol}}{\text{eq}} \times \frac{\text{g}}{\text{mol}}}{\text{mL} \times \frac{\text{L}}{1000\text{mL}}}$$

Select the variables to change value.
U = volume dispensed in liters.

1.000 eq/L -> titrant conc.

1.000 mol/eq -> (sample/titrant)

1.000 g/mol -> mw of sample

1.000 mL -> sample volume

Select	Escape	Save / Exit		
--------	--------	----------------	--	--

Variables can be set according to the amount of sample and titrant used.

METHODS

5.5.13.5 Standardize Titrant by Weight

This calculation is used when the concentration of the titrant is determined using a solid standard. Determination of the titrant concentration is based on the primary standard weight (in grams).

The calculation is based on the selected titrant unit. If the titrant unit is M (mol/L), the formula used to calculate the result is displayed below:

Titrant Units				
Select the titrant unit.				
<div style="border: 1px solid black; padding: 2px;">M (mol/L) N (eq/L) g/L mg/L</div>				
Select	Escape			

Calculating Titrant Concentration				
The titrant concentration unit is M (mol/L).				
The calculation is:				
$\frac{g \times \frac{mol}{g} \times \frac{mol}{mol}}{U}$				
Select the variables to change value. U = volume dispensed in liters.				
<div style="border: 1px solid black; padding: 2px;">0.200 g -> standard weight 204.23 g/mol -> mw of standard 1.000 mol/mol -> (titrant/standard)</div>				
Select	Escape	Save / Exit		

5.5.13.6 Standardize Titrant by Volume

This calculation is used when the concentration of the titrant is determined using a primary standard solution. Determination of the titrant concentration is based on the primary standard volume (in milliliters).

The titrator will perform the calculation based on the titrant unit selected.

The calculation is based on the selected titrant unit. If the titrant unit is N (eq/L) the formula used to calculate the result is displayed below:

Titrant Units				
Select the titrant unit.				
<div style="border: 1px solid black; padding: 2px;">M (mol/L) N (eq/L) g/L mg/L</div>				
Select	Escape			

Calculating Titrant Concentration				
The titrant concentration unit is N (eq/L).				
The calculation is:				
$\frac{mL \times \frac{L}{1000mL} \times \frac{eq}{L}}{U}$				
Select the variables to change value. U = volume dispensed in liters.				
<div style="border: 1px solid black; padding: 2px;">1.684 mL -> standard volume 2.375 eq/L -> standard conc.</div>				
Select	Escape	Save / Exit		



5.5.13.7 Generic Formula

The user can define their own calculation formula based on the final result units in a solid or liquid sample.

Final Result Units:

ppt (g/kg)	parts per thousand (grams/kilogram)
ppt (g/L)	parts per thousand (grams/liter)
ppm	parts per million (milligrams/kilogram)
ppm (mg/L)	parts per million (milligrams/liter)
ppb	parts per billion (micrograms/kilogram)
ppb (µg/L)	parts per billion (micrograms/liter)
% (g/100 g)	percentage in weight (grams/100 grams)
M (mol/L)	Molarity (moles/liter)
mg/g	milligrams/gram
N (eq/L)	Normality (equivalents/liter)
mg/kg	milligrams/kilogram
mg/L	milligrams/liter
mol/kg	moles/kilogram
µg/L	micrograms/liter
mmol/g	millimoles/gram
eq/kg	equivalents/kilogram
mmol/L	millimoles/liter
meq/kg	milliequivalents/kilogram
mg/mL	milligrams/milliliter
mg/100 mL	milligrams/100 milliliters
g/100 mL	grams/100 milliliters
eq/L	equivalents/liter
meq/L	milliequivalents/liter
No Unit	No result unit

The titrator will calculate the results based on the selected unit.
 The formula can be either for titrant standardization or sample analysis.



METHODS

Where:

C = the concentration of the titrant

F1 = general factor

F2 = general factor

F3 = general factor

S = sample size, in grams or milliliters

V = the volume delivered, in liters, to reach the preset or equivalence end point
(determined by the titrator)

General factors:

Weight Conversion:

One of the general factors should be a weight conversion factor.

Examples of concentration units:

mol/L	moles/Liter
eq/L	equivalents/Liter
g/L	grams/Liter
mg/L	milligram/Liter

Reaction Ratio:

The reaction ratio is the ratio between the analyte and titrant or standard and titrant.

Examples of ratios:

mol/mol	moles of sample/moles of titrant
mol/eq	moles of sample/equivalents of titrant
eq/mol	equivalents of sample/moles of titrant
mol/mol	moles of titrant/moles of standard
eq/mol	equivalents of titrant/moles of standard

Example: 2 moles of NaOH react with 1 mole of H₂SO₄

Unit Conversion factor:

Used to convert between various measurement units.

Examples: L/1000 → mL
g/1000 → mg

Weight Conversion factor:

Used to convert between weight measurement bases (kg, g, mg, μg or mole, mmole).

Example: g → mol

5.5.13.8 Back Titrations

Calculations							
Select either the calculation to be performed or modify the variables.							
<table border="1"> <tr> <td>Sample Calc. by Weight</td> </tr> <tr> <td>Sample Calc. by Volume</td> </tr> <tr> <td>Generic Formula</td> </tr> </table>					Sample Calc. by Weight	Sample Calc. by Volume	Generic Formula
Sample Calc. by Weight							
Sample Calc. by Volume							
Generic Formula							
Select	Escape						

5.5.13.8.1 Sample Calculations by Weight

Select the titrant 1 unit, the titrant 2 unit, and the final result unit.

Titrant 1 Units								
Select the titrant 1 unit.								
<table border="1"> <tr> <td>M (mol/L)</td> </tr> <tr> <td>N (eq/L)</td> </tr> <tr> <td>g/L</td> </tr> <tr> <td>mg/L</td> </tr> </table>					M (mol/L)	N (eq/L)	g/L	mg/L
M (mol/L)								
N (eq/L)								
g/L								
mg/L								
Select	Escape							

Titrant 2 Units								
Select the titrant 2 unit.								
<table border="1"> <tr> <td>M (mol/L)</td> </tr> <tr> <td>N (eq/L)</td> </tr> <tr> <td>g/L</td> </tr> <tr> <td>mg/L</td> </tr> </table>					M (mol/L)	N (eq/L)	g/L	mg/L
M (mol/L)								
N (eq/L)								
g/L								
mg/L								
Select	Escape							

Final Result Units														
Select the unit for your results.														
<table border="1"> <tr> <td>ppt (g/Kg)</td> </tr> <tr> <td>ppm (mg/Kg)</td> </tr> <tr> <td>ppb (µg/Kg)</td> </tr> <tr> <td>% = (g/100g)</td> </tr> <tr> <td>mg/g</td> </tr> <tr> <td>mg/Kg</td> </tr> <tr> <td>mol/Kg</td> </tr> <tr> <td>mmol/g</td> </tr> <tr> <td>eq/Kg</td> </tr> <tr> <td>meq/Kg</td> </tr> </table>					ppt (g/Kg)	ppm (mg/Kg)	ppb (µg/Kg)	% = (g/100g)	mg/g	mg/Kg	mol/Kg	mmol/g	eq/Kg	meq/Kg
ppt (g/Kg)														
ppm (mg/Kg)														
ppb (µg/Kg)														
% = (g/100g)														
mg/g														
mg/Kg														
mol/Kg														
mmol/g														
eq/Kg														
meq/Kg														
Select	Escape													

METHODS

A formula example is shown below using M (mol/L) as the titrant 1 units, M (mol/L) as the titrant 2 units, mg/g and the final result units. This formula is used to calculate the amount of titrant 1 to dispense:

```

Calc. Direct Titr. Volume

Titr1 Unit: M (mol/L)-->Result Unit: L

The calculation is:

      g x mol x f
      g
-----
      mol x mol
      L x mol

Select the variables to change value.

1.000 g -> sample weight
1.000 g/mol -> mw of sample
1.000 f ->(excess factor)
1.000 mol/L -> titrant 1 conc.
1.000 mol/mol -> (sample/titrant 1)
  
```

Select	Escape		Next	
--------	--------	--	------	--

The formula is based on the assumption that the sample concentration is 100% w/w. The titrator will calculate the volume of titrant 1 needed to consume the sample and multiply it with the excess factor in order to raise or lower the amount of titrant 1 dispensed. Variables can be set according to the amount of sample and titrant used.

Press to proceed to the next formula.

If you do not want the titrator to calculate the volume of titrant 1 to add, see Section 5.5.18 *Titrant 1 Entry*.

The remaining volume of titrant 1 needs to be calculated.

The following formula is used to calculate the remaining volume of titrant 1 after the reaction with the sample:

```

Calc. Excess Volume Of Titr1.

M (mol/L) + M (mol/L) --> L

The calculation is:

      U2 x mol x mol
      L x mol
-----
      mol
      L

U1 = U1tot - U1 excess
Select the variables to change value.
U2 = backtitr. vol. dispensed in liters.

1.000 mol/L -> titrant2 conc.
1.000 mol/mol -> (titrant1/titrant2)
1.000 mol/L -> titrant1 conc.
  
```

Select	Escape		Next	
--------	--------	--	------	--

When all of the variables are set, press Next to proceed with the "Calculating Sample Concentration" formula:

Calculating Sample Concentration

Final unit is mg/g:
The calculation is:

$$\frac{U1 \times \frac{\text{mol}}{\text{L}} \times \frac{\text{mol}}{\text{mol}}}{\frac{\text{g} \times \frac{\text{mol}}{\text{g}} \times \frac{\text{g}}{1000\text{mg}}}$$

Select the variables to change value.
U1 = volume dispensed in liters

1.000 mol/L -> titrant1 conc.
1.000 mol/mol -> (sample/titrant1)
1.000 g -> sample weight
1.000 g/mol -> mw of sample

Select	Escape	Save / Exit	
--------	--------	----------------	--

5.5.13.8.2 Sample Calculations by Volume

Select the titrant 1 unit, titrant 2 unit, and the final result unit.

Titrant 1 Units

Select the titrant 1 unit.

M (mol/L)
N (eq/L)
g/L
mg/L

Select	Escape		
--------	--------	--	--

Titrant 2 Units

Select the titrant 2 unit.

M (mol/L)
N (eq/L)
g/L
mg/L

Select	Escape		
--------	--------	--	--

Final Result Units

Select the unit for your results.

ppt (g/L)
ppm (mg/L)
ppb (ug/L)
M (mol/L)
N (eq/L)
g/L
mg/L
ug/L
mol/L
mmol/L
mg/mL
mg/100mL
g/100 mL
eq/L

Select	Escape	Page Up	Page Down
--------	--------	------------	--------------

METHODS

After you have selected the titrant 1, titrant 2, and the final result units, the titrator will display a screen with a formula used to calculate the amount of titrant 1 (used in the first stage of back titration) to be dispensed.

Calc. Direct Titr. Volume

Titr1 Unit: M (mol/L)-->Result Unit: L

The calculation is:

$$\frac{\text{mL} \times \frac{\text{L}}{1000\text{mL}} \times \frac{\text{g}}{\text{L}} \times \frac{\text{mol}}{\text{g}} \times f}{\frac{\text{mol}}{\text{L}} \times \frac{\text{mol}}{\text{mol}}}$$

Select the variables to change value.

1.000 mL -> sample volume

1.000 g/L -> sample max conc.

1.000 g/mol -> mw of sample

1.000 f ->(excess factor)

1.000 mol/L -> titrant1 conc.

Select	Escape		Next
--------	--------	--	------

The formula is based on the assumption that the sample concentration is 100% v/v. The titrator will calculate the volume of titrant 1 needed to consume the sample and multiply it with the excess factor in order to raise or lower the amount of titrant 1 dispensed.

Variables can be set according to the amount of sample and titrant used.

Press Next to proceed to the next formula.

If you do not want the titrator to calculate the volume of titrant 1 to add, see Section 5.5.18 *Titrant Entry*.

The following formula is used to calculate the remaining volume of titrant 1 after the reaction with the sample.

Calculating Sample Concentration

Final unit is mg/g:

The calculation is:

$$\frac{U1 \times \frac{\text{mol}}{\text{L}} \times \frac{\text{mol}}{\text{mol}}}{\frac{\text{g}}{\text{g}} \times \frac{\text{mol}}{\text{g}} \times 1000\text{mg}}$$

Select the variables to change value.

U1 = volume dispensed in liters

1.000 mol/L -> titrant1 conc.

1.000 mol/mol -> (sample/titrant1)

1.000 g -> sample weight

1.000 g/mol -> mw of sample

Select	Escape	Save / Exit	
--------	--------	----------------	--

When all the variables are set, press Next to proceed with the "Calculating Sample Concentration" formula:

Calculating Sample Concentration

Final unit is g/L:
The calculation is:

$$\frac{U1 \times \frac{\text{mol}}{\text{L}} \times \frac{\text{mol}}{\text{mol}}}{\text{mL} \times \frac{\text{L}}{1000\text{mL}} \times \frac{\text{mol}}{\text{g}}}$$

Select the variables to change value.
U1 = volume dispensed in liters

1.000 mol/L -> titrant1 conc.

1.000 mol/mol -> (sample/titrant1)

1.000 mL -> sample volume

1.000 g/mol -> mw of sample

Select	Escape	Save / Exit	
--------	--------	----------------	--

5.5.13.8.3 Generic Formula

This option allows the user to define their calculation formula for the "Direct Titration Volume", "Calculating Excess Volume of Titrant 1" and "Final Sample Concentration" in a solid or liquid sample.

5.5.14 Dilution Option

When the initial sample is diluted, a titration is made with an aliquot of the diluted sample, dilution calculations can be used.

The calculations are based on the original sample weight (volume) in order to express the results for the initial sample.

Dilution Parameters

Select the option.

Final Dilution Volume: 100.000 mL

Aliquot Volume: 10.000 mL

Analyte size to be diluted: 1.000 g

Select	Escape		
--------	--------	--	--

Final Dilution Volume: The volume of the sample after dilution

Aliquot Volume: Sample volume used for the titration

Analyte size to be diluted: The initial sample weight (volume)

METHODS

The sample size used in the calculations:

$$\frac{\text{Analyte size to be diluted} * \text{Aliquot Volume}}{\text{Final Dilution Volume}}$$

5.5.15 Titrant Name

Enter the name of the titrant (up to 24 characters). This name will appear in the titration report.

Titrant Name

Select the highlighted letter by using the arrow keys then press "Enter".
Select the empty field for a space.
Press Accept to save the entered text.

A	B	C	D	E	F	G	H	I	J	K	L	M											
N	O	P	Q	R	S	T	U	V	W	X	Y	Z											
a	b	c	d	e	f	g	h	i	j	k	l	m											
n	o	p	q	r	s	t	u	v	w	x	y	z											
À	Á	Â	Ã	Ä	Å	Ç	È	É	Ê	Ë	Ì	Í	Î	Ï	Ñ	Ò	Ó	Ô	Õ	Ö	×	÷	\$
0	1	2	3	4	5	6	7	8	9	%	#	.											
,	?	!	()	[]	<	>	=	/	+	-											
*	\	_	`	^	'	:																	

0.1N NaOH

Accept	Escape	Delete Letter	Cursor Left	Cursor Right
--------	--------	---------------	-------------	--------------

5.5.16 Titrant Concentration

Enter the concentration of the titrant to be used. When determining the titrant concentration only the concentration unit is displayed. The titrant concentration can not be set.

5.5.17 Analyte Size

Enter the size of the sample (for sample concentration determinations) or standard (for titrant concentration determination).

5.5.18 Analyte Entry

Select the analyte entry mode.

Analyte Entry

Select the entry mode of analyte.

Fixed Weight or Volume

Manual Weight or Volume

Verify the correct formula is being used,
I.E. weight or volume analyte type.

Select	Escape			
--------	--------	--	--	--

5.5.18.1 Fixed Weight or Volume

Each titration will use a set weight or volume in the calculations.

5.5.18.2 Manual Weight or Volume

Each titration the exact weight or volume can be entered. The titrator will prompt for the analyte weight volume at the beginning of each titration.

5.5.18.3 Same as previous (Linked Method Only)

The same weight or volume is used for both methods.

5.5.19 Titrant 1 Entry (Back Titration Only)

Select the mode for evaluating the necessary quantity of titrant 1 used in the back titration process (phase 1).

Titrant 1 Entry						
Select the entry mode of titrant 1.						
<table border="1"> <tr> <td>Calculated By Formula</td> </tr> <tr> <td>Fixed By User</td> </tr> </table>					Calculated By Formula	Fixed By User
Calculated By Formula						
Fixed By User						
Select	Escape					

5.5.19.1 Calculated by Formula

The volume of titrant 1 to be dispensed in the phase 1 of back titration will be calculated by the titrator (see section 5.5.12.2 *Back Titrations*).

5.5.19.2 Fixed by User

A fixed volume of titrant 1 will be used during the first phase of back titration process (direct titration).

Direct Titration Volume						
Enter the volume of titrant which will be dispensed during direct titration.						
<table border="1"> <tr> <td>10.000</td> <td>mL</td> </tr> </table>					10.000	mL
10.000	mL					
This volume will be dispensed when Fixed By User option is selected.						
ACCEPT	Escape	Delete Digit				

METHODS

5.5.20 Maximum Titrant Volume

The maximum titrant volume used in the titration must be set according to the analysis. If the titration end point (fixed or equivalence End Point) is not reached, the titration will be terminated after the maximum titrant volume has been dispensed. The error message ("Limits Exceeded") will appear on the display.

Maximum Titrant 2 Volume				
Enter the maximum titrant volume to be dispensed.				
25.000 mL				
Recommend the total volume of the burette.				
Accept	Escape	Delete Digit		

Range is from 0.100 to 100.000 mL.

5.5.21 Stirring Speed

The stirring speed can be set between 100 and 2500 RPM with 100 RPM resolution.

Stirring Speed				
Enter the speed of the stirrer during the titration.				
1000 RPM				
The range is from 0 to 2500 RPM.				
Accept	Escape	Delete Digit		

The stir speed from the current method is used during the entire titration.

The speed can be adjusted at any time using the \triangle and ∇ keys when the stirrer is on.

5.5.22 Potential Range

The input potential range can be set by the user. The titration will be terminated and an error message will appear if the potential is outside these limits.

These limits provide protection against a titration that does not generate an end point due to potential over-range.

Potential Range				
Enter the upper and lower potential.				
<div style="background-color: black; color: white; display: inline-block; padding: 2px;">2000.0</div> mV - Upper Limit -2000.0 mV - Lower Limit				
Press Next to move to the next entry.				
ACCEPT	Escape	Delete Digit	Next	

The ranges can be set between -2000.0 to 2000.0 mV.

5.5.23 Volume/Flow Rate

The flow rate for the dosing system can be set by the user in an interval of 0.1 to two times the burette volume:

- 0.1 to 10 mL/min for a 5 mL burette
- 0.1 to 20 mL/min for a 10 mL burette
- 0.1 to 50 mL/min for a 25 mL burette

Flow Rate				
Enter the titrant flow rate.				
<div style="background-color: black; color: white; display: inline-block; padding: 2px;">50.0</div> mL/min				
The range is from 0.1 to twice the total volume of the burette.				
ACCEPT	Escape	Delete Digit		

Note: The titrator will automatically detect the burette size, and display the correct high limit volume.

The flow rate is set for all burette operations.

METHODS

5.5.24 Signal Averaging

This option enables filtering on the mV/pH reading.

If *1 Reading* is selected, the filtering is disabled. The titrator will take the last reading and place it into a "moving window" along with the last 2, 3 or 4 readings (depending on the selected option). The average of those readings is displayed and used for calculations.

Signal Averaging																								
Select the number of readings to be averaged.																								
<table border="1"><tr><td>1 Reading</td><td colspan="4" style="background-color: black;"></td></tr><tr><td>2 Readings</td><td colspan="4"></td></tr><tr><td>3 Readings</td><td colspan="4"></td></tr><tr><td>4 Readings</td><td colspan="4"></td></tr></table>					1 Reading					2 Readings					3 Readings					4 Readings				
1 Reading																								
2 Readings																								
3 Readings																								
4 Readings																								
Select	Escape																							

Averaging more readings is helpful when a noisy signal is received from the electrode.


5.5.25 Final Result Format

This option allows you to set the format for displaying the final titration result.

Significant Figures																								
Select the desired format for displaying the final titration result.																								
<table border="1"><tr><td>XX</td><td colspan="4"></td></tr><tr><td>XXX</td><td colspan="4"></td></tr><tr><td>XXXX</td><td colspan="4"></td></tr><tr><td>XXXXX</td><td colspan="4" style="background-color: black;"></td></tr></table>					XX					XXX					XXXX					XXXXX				
XX																								
XXX																								
XXXX																								
XXXXX																								
Select	Escape																							

5.5.26 Linked Method

This option allows the user to link two titration methods. If *No Link* is selected only the current method will run. If a method is selected it will run after the current method. See Appendix 4 for additional information.



5.6 Printing

To print method parameters, press Method Options from the main screen.

Press Print Method and wait a few seconds until the printer completes the job.

If no printer is connected to the dedicated socket, or if the printer is offline, an error message will appear on the display (see 9.3.3 *Connecting a Printer* section, for information about connecting a printer to the titrator).

Chapter 6. Contents

6	TITRATION MODE	6 - 3
6.1	Running a Titration	6 - 3
6.1.1	Starting a Titration	6 - 3
6.1.2	Suspending a Titration	6 - 3
6.1.3	Viewing the Titration Curve	6 - 3
6.2	Stopping a Titration	6 - 5

6 TITRATION MODE

6.1 Running a Titration

Before beginning a titration, make sure that the following conditions are met:

- At least one pump is properly installed.
- A burette is inserted in the pump and filled with titrant.
- The aspiration tube is inserted in the titrant bottle and primed. The dispensing tube is over the titration beaker.
- The standard or sample has been carefully weighed / measured into the titration beaker.
- The electrode(s) and the temperature probe are inserted in the titration beaker.
- The desired method is selected and the parameters are set to the optimal values.

6.1.1 Starting a Titration

To start a new titration, press  from the main screen.

When a titration begins:

- The stirrer will turn on (if detected and enabled).
- If the pre-stirring time option is enabled, the sample will be stirred until the set time elapses (see section 5.5.9 *Pre-Titration Stir Time*).
- If the pre-titration volume option is enabled, the set volume will be dispensed (see section 5.5.8 *Pre-Titration Volume*).
- The titrator will start the analysis and continue to deliver titrant until the end point is detected or the titration is terminated.

6.1.2 Suspending a Titration

While a titration is in progress, you can temporarily stop it by pressing . The burette will stop dispensing titrant.

To continue the titration press .

6.1.3 Viewing the Titration Curve

During a titration, the potentiometric curve and the derivative curve (equivalence point only) can be displayed on the **Titration Graph** screen, by pressing . The titration report ID is also displayed inside the graph window.

TITRATION MODE

The potentiometric curve and the derivative curve are scaled to fit simultaneously inside the display.

When a titration end point is successfully detected, the volume is displayed on the graph and marked with an "x":

The contents of the graph as related to an end point type are as follows:

Equivalence End Point (pH) - the pH readings and the selected derivative vs. volume of titrant are displayed (see Figure 1).

Equivalence End Point (mV) - the mV readings and the selected derivative vs. volume of titrant are displayed (see Figure 2).

Fixed End Point (pH) - the pH readings vs. volume of titrant are displayed (see Figure 3).

Fixed End Point (mV) - only the mV vs volume curve is displayed (see Figure 4).

Figure 1




Figure 2




Figure 3





Figure 4



Select Trace

- changes the y-axis scale to either the mV (or pH) readings or the selected derivative values (of mV or pH). Available only for titrations with equivalence endpoints.

Save as Bitmap

- allows you to save the graph as a bitmap file. Available only when the titration is finished.

6.2 Stopping a Titration

The titration is terminated when one of the following conditions is met:

- **Titration Completed.** This is the only mode with valid final result values. The endpoint was successfully detected. The final results will be displayed.
- **Manually Terminated.** The current titration terminated by the user before the endpoint was detected.
- **Limits Exceeded.** The maximum titrant volume was delivered without reaching the endpoint. An error message is displayed on the screen.
- **Critical Error.** A critical error occurred and the titration was stopped. These errors are typically related to the dosing system. An error message is displayed on the screen.
- **Potential Out of Range.** The measured values from the electrode are outside the potential range. An error message is displayed on the screen.

Chapter 7. Contents


7	pH, mV & ISE MODE	7 - 3
7.1	pH Mode	7 - 4
7.1.1	Display	7 - 4
7.1.2	pH Setup	7 - 5
7.1.2.1	Buffer Entry Type	7 - 5
7.1.2.2	First Calibration Point	7 - 6
7.1.2.3	Edit Custom Buffers	7 - 6
7.1.2.4	Edit Buffer Group	7 - 7
7.1.2.5	Calibration Reminder	7 - 8
7.1.2.6	Set Reminder Period	7 - 8
7.1.2.7	Clear Calibration	7 - 9
7.1.2.8	pH GLP Data	7 - 9
7.1.2.9	Logging Interval	7 - 10
7.1.2.10	Stability Criteria	7 - 10
7.1.2.11	pH Resolution	7 - 11
7.1.2.12	Stirrer Configuration	7 - 11
7.1.2.13	Stirring Speed	7 - 12
7.1.3	pH Calibration	7 - 12
7.2	mV Mode	7-15
7.2.1	Display	7 - 15
7.2.2	mV Setup	7 - 15
7.2.2.1	Clear Relative mV Offset	7 - 16
7.2.2.2	Logging Interval	7 - 16
7.2.2.3	Stability Criteria	7 - 17
7.2.2.4	Stirrer Configuration	7 - 17
7.2.2.5	Stirring Speed	7 - 18
7.2.3	mV Rel Calibration	7 - 18
7.3	ISE Mode	7-19
7.3.1	Display	7 - 19
7.3.2	ISE Setup	7 - 20
7.3.2.1	Calibration Group	7 - 20
7.3.2.2	Temperature Compensation	7 - 21
7.3.2.3	Isopotential Point	7 - 21

7.3.2.4	Edit Custom Standards	7 - 22
7.3.2.5	Edit Standard Group	7 - 22
7.3.2.6	Calibration Reminder	7 - 23
7.3.2.7	Set Reminder Period	7 - 23
7.3.2.8	Clear Calibration	7 - 24
7.3.2.9	ISE GLP Data	7 - 24
7.3.2.10	Electrode Type	7 - 25
7.3.2.11	Concentration Unit	7 - 26
7.3.2.12	Logging Interval	7 - 26
7.3.2.13	Stability Criteria	7 - 27
7.3.2.14	ISE Significant Digits	7 - 27
7.3.2.15	Stirrer Configuration	7 - 28
7.3.2.16	Stirring Speed	7 - 28
7.3.3	ISE Calibration	7 - 28
7.4	Logging	7-30
7.4.1	Automatic Logging	7 - 30
7.4.2	Manual Logging	7 - 30

7 pH, mV & ISE MODE

By pressing from the main screen, the titrator can be switched to **Titrator**, **pH**, **mV** or **ISE** modes.

One Analog Board




Switches to **Titrator** mode.

Switches to **pH** mode.

Switches to **mV** mode.

Switches to **ISE** mode.

Two Analog Boards



Switches to **Titrator** mode.

or

Switches the Analog Input for pH, mV and ISE mode.

Switches to **pH** mode.

Switches to **mV** mode.


Switches to **ISE** mode.

pH MODE ISE MODE

7.1 pH Mode

7.1.1 Display

The **pH** screen is shown below with short explanations of the screen segments.




pH Mode Option keys:

- General Options** The General Options screen gives access to options that are not directly related to the measurement process (see chapter 4 *General Options* for more information).
- Save Reading** Stores the current pH reading (see section 7.4.2 *Manual Logging*).
- OR**
- Start Log** Starts the pH automatic log (see section 7.4.1 *Automatic Logging*).
- pH Calibr.** Enter the pH calibration screen (see section 7.1.3 *pH Calibration*).
- pH Setup** Enter the pH setup screen, parameters are associated with pH measurements and calibration (see section 7.1.2 *pH Setup*).
- Mode** Allow the user to switch between the available measurement modes: **Titration**, **pH**, **mV** or **ISE** mode.

7.1.2 pH Setup

To access the pH Setup, press  option key while in pH mode.



Use  and  keys to highlight the desired option.

Press  or  to access the selected option.


7.1.2.1 Buffer Entry Type

Select the pH buffer entry mode used for calibration:

Automatic - the instrument automatically selects the pH calibration point as the closest buffer from the predefined Buffer Group (see section 7.1.2.4 Edit Buffer Group).

Semiautomatic - the instrument automatically selects the closest buffer from Available Buffers (standard and custom buffers).

Manual - the calibration buffer must be manually selected by the user during calibration from the available buffer list (standard and custom buffers).



pH, mV & ISE MODE

7.1.2.2 First Calibration Point

Two options are available for the First Calibration Point: *Point* and *Offset*.

If *Point* option is selected, the slope values adjacent to the calibration points will be reevaluated (normal calibration).

If at least a two-point calibration has been performed and an offset correction is needed, perform a one-point calibration using the *Offset* option. The existing slope values will not be changed.

First Cal Point						
Select a menu option.						
<table border="1"><tr><td>Point</td></tr><tr><td>Offset</td></tr></table>					Point	Offset
Point						
Offset						
Select	Escape					

7.1.2.3 Edit Custom Buffers

If you wish to use buffers other than the standard ones, the Edit Custom Buffers option is available, allowing you to set the desired pH buffers. Up to five pH custom buffers can be set.

Note: Custom buffers are not temperature compensated. The value of the buffer at the calibration temperature should be entered. The standard buffers are automatically temperature compensated.

Edit Custom Buffers														
Press <Edit> to edit selected buffer.														
Press <Invalid. Buffer> to invalidate the custom buffer.														
<table border="1"><tr><td>Cust</td><td>Cust</td><td>Cust</td><td>Cust</td><td>Cust</td></tr><tr><td>6.870</td><td>9.230</td><td>12.750</td><td>----</td><td>----</td></tr></table>					Cust	Cust	Cust	Cust	Cust	6.870	9.230	12.750	----	----
Cust	Cust	Cust	Cust	Cust										
6.870	9.230	12.750	----	----										
Use arrows keys to select the buffer.														
Invalid. Buffer	ESCAPE	Edit	◀	▶										

- Use ◀ and ▶ keys to select the desired buffer.
- Press to invalidate the selected buffer.
- Press to edit the selected buffer; use the numeric keys to edit the buffer values.
- Press to save the value.
- Press to return to pH Setup menu.

7.1.2.4 Edit Buffer Group

Select up to five buffers from the available buffers (Hanna or Custom) to be used for automatic buffer recognition (Automatic Buffer Entry Type).

Within the Buffer Group, pH values must be at least 1.5 pH far apart.

If the Buffer Group already contains five pH buffers, at least one pH buffer has to be removed in order to add another buffer.

The screenshot shows the 'Edit Buffer Group' menu. It is divided into two main sections: 'Available Buffers' and 'Buffer Group'. The 'Available Buffers' section contains five buffer slots: three 'Hanna' buffers with values 9.177, 10.010, and 12.450; one 'Cust' buffer with value 6.870; and one 'Cust' buffer with value 9.230. Below these is another 'Cust' buffer with value 12.750. The 'Buffer Group' section contains five slots: three 'Hanna' buffers with values 1.679, 4.010, and 7.010; one 'Cust' buffer with value 9.230; and one 'Hanna' buffer with value 12.450. At the bottom of the screen are five buttons: 'Remove', 'Escape', a right-pointing arrow, an up-pointing triangle, and a down-pointing triangle.


- Use the arrow keys to select the pH buffer to be included/removed in/from the buffer group.
- Press or to add/remove the selected pH Buffer to/from buffer group.
- Press to return to pH Setup menu.

pH, mV & ISE MODE

7.1.2.5 Calibration reminder

In order to have accurate readings, the electrode must be calibrated frequently. Three options are available for calibration reminder:

- Daily* - the calibration reminder will appear daily at specified time.
- Periodic* - the calibration reminder will appear after the set time has elapsed since the last calibration.
- Disable* - the calibration reminder will not appear.




7.1.2.6 Set Reminder Period

If Daily or Periodic option was selected for the Calibration Reminder, the reminder period must also be set.

For a daily reminder period, the time of day can be set.


For a periodic reminder period, the number of days, hours and minutes can be set.



- Press to move the cursor to the next field.
- Press to save the changes or to return to the previous screen.
- Press to disable the calibration reminder and return to pH setup.

7.1.2.7 Clear Calibration


This option clears the existing pH calibration for the selected channel. If the calibration is cleared, another calibration has to be performed.



- Press to clear the previous calibration or to return to the previous screen without clearing the calibration.

7.1.2.8 pH GLP Data

Displays the pH calibration data.



pH, mV & ISE MODE

7.1.2.9 Logging Interval

Set the logging interval to be used for automatic logging.

Logging Interval				
Enter the data logging interval.				
0	0	2		
hours	minutes	seconds		
Press Next to move to the next entry.				
Accept	Escape	Delete Digit	Next	Off

7.1.2.10 Stability Criteria

Select the signal stability criteria:

- Fast quicker results with less accuracy
- Medium medium speed results with medium accuracy
- Accurate slower results with high accuracy

Stability Criteria				
Select a menu option.				
Fast				
Medium				
Accurate				
Select	Escape			

7.1.2.11 pH Resolution

Set the desired pH resolution: one(x.x), two (x.xx) or three (x.xxx) decimals.

pH Resolution							
Select the pH resolution for direct reading.							
<table border="1"><tr><td>X.X</td></tr><tr><td>X.XX</td></tr><tr><td>X.XXX</td></tr></table>					X.X	X.XX	X.XXX
X.X							
X.XX							
X.XXX							
Select	Escape						

7.1.2.12 Stirrer Configuration

Set the stirrer configuration: Stirrer 1, Stirrer 2, or Disabled.

Stirrer Configuration							
Select the option.							
<table border="1"><tr><td>Disabled</td></tr><tr><td>Stirrer 1</td></tr><tr><td>Stirrer 2</td></tr></table>					Disabled	Stirrer 1	Stirrer 2
Disabled							
Stirrer 1							
Stirrer 2							
Select	Escape						

pH, mV & ISE MODE

7.1.2.13 Stirring Speed

The stirring speed for the selected stirrer can be set.

Stirring Speed				
Enter the speed of the stirrer within below range.				
██████████ 1100 RPM				
The range is from 0 to 2500 RPM.				
Accept	Escape	Delete Digit		

7.1.3 pH Calibration

Calibrate the instrument often, especially if high accuracy is required.
The instrument should be recalibrated:

- Whenever the pH electrode is replaced.
- At least once a week.
- After testing aggressive chemicals.
- When "No pH Calibration" or "pH Calibration Expired" message appears on the LCD, in the Reminder messages area.

pH Calibration				
1.674 _{pH}				
Manual	25.0 °C		Unstable mV	318.0
Calibrated Buffers				
Hanna 1.679	Hanna 4.010	Hanna 7.010	Cust 9.230	Hanna 12.450
Last Calibration: 12:16:09 Aug 05, 2009				
Press <Clear Cal> to clear old calibr. Press <Edit> to edit manual temperature.				
Clear Cal	Escape	Edit		

PREPARATION

Pour small quantities of the buffer solutions into clean beakers. If possible, use plastic beakers to minimize any EMC interferences.

For accurate calibration and to minimize cross-contamination, use two beakers for each

buffer solution: one for rinsing the electrode and one for calibration. If you are measuring in the acidic range, use pH 7.01 or 6.86 as the first buffer and pH 4.01/3.00 or 1.68 as the second buffer. If you are measuring in the alkaline range, use pH 7.01 or 6.86 as the first buffer and pH 10.01/9.18 or 12.45 as the second buffer.

For extended range measurements (acidic and alkaline), perform a five-point calibration by selecting five buffers across the entire pH range.

CALIBRATION PROCEDURE

During calibration, the user has a choice of 8 standard buffers: pH 1.68, 3.00, 4.01, 6.86, 7.01, 9.18, 10.01, 12.45 and up to 5 custom buffers.

For accurate measurements it is recommended to perform a five-point calibration. However, at least a two-point calibration is suggested. For pH titrations the selected buffers should bracket your endpoint (e.g.: if your endpoint value is at 8.5, use 7.01 and 9.18 for calibration).

Three buffer entry types are available: Automatic, Semiautomatic and Manual Selection (see section 7.1.2.1 *Buffer Entry Type*).

To begin calibration:

- Press . If the instrument was calibrated before and the calibration was not cleared, the old calibration can be cleared by pressing .

Note: *It is very important to clear calibration history when a new electrode is used.*

- Immerse the pH electrode and the temperature probe approximately 4 cm (1.5") into a buffer solution and stir gently. The temperature probe should be close to the pH electrode.
- Select the pH calibration buffer value with or .
- Press to update the calibration. Once the reading has stabilized, the calibration buffer will be added to the Calibrated Buffers section.
- Rinse the pH electrode and the temperature probe, then immerse them into the next buffer solution and follow the above procedure or press to exit the procedure.

Notes: *The new calibration points will replace old ones if the difference between them is +/- 0.2 pH. Buffers used in older calibrations will not have a solid background.*

If calibrating with a standard buffer in MTC mode, the pH value and temperature can be modified by pressing . The values can be adjusted using the numeric keys.

Press to save the new values.

pH, mV & ISE MODE

Manual Edit				
Edit pH buffer and manual temperature.				
Buffer: ████████ 7.005 pH				
Temperature: 25.0 °C				
Low limit: 6.990 pH				
High limit: 7.030 pH				
Press Next to move to the next entry.				
Accept	Escape	Delete Digit	Next	

In ATC mode, the pH value can be modified by pressing Edit.

If the Automatic buffer entry type was selected for the calibration procedure, the instrument will automatically select the closest buffer to the measured pH value from the edit buffer group (see section 7.1.2.4 Buffer Group Edit).

If the Semiautomatic buffer entry type was selected for the calibration procedure, the instrument will automatically select the closest buffers to the measured pH value from all the available buffers and the buffer value can be selected with Previous Buffer or Next Buffer.


CALIBRATION MESSAGES:

- **Wrong Buffer. Please check the buffer:** This message appears when the difference between the pH reading and the value of the selected calibration buffer is significant. If the message is displayed, check if you have selected the appropriate calibration buffer.
- **Wrong buffer temperature:** This message appears if the buffer temperature is out of the defined temperature range.
- **Clean the electrode or check the buffer. Press Accept to update calibration:** This message alerts the user that some dirt or deposits could be on the electrode.
- **Slope too low. Please check the buffer:** This message appears if the current slope is under 80% or over 110% of default slope. Recalibrate the instrument using fresh buffers.
- **Slope too high. Press Clear Cal to clear the old calibration:** This message appears as a result of an erroneous slope condition. Follow displayed instructions.

7.2 mV Mode

7.2.1 Display

The **mV** screen is shown below.



mV Mode Option Keys:



The General Options screen gives you access to options that are not directly related to the measurement process (See chapter 4 *General Options* for more information).



Stores the current mV reading (see section 7.4.2 *Manual Logging*).

OR



Starts the mV automatic log (see section 7.4.1 *Automatic Logging*).



Enter the relative mV calibration screen (see section 7.2.3 *Relative mV Calibration*).




Enter the mV setup screen. Parameters are associated with mV measurement and calibration.



Allows the user to switch between the available measurement modes: Titrator, pH, mV or ISE mode.

7.2.2 mV Setup



pH, mV & ISE MODE

7.2.2.1 Clear Relative mV Offset

Clear the relative mV offset and return to absolute mV measurement.

Clear Relative mV Offset				
Press Clear to clear the relative mV offset.				
Press Escape to return without clearing the relative mV offset.				
Clear	Escape			

- Press to clear the relative mV offset or to return to the previous screen.

7.2.2.2 Logging Interval

Set the logging interval.

Logging Interval				
Enter the data logging interval.				
0 0 2				
hours minutes seconds				
Press Next to move to the next entry.				
Accept	Escape	Delete Digit	Next	Off

7.2.2.3 Stability Criteria

Select the signal stability criteria:

- Fast quicker results with less accuracy
- Medium medium speed results with medium accuracy
- Accurate slower results with high accuracy

Stability Criteria				
Select a menu option.				
<div style="border: 1px solid black; padding: 2px;">Fast Medium Accurate</div>				
Select	Escape			

7.2.2.4 Stirrer Configuration

Set the stirrer configuration: Stirrer 1, Stirrer 2, or Disabled.

Stirrer Configuration				
Select the option.				
<div style="border: 1px solid black; padding: 2px;">Disabled Stirrer 1 Stirrer 2</div>				
Select	Escape			

pH, mV & ISE MODE

7.2.2.5 Stirring Speed

The stirring speed for the selected stirrer can be set.

Stirring Speed				
Enter the speed of the stirrer within below range.				
████████ 1100 RPM				
The range is from 0 to 2500 RPM.				
Accept	Escape	Delete Digit		

7.2.3 mV Rel Calibration


Relative mU				
Set the value for the relative mU offset.				
Absolute mU: -212.6 mU				
Relative mU: ████████ 12 mU				
Low limit: -2212.6 mU				
High limit: 1787.4 mU				
Accept	Escape	Delete Digit		

- Press to accept the value.
- Press to delete the last digit.
- Press to cancel this operation and return to the previous screen.

7.3 ISE Mode

7.3.1 Display

The **ISE** screen is shown below.



ISE Mode option keys:

- General Options

The General Options screen gives access to options that are not directly related to the measurement process (see chapter 4 *General Options* for more information).
- Save Reading

Stores the current concentration reading (see section 7.4.2 *Manual Logging*).
- or
- Start Log

Starts the ISE automatic log (see section 7.4.1 *Automatic Logging*).
- ISE Calibration

Enter the ISE calibration screen (see section 7.3.3 *ISE Calibration*).
- ISE Setup


Enter the ISE setup screen. Parameters are associated with ISE measurements and calibration.
- Mode

Allows the user to switch between the available measurement modes: Titrator, pH, mV and ISE mode.

pH, mV & ISE MODE

7.3.2 ISE Setup

To access the ISE Setup, press  option key in ISE mode.




7.3.2.1 Calibration Group

Selecting the set of available standards to be used in calibration:

All standards: the set of available standards includes the Standard Solutions and Custom Solutions.

Standards Group: the set of available standards includes the standards selected by the user.



7.3.2.2 Temperature Compensation

Enable or disable temperature compensation for ISE measurements.

Temperature Compensation						
Select the option.						
<table border="1"><tr><td>Disabled</td></tr><tr><td>Enabled</td></tr></table>					Disabled	Enabled
Disabled						
Enabled						
Select	Escape					

Note: If you enabled Temperature Compensation, then the isopotential point must be set.

7.3.2.3 Isopotential Point

This option is available only if temperature compensation is enabled.


This option allows the user to set an isopotential point for the selected electrode. Ion selective electrodes have different isopotential points. The isopotential point is edited in ppm units only. The isopotential point should be entered if it is known and if measurements are going to be made at several temperatures.

Isopotential Point					
Enter the value for isopotential point.					
<table border="1"><tr><td>20</td></tr></table> ppm					20
20					
Low limit: 1E-2 ppm					
High limit: 1E+5 ppm					
Accept	Escape	Delete Digit			

pH, mV & ISE MODE

7.3.2.4 Edit Custom Standards


Edit the custom standard list. Up to five can be used in calibration.



- Use the ◀ and ▶ keys to select the standard.
- Press to invalidate the standard.
- Press to edit the selected custom standard; use the numeric keys to edit the standard.

7.3.2.5. Edit Standard Group

Select up to 5 standards from the available standards (Predefined and Custom) to be used during calibration.




- Use the arrow keys to select the standard to be included/removed in/from the standard group.

- Press or to add/remove the selected standard to/from standard Group.
- Press to return to ISE Setup menu.

7.3.2.6. Calibration Reminder

In order to have accurate readings, the electrode must be calibrated frequently. Three options are available for the calibration reminder:

- Daily* - the calibration reminder will appear daily at specified time.
- Periodic* - the calibration reminder will appear after the set time has elapsed since the last calibration.
- Disable* - the calibration reminder will not appear.




7.3.2.7 Set Reminder Period

If Daily or Periodic option was selected for the Calibration Reminder, the reminder period must also be set.

For a daily reminder period the time of day can be set.

For a periodic reminder period the number of days, hours and minutes can be set.



pH, mV & ISE MODE

- Press to move the cursor to the next field.
- Press to save the changes or to return to the previous screen.
- Press to disable the calibration reminder and return to pH setup.

7.3.2.8 Clear Calibration

This option clears the existing ISE calibration. If the calibration is cleared, a new calibration must be done in order to take measurements.

Clear Calibration				
Press <Clear> to clear all calibration points.				
Press <Escape> to return without clearing the calibration points.				
Clear	Escape			

- Press to clear the previous calibration or to return to the previous screen.


7.3.2.9 ISE GLP Data

Displays the ISE calibration data.

ISE GLP Data				
Analog 2				
Last Calibration: 02:36 PM Feb 24, 2010				
Slope: 101.0% ISE: Calcium				
10.0 ppm, 1.3mV 21.0°C A				
02:31:20 PM Feb 24, 2010				
100 ppm, 31.2mV 21.0°C A				
02:31:37 PM Feb 24, 2010				
1000 ppm, 61.1mV 21.1°C A				
02:36:17 PM Feb 24, 2010				
	Escape			

7.3.2.10 Electrode Type

Select the Ion Selective Electrode used for measurements from a list: Ammonia, Bromide, Cadmium, Calcium, Carbon Dioxide, Chloride, Cupric, Cyanide, Fluoride, Iodide, Lead, Nitrate, Potassium, Silver, Sodium, Sulfate, Sulfide or five custom ISEs. For the standard ISE it is possible to view the ion constants (Name, Molar Weight and Electric Charge/Slope), while for the custom ISE all these constants must be manually set.



For Standard ISE:

- Press to see the ion constants, press at any time to exit Ion Constants view.

For Custom ISE:

- Press to edit the ion constants for the selected custom ISE. Use the and keys to select the desired Ion constant and press to edit the value or to cancel operation.
- Set the Ion Name (up to 10 characters can be entered).
- Set the appropriate molecular weight (in g / mol) using the numeric keys. Press to save the value or press to return to the previous screen.
- Select the appropriate Electric Charge / Slope. Use the and keys to select the value and then press . If the Ion electric charge is None, its slope can be manually set by pressing . Press to save the value or press to return to the previous screen.

pH, mV & ISE MODE

7.3.2.11 Concentration Unit

Select the desired concentration unit for the measured ion or chemical compound. The available concentration units are: ppt (g/L), ppm (mg/L), ppb ($\mu\text{g/L}$), mg/mL, M (mol/L), mmol/L, %w/v or user defined.

Concentration Unit				
Select a menu option.				
ppt				
g/L				
ppm				
mg/L				
ppb				
$\mu\text{g/L}$				
mg/mL				
M				
mol/L				
mmol/L				
%w/v				
User				
Select	Escape			

7.3.3.12 Logging Interval


Set the logging interval to be used.

Logging Interval				
Enter the data logging interval.				
0				
hours	minutes	seconds		
Press Next to move to the next entry.				
Accept	Escape	Delete Digit	Next	Off

7.3.3.13 Stability Criteria


This option allows the user to select the signal stability criteria for the measured parameters:

- Fast quicker results with less accuracy
- Medium medium speed results with medium accuracy
- Accurate slower results with high accuracy



7.3.3.14 ISE Significant Digits

Select the number of significant digits to be displayed: one(x), two(xx) or three(xxx).



pH, mV & ISE MODE

7.3.2.15 Stirrer Configuration

Set the stirrer configuration: Stirrer 1, Stirrer 2 or Disabled.

Stirrer Configuration							
Select the option.							
<table border="1"><tr><td>Disabled</td></tr><tr><td>Stirrer 1</td></tr><tr><td>Stirrer 2</td></tr></table>					Disabled	Stirrer 1	Stirrer 2
Disabled							
Stirrer 1							
Stirrer 2							
Select	Escape						

7.3.2.16 Stirring Speed

The stirring speed for the selected stirrer can be set.

Stirring Speed					
Enter the speed of the stirrer within below range.					
<table border="1"><tr><td>1100</td></tr></table> RPM					1100
1100					
The range is from 0 to 2500 RPM.					
Accept	Escape	Delete Digit			

7.3.3 ISE Calibration

It is recommended to calibrate the instruments frequently if high accuracy is required. The instrument should also be recalibrated whenever the "Calibrate Electrode" message appears on the LCD.

Due to electrode conditioning time, the electrode must be immersed for several seconds to stabilize. The user will be guided step by step during calibration with easy-to-follow messages on the display. This will make the calibration a simple and error-free procedure.

PREPARATION:

Pour small quantities of the standard solution into clean beakers. If possible, use plastic beakers to minimize any EMC interferences.

For accurate calibration and to minimize cross-contamination, use two beakers for each standard solution: one for rinsing the electrode and one for calibration.

Note: For accurate measurements, add the appropriate ISA (Ionic Strength Adjustment) to the calibration standards.

CALIBRATION PROCEDURE:

Before calibrating, make sure that the appropriate Electrode Type and concentration unit has been selected in ISE Setup.

Up to a five point calibration is possible using any combination of five memorized standard solutions and five custom solutions.

The ISE calibration and measurement can be performed with or without temperature compensation. If the temperature compensation option is enabled, the isopotential point of the electrode must be set in ISE Setup.

The current standard will be manually selected by the user from the available standards list. The list of available standards depends of the Manual Entry setting.

To calibrate the instrument using Manual Entry:

- Press from the main screen. If the instrument was calibrated before and the calibration was not cleared, the old calibration can be cleared by pressing .

The screenshot shows the 'ISE Calibration' screen. At the top, it displays 'Analogs 2' and a large reading of '1010 PPM'. Below this, it shows 'ISE: Calcium' and 'Stirrer 2' with a speed of '1200/1200 RPM'. The temperature is 'Temp 21.0 °C' and the mV reading is '61.1 mV'. A 'Calibrated Standards' list shows two beakers with '10.0' and '100'. The 'Last Calibration' is '02:31 PM Feb 24, 2010'. A prompt says 'Press <Accept> to update calibration.' At the bottom, there are five buttons: 'Accept', 'Escape', 'Edit', 'Next Standard', and 'Previous Standard'.

- Immerse the Ion Selective Electrode and the temperature probe approximately 2 cm into the lowest concentrated standard solution.
- Select the concentration with or .
- When the reading has stabilized, press to update the calibration. The calibration point value will be added to the Calibrated Standard list.
- Select and repeat the procedure with all of the available standards.

pH, mV & ISE MODE

7.4 Logging

Data logging is available in pH, mV or ISE mode. It can be a single data point or a continuous timed log.

To customize the logging report:

- Press  to display the **Data Parameters** screen.
- Highlight the *Setup pH/mV/ISE Report* option and press  to display the **Setup pH/mV/ISE Report** screen.
- Use the  and  keys to highlight the data field that you want to show/hide in the pH/mV/ISE report and then press  to activate/deactivate it.
- Each field marked by "*" is an active field selected for the report.
- Press  to save the customized report.

7.4.1 Automatic Logging

The logging interval is set in the pH / mV / ISE setup screen.

Press  to start the log.

The logging interval and name of logging file will be also displayed on the measure screen.

To stop the automatic logging, press  again.

7.4.2 Manual Logging

To manually log pH, mV or ISE readings, press  from the **pH, mV** or **ISE** screen.

A new record will be added to the report every time  is pressed.

Chapter 8. Contents

8	AUXILIARY FUNCTIONS	8 - 3
8.1	Burette	8 - 3
8.1.1	Prime Burette	8 - 4
8.1.2	Rinse Tip	8 - 5
8.1.3	Manual Dispense	8 - 5
8.1.4	Purge Burette	8 - 5
8.2	Stirrer	8 - 6
8.3	Results	8 - 6
8.3.1	Review Last Titration Report	8 - 7
8.3.2	Review Available Reports	8 - 8
8.3.3	GLP Data	8 - 8
8.3.4	Instrumentation Data	8 - 9
8.3.5	Logging Interval – pH/mV/ISE	8-10
8.3.6	Setup pH/mV/ISE Report.....	8-10
8.3.7	Setup Titration Report.....	8-10

8 AUXILIARY FUNCTIONS

8.1 Burette

To access the **Burette** screen, press Burette from the main titration screen. Highlight the desired option and then press Select.

Burette								
Select a menu option.								
<table><tr><td style="background-color: black; color: white;">Prime Burette</td></tr><tr><td>Rinse Tip</td></tr><tr><td>Manual Dispense</td></tr><tr><td>Purge Burette</td></tr></table>					Prime Burette	Rinse Tip	Manual Dispense	Purge Burette
Prime Burette								
Rinse Tip								
Manual Dispense								
Purge Burette								
The current pump is: Pump 1 Current burette volume is 25 mL.								
Select	Escape	Choose Pump						

Choose Pump allows you to select the desired pump for burette operations (it is only active if two pumps are connected).


Pump Setting						
Select the current pump.						
<table><tr><td style="background-color: black; color: white;">Pump 1</td></tr><tr><td>Pump 2</td></tr></table>					Pump 1	Pump 2
Pump 1						
Pump 2						
Select	Escape					

AUXILIARY FUNCTIONS

8.1.1 Prime Burette

The *Prime Burette* option is used to mechanically fill the burette before starting a set of titrations. The priming process consists of several cycles of filling and emptying the burette with titrant.

Two rinse cycles of burette are shown in the figure below. The dispensing tube is connected on the right side and the aspiration tube on the left side.



Note: Before starting this operation, the aspiration tube must be inserted in the titrant bottle. A waste container should be placed under the dispensing tip to collect the waste solution.

To prime the burette, select *Prime Burette* from the **Burette** screen. Enter the number of rinses and press .

The number of burette rinses can be set between 1 and 5 (we recommend at least three rinses to assure that the air bubbles are completely removed).

Total Burette Rinses				
Enter the total number of burette rinses.				
3				
A minimum of three rinses is recommended.				
Accept	Escape	Delete Digit		

8.1.2 Rinse Tip

A 2 mL dose of titrant will be dispensed from the burette when this operation is selected. This operation will eliminate the air from the dispensing tip.

8.1.3 Manual Dispense

Manual Dispense option allows a defined titrant volume to be dosed. Select the *Manual Dispense* option and press . The **Manual Volume Dispense** screen will become active and the display will prompt you to enter the desired volume to be dispensed.

Manual Volume Dispense				
Enter the amount of volume to be dispensed.				
1.000 mL				
Current burette volume is 10 mL.				
ACCEPT	Escape	Delete Digit		

The manual dispense volume must be between the limits shown below:

- 0.001 to 4.500 mL for a 5 mL burette
- 0.001 to 9.000 mL for a 10 mL burette
- 0.005 to 22.500 mL for a 25 mL burette


8.1.4 Purge Burette

This option allows the burette to be emptied before cleaning and/or storing the burette. The burette is flushed twice.

Note: Before starting this operation, remove the aspiration tube from the titrant bottle.

AUXILIARY FUNCTIONS

The figures below show the steps in a purge burette operation.



8.2 Stirrer


The stirrer can be turned on and off by pressing .

The stirring speed is set within the method parameters (see section 5.5.17 *Stirring Speed*).

During the titration process, the stirring speed can be manually adjusted by using the and keys.

8.3 Results

From the **Data Parameters** screen you can access the following options:



8.3.1 Review Last Titration Report

The last titration report can be reviewed.

The titration graph can be reviewed by selecting .

Review Result				
HI902 - Titration Report				
Method Name: copy of Titr. Conc. 0.1N				
Time & Date: 02:06 PM Feb 24, 2010				
Titration ID: Ti_00022				
Titration Results				
Method Name: copy of Titr. Conc. 0.1N				
Time & Date: 02:06 PM Feb 24, 2010				
Analyte Size: 0.200 g				
End Point Volume: 0.985 mL				
pH Equivalence Point: 6.272				
Results: 9.9395E-1 N (eq/L)				
View Graph	ESCAPE	Print Report	Page Up	Page Down

The information seen in the report is based on the selections made in the **Setup Titration Report** screen.

The following option keys are available:

 Review the titration graph. The potentiometric titration curve is displayed. If the *Equivalence End Point* option was selected, the derivative curve (1st derivative, 2nd derivative) is simultaneously displayed. Pressing  will change the vertical axes scale units.

 Print the titration report.

 Return to the previous screen.


  Keys can be used to scroll through the pages.

AUXILIARY FUNCTIONS

8.3.2 Review Available Reports

Up to 100 reports can be saved on the titrator. To view one of the saved reports, highlight a report and then press .

All of the saved reports can be reviewed and printed.




The report contains only the information selected in the **Setup Titration Report** and **Setup pH/mV/ISE Report** screens during report configuration.

The following option keys are available:

-  Review the selected graph.
-  Review the selected report.
-  Print the selected report.
-  Delete the selected report.
-  Return to the previous screen.

8.3.3 GLP Data

Enter up to 20 alphanumeric characters for each option from **GLP Data** screen.



- Sample Name* Allows the sample name to be recorded in each report. The sample name will increase by one, with each new titration or logging report, if the last character is a number.
- Company Name* Allows the company name to be recorded in each report.
- Operator Name* Allows the operator name to be recorded in each report.
- Electrode Name* Allows the electrode name to be recorded in each report.
- Fields 1, 2, 3* Allows any additional information to be recorded in each report. The fields must be selected from **Setup Titration Report** screen (see section 8.3.6 *Setup pH/mV/ISE* and section 8.3.7 *Setup Titration Report*) in order to be displayed in the titration report.

8.3.4 Instrumentation Data

Displays titrator configuration data.

```

Meter Information
HI 902C Titrator

SERIAL NUMBER
Titrator Serial Number:      02103805
Analog Board 1 Serial Number: 30103805
Analog Board 2 Serial Number: 30103806
Pump 1 Serial Number:        70103805

SOFTWARE VERSION
Titrator Software Version:    v2.00
Base Board Software Version:  v2.01
Pump 1 Software Version:      v1.4

Analog 1 Calibration Date:    Oct 18, 2010
Analog 2 Calibration Date:    Oct 18, 2010
    
```

	Escape	Print		
--	--------	-------	--	--

Titrator Serial Number: The serial number of the titrator base board.

Analog Board 1 (and/or 2) Serial Number: The serial number of the analog board.

Pump 1 (and/or 2) Serial Number: The serial number of the connected pump.

Titrator Software Version: The current software version installed on the titrator.

Base Board Software Version: The current software version present on the base board of the titrator.

Pump 1 (and/or 2) Software Version: The current software version for the pump.

Analog 1 (and/or 2) Calibration Date: Manufacturer calibration date of the analog board.

Note: If more than 1 year elapsed from the calibration date of the analog board 1 and/or 2, the message **Analog 1 Calibration Due** and/or **Analog 2 Calibration Due** will appear on the main screen. The analog board(s) need to be recalibrated.

AUXILIARY FUNCTIONS

8.3.5 Logging Interval – pH/mV/ISE

Allows the user to setup a time interval for auto-logging while in pH, mV, or ISE mode (see section 7.3 Logging).

8.3.6 Setup pH/mV/ISE Report

Customize a unique report to record the pH, mV, and ISE measurements. An asterisk means that it will be included in the report.

Setup pH/mV/ISE Report				
Select fields to be saved in the report.				
<ul style="list-style-type: none">* Result and Units* Potential* Temperature and Units* Date and Time* Calibration DataSample NameCompany NameOperator NameElectrode NameField 1Field 2Field 3Software VersionsSerial Numbers				
Select	Escape	Save Report		

8.3.7 Setup Titration Report

Customize a unique report to record the titration results. An asterisk means that it will be included in the titration report.

Setup Titration Report				
Select fields to be saved in the report.				
<ul style="list-style-type: none">* Result and Units* Titration Method* Initial and Final Readings* Analyte Size* End Point Volume* Titration Duration* Date and Time* Titration Ended ByAll Data PointsMethod ParametersStandardization DataSample NameCompany NameOperator Name				
Select	Escape	Save Report	Page Up	Page Down

Chapter 9. Contents

9	MAINTENANCE, PERIPHERALS	9 - 3
9.1	Burette Maintenance	9 - 3
9.1.1	Burette Assembly	9 - 3
9.1.2	Changing the Burette	9 - 3
9.1.3	Disassembling the Burette	9 - 4
9.1.4	Assembling the Burette	9 - 4
9.1.5	Cleaning the Burette.....	9 - 5
9.1.6	Burette Preparation (Filling with Titrant)	9 - 7
9.2	Fuses Replacement	9 - 8
9.3	Peripherals.....	9 - 9
9.3.1	Connecting an External Display	9 - 9
9.3.2	Connecting an External PC Keyboard	9 - 9
9.3.3	Connecting a Printer	9 - 11
9.3.4	Connecting to a Computer	9 - 11

9 MAINTENANCE, PERIPHERALS

The 25 mL burette included with the titrator exceeds the ISO 8655 standard for accurate delivery of liquids by a motor-driven piston burette.


9.1 Burette Maintenance

9.1.1 Burette Assembly

The burette is delivered with a 25 mL syringe inside and with all of the accessories mounted (see section *2.1 Unpacking* for burette assembly details). The burette assembly consists of a rigid housing which holds the glass syringe, a 3/2-way valve, and titrant tubing.

9.1.2 Changing the Burette

Remove the burette from the pump assembly by sliding it forward and then slide the new burette into place (see the picture below).




MAINTENANCE, PERIPHERALS

9.1.3 Disassembling the Burette

The aspiration and the dispensing tubes have fittings and tube protectors. The aspiration tube will be mounted in the left side and the dispensing tube will be mounted in the right side of the burette.

To remove the dispensing tube and the aspiration tube follow these steps:


- Slide up the tube protector (3);
- Remove the tube lock (4) from the burette holder.
- Unscrew the fitting (2).
- Remove the tube.



9.1.4 Assembling the Burette

To attach the dispensing tube and the aspiration tube, follow these steps:

- Insert the flat-shaped end of the dispensing tube into the valve outlet, and screw in the fitting so that the highest of its 9 cuts stays vertically in the final position (2).
- Bend the tube up into the vertical position to enter the highest cut of the fitting (5).
- Put on the tube lock on the fitting (4).
- Slide down tightly the tube protector (3) into the dedicated gap of the tube lock.



9.1.5 Cleaning the Burette

To clean the burette, follow these steps:

- If the burette is filled with titrant, remove the aspiration tube from the titrant bottle and purge burette (see section 8.1.4 *Purge Burette*).
- Insert the aspiration tube into cleaning solution, deionized water or titrant solvent.
- Prime burette to fill the burette (use 2 rinses) (see section 8.1.1 *Prime Burette*).
- During second refilling of the burette remove the aspiration tube out of the cleaning solution, deionized water, or solvent and allow the air to replace the liquid in the burette. This will clean the aspiration tube.

If this simple cleaning procedure is not adequate, continue with these steps:

- Slide the burette out from the pump assembly.
- Remove the dispensing and aspiration tubes. Clean them separately or insert new ones.
- Remove the protective cap from the bottom of the burette assembly by using the special tool.
- Remove the syringe from the burette assembly by unscrewing it with your fingers.
- Extract the piston from the syringe.
- Clean both the piston and the syringe with appropriate cleaning solution. Rinse with deionized water.
- Remove the excess liquid.

MAINTENANCE, PERIPHERALS


Warning: Avoid contacting the titrant with bare hands.

Avoid spilling titrant.

Clean the external side of the syringe and piston to remove aggressive chemicals.

Do not touch the Teflon part of the piston or internal walls of the burette with bare hands or greasy materials.

- Reinsert the piston into the syringe.
- Reinsert the syringe by screwing it in the valve with your fingers.
- Reinsert the protective cap to the bottom of the burette assembly. Carefully position the cap into the burette.
- Slide the burette into the burette stand. Notice the position of the piston shaft to the pump couple.
- Priming the burette three times with new titrant is recommended.



9.1.6 Burette Preparation (Titrant Filling)

Before starting a titration, the burette must be properly filled with titrant in order to obtain an accurate and repeatable result. To fill the burette, follow the next steps and recommendations:

- If necessary, clean the burette and make sure it is empty.
- From the main screen press .
- Highlight *Prime Burette* option and press .
- Enter the number of times the burette needs to be rinsed (minimum three rinses are recommend allowing air bubbles to be evacuated).
- Press .
- Insert the aspiration tube into the titrant bottle only when the piston is going down and has reached about $\frac{1}{4}$ from the upper end.

To avoid the presence of the air bubbles inside the burette, make sure to have a continuous liquid flow inside the burette. A little air just above the liquid level at the first filling is normal. The next filling will evacuate all of the air; no air will be left in the valve.

Sometimes during this process, slight finger tapping on the tubes is helpful to remove any residual air bubbles from the tubes.

If air bubbles are still present:


- Remove the aspiration tube from the titrant bottle.
- Repeat burette preparation procedure.
- If no success, clean the burette again.

MAINTENANCE, PERIPHERALS

9.2 Fuse Replacement

To replace the fuses, follow these steps:

- Turn off the titrator.
- Remove the power cord from the power connector located on the rear side of the titrator case.
- With a screw driver open the fuse holder lid.



- Extract the fuse holder.
- Replace the fuses (for recommended fuses replacement please see the label located above the power switch).
- Close the fuse holder lid.
- Connect the power cord.

Note: For other maintenance operations please contact your dealer or the nearest Hanna Customer Service Center.


9.3 Peripherals

Warning! Connection/disconnection of POWER CORD, PUMP ASSEMBLY, EXTERNAL PC DISPLAY, PRINTER, RS232 INTERFACE, EXPANSION DEVICE must only be done when titrator and external devices are turned off.

MAINTENANCE, PERIPHERALS

9.3.1 Connecting an External Display


The information shown on the titrator display can be viewed also on a Standard VGA display connected with a 15-pin cable, as presented below.



Connect the external display to the display socket.
Turn on the titrator and then the external display.

9.3.2 Connecting an External PC Keyboard



This connection allows you to use an external PS/2 PC Keyboard in addition to the titrator's keypad.



Connect an external PC Keyboard (PS/2 connector).


MAINTENANCE, PERIPHERALS

The correspondence between the Titrator's Keypad and the United States 101-type external keyboard are:

External PC Keyboard (United States 101)	Titration Keypad
Function Key F-1	
Function Key F-2	
Function Key F-3	
Function Key F-4	
Function Key F-5	Option Key 1 (from left to right)
Function Key F-6	Option Key 2 (from left to right)
Function Key F-7	Option Key 3 (from left to right)
Function Key F-8	Option Key 4 (from left to right)
Function Key F-9	Option Key 5 (from left to right)
Function Key F-10	
Arrow Key: Up	
Arrow Key: Down	
Arrow Key: Left	
Arrow Key: Right	
Page Up	
Page Down	
Numeric Keys: 0 to 9	
Tab	
Enter	
Alphanumeric Keys	Allow alphanumeric entries.

9.3.3 Connecting a Printer

A variety of parallel printers can be connected to the parallel port of the titrator using a standard DB25-pin cable.




Warning: The titrator and the external printer must be both turned OFF before they are connected.

Connect the external printer to the standard 25-pin Socket.
Turn on the titrator and then the printer.

9.3.4 Connecting to a Computer

The titrator can be connected to a computer using a USB cable. HI 900 PC application needs to be installed on the PC.



MAINTENANCE, PERIPHERALS


Connect the cable to the USB port on the rear panel of the titrator.

Connect the cable to the USB port on the PC.

Select the **USB Communication** screen on the titrator following the path:

General Options - USB Link with PC

Launch the HI 900 PC application and then select the appropriate USB Port on the PC.



The HI 900 PC application allows the transfer of methods and reports between the titrator and PC.

Appendix 1. Contents

A1 TECHNICAL SPECIFICATIONS A1-3

A1 TECHNICAL SPECIFICATIONS

mV	Range	- 2000.0 to 2000.0 mV
	Resolution	0.1 mV
	Accuracy	±0.1 mV
pH	Range	- 2.000 to 20.000 pH
	Resolution	0.1 / 0.01 / 0.001 pH
	Accuracy	±0.001 pH
ISE	Range	1x10 ⁻⁶ to 9.99x10 ¹⁰
	Resolution	1 / 0.1 / 0.01
	Accuracy	±0.5% (monovalent ion) ±1.0% (divalent ion)
Temperature	Range	- 5.0 to 105.0 °C 23.0 to 221.0 °F 268.2 to 378.2 K
	Resolution	0.1 °C / 0.1 °F / 0.1 K
	Accuracy	±0.1 °C / ±0.4 °F / ±0.2 K
Burette Sizes	Resolution	0.001 mL
	Accuracy	±0.005 mL (5 mL Burette) ±0.010 mL (10 mL Burette) ±0.025 mL (25 mL Burette)

Graphic Display 7.5" graphical color display with backlight.

Languages English, Portuguese, Spanish.

Titration Methods up to 100 (standard and user methods)

Burette size auto-detection and interchangeable burettes. The titrator automatically detects the size of the burette when it is slid in the pump assembly.

Propeller Stirrer with Programmable Stir Speed. The stirring speed can be set between 100 and 2500 RPM with 100 RPM resolution.

Flow Rate: user-selectable (see Section 5.5.22 *Volume/Flow Rate*).

mV / pH / ISE Measurement modes.

Automatically Temperature Compensated pH Measurements.

pH Calibration with up to 5 buffers using *Auto-Entry* or *Manual-Entry* options; temperature compensated buffers are stored internally for *Auto-Entry* option.

Relative mV calibration: single point offset.

ISE Calibration: with up to 5 standards.

APPENDIX 1

Potentiometric Titrations: Acid-Base (pH or mV-Mode), Redox, Precipitation, Complexometric, Non-Aqueous, Ion-Selective, Argentometric.

Titer Determination.

Fixed mV or pH End Point Detection.

Single Equivalence Point Detection, with the 1st or 2nd Derivatives of the titration curve.

Multiple Equivalence Point Detection

Flexible Concentration Calculations, with many concentration units.

Graph Display during titration, graphs of the stored titration data (mV-Volume or pH-Volume titration curve, 1st derivative curve or 2nd derivative curve, in pH-mode or mV-mode) and pH/mV values versus time-data logging results.

Data Storage: up to 100 complete titration report and pH/mV/ISE reports.

Files Copied to and Restored from USB Storage Device: Standard Methods, User Methods, Titration and pH/mV/ISE Logging Reports, and Bitmap Files can be transferred to a PC using a USB storage device.

Peripheral Units:

External VGA Display

External PC Keyboard

Printer

PC via USB Interface

GLP Conformity: Good Laboratory Practice and Instrumentation Data storage and printing capabilities.

Mains: 110/220 Vac ; 50-60 Hz

Power Draw: 1.3 Amps

2 Exchangeable Fuses

Enclosure Material: ABS plastic and Steel

Keypad: Polycarbonate

Dimensions: Width x Depth x Height = 390 x 350 x 380 mm

Weight: approx. 22 lbs. (10 Kg) (with 1 pump, stirrer and sensors)

Operating Environment: 10 to 40 °C, up to 95% relative humidity

Storage Environment: -20 to 70 °C, up to 95% relative humidity

Appendix 2. Back Titration

A2	BACK TITRATION	A2-3
A2.1	Applicability Domain	A2-3
A2.2	Method Principles	A2-3
A2.3	Example of a Back Titration	A2-3
A2.3.1	Introduction	A2-3
A2.3.2	Setting Up the Method Parameters	A2-4
A2.3.3	Preparing the Sample	A2-5
A2.3.4	Perform the Titration	A2-5

A2 BACK TITRATION

A2.1 Applicability Domain

Back titrations are generally used for one of the following reasons:

- Reaction kinetics are too slow for the direct titration of the analyte;
- The metal precipitates at the desired pH (complexometric titrations);
- The reaction between titrant and analyte produces some auxiliary compounds (e.g. CO₂) that can affect the mV signal and also the equivalence point detection;
- Titrations with flat first derivative curves.

A2.2 Method Principles


An excess of reagent is added to the sample solution, helping a slow reaction go to completion. The unreacted excess reagent is then titrated. The difference in the total volume of the first reagent added and the amount determined from the second titration is the quantity of reagent required to complete the reaction. The sample concentration is calculated using this value.

A2.3 Example of a Back Titration


An example of a back titration is determining the neutralizing capacity of antacid.

A2.3.1 Introduction

An excess of stomach acid (primarily HCl) causes heartburn and acid indigestion. Commercial antacids active ingredient is a basic salt such as Mg(OH)₂ (milk of magnesia), NaHCO₃ (sodium bicarbonate), CaCO₃ (calcium carbonate) or Al(OH)₃ (aluminum hydroxide). In this example we will analyze a typical antacid containing CaCO₃, which reacts with the acid to form a salt, water and gas:



A conventional acid/base titration is very difficult because the active ingredient is sparingly soluble in water, and the CO₂ gas has an influence in pH measurements. In order to overcome this the tablets will be completely dissolved in excess acid:



The excess acid will be titrated with NaOH to determine the amount of CaCO₃ present.



The final concentration is expressed as g CaCO₃/g.

APPENDIX 2

A2.3.2 Setting Up the Method Parameters

Use the following settings to program the method (see section 5.5 *Method Options* section):

Name:	Antacid Power
Method Revision:	1.0
Titration Type:	Back-Titration
Break at titrant change:	YES
Analog Board:	Analog 1
Stirrer Configuration:	Stirrer 1
Pump Configuration:	
Titrant 1 pump:	Pump 1
Titrant 2 pump:	Pump 2
Dosing Type:	Dynamic
min Vol:	0.010 mL
max Vol:	0.500 mL
delta E:	4.500 mL
End Point Mode:	pH 1EQ point, 1st Der
Recognition Option:	
Threshold:	500 mV / mL
Range:	NO
Filtered Derivatives:	NO
Pre-Titration Volume:	0.000 mL
Pre-Titration Stir Time:	5 Sec
Measurement Mode:	Signal Stability
delta E:	0.3 mV
delta t:	1.5 Sec
t min wait:	5 Sec
t max wait:	15 Sec
Electrode Type:	pH
Calculations:	Sample Calc. by Weight
Titrant 1 Units	M (mol/L)
Titrant 2 Units	M (mol/L)
Final Result Units	g/g
Titrant 1 Name:	HCl
Titrant 1 Conc.:	1 M(mol/L)
Titrant 2 Name:	NaOH
Titrant 2 Conc.:	1 M(mol/L)
Analyte Size:	0.5 g
Analyte Entry:	Manual
Titrant 1 Entry:	Calculated
Maximum Titrant 2 Volume:	25.000 mL
Stirring Speed:	1500 RPM
Potential Range:	-2000.0 to 2000.0 mV
Flow Rate:	30.0 mL/min
Signal Averaging:	1 Reading
Significant Figures:	XXX
Linked to:	No link

In order to start the analysis, both burettes must be filled with the appropriate reagents:

- Pump 1 - a burette filled with 1 M HCl.

- Pump 2 - a burette filled with 1 M NaOH.


A2.3.3 Preparing the Sample

For this example commercial antacid pills were used.

- Crush several pills with a mortar and pestle. From the crushed pills weight roughly 0.5 g into a 150 mL sample beaker.
- Pour about 20 mL distilled and deionized water into the beaker.
- Raise the stirrer assembly.
- Place the beaker under the stirrer assembly.
- Lower the stirrer assembly until it rests on its positioning collar.
- Adjust the height of the stirrer so it is as close as possible to the bottom of the beaker.
- If necessary add more distilled or deionized water so that the pH electrode bulb is completely immersed in the sample and the reference junction of the electrode is 5-6 mm below the surface.

A2.3.4 Perform the Titration

- When the method is selected (the main screen displays "Antacid Power") press .
- The titrator will prompt for the sample weight. Enter the exact weight value (with 4 digits) and press .
- The calculated volume of titrant 1 is displayed (this value can be modified by the user). The displayed volume will be dispensed during the first phase of back-titration.
- Press to proceed with the next step.
- The titrator will start to dispense titrant 1 (1 M HCl).



APPENDIX 2

- When this phase is completed, the titrator will stop and the ***Continue Back Titration*** screen is displayed.
- Slide the stirrer assembly up.
- Put the beaker on a hotplate.
- Heat gently until all of the effervescence has ceased, then boil it for 1-2 minutes. Some of the inactive tablet material may not dissolve; however, this should not interfere with the titration.
- Cool the solution to room temperature and put the beaker under the stirrer assembly.
- Lower the stirrer assembly until it rest on its positioning collar.
- Press to proceed to the second phase of the analysis.
- The results will be displayed on the screen.

Appendix 3. Multiple Equivalence Point Titration

A3	MULTIPLE EQUIVALENCE POINTS	A3-3
A3.1	Applicability Domain	A3-3
A3.2	Performing a Titration	A3-3
A3.3	Example of a Multiple Equivalence Point Titration	A3-4

A3 MULTIPLE EQUIVALENCE POINTS

A3.1 Applicability Domain

Multiple equivalence point titrations can be used to determine:

- the component concentrations from a synthetic mixture (e.g.: Mixture of HCl, CH₃COOH, NH₄Cl);
- the concentration of a polyprotic acid in its titratable ionization stages (e.g. H₃PO₄, two equivalence points in aqueous medium).

A3.2 Performing a Titration

Set up the method parameters as follows:

- Set the end point mode to equivalence end point (mV or pH).
- Set the number of equivalence points, up to 5 points can be detected.

Number of Equivalence Points				
Enter the number of equivalence points to be found.				
3 points				
The range is between 1 and 5 equivalence points.				
ACCEPT	Escape	Delete Digit		

- Select the end point determination, first or second derivative.
- The rest of the method parameters can be modified (if needed) and save the method.

Note: In the method calculation screen the molecular weight and reaction ratio can be entered for each point.

APPENDIX 3

A3.3 Example of a Multiple Equivalence Point Titration

An example of a multi-equivalence point titrations is titrating a mixture of HCl, CH₃COOH and NH₄Cl with Sodium Hydroxide.

During this type of titration the user is able to view the following information:

- After first equivalence point is detected, the titration screen will show  and the


16:17:25 Jun 21 2006			
HCl + NH4Cl + CH3COOH			
3.1328E-2 N (eq/L)			
EQ points detected: 1		1500/1500 RPM	
In Progress			
Pump 1 Running			
Burette: 25 mL			
Unstable			
ATC	Volume Delivered	mV	
35.2 °C	0.727 mL	181.6	
View Details	View Curve	Suspend	Stop

number of equivalence points detected.

EQ Points Report			
	mV	Volume [mL]	Results [N (eq/L)]
EQ1:	228.4	0.393	3.6710E-2
Escape			

By pressing  the titrator will display **EQ Points Report**. This screen contains information about the detected equivalence points.

- On the graph, the equivalence points are marked with an "x" and titrant volume is displayed.
- For the first equivalence point the volume of titrant is equal to the volume of titrant dispensed.



- The volumes of each additional point (marked with an "x") are calculated as the difference between the total volume dispensed to reach the current equivalence point and the sum of the previously detected equivalence points.

The concentration is calculated with the formula selected in the *Calculations* option from the **View / Modify Method** screen.

Appendix 4. Linked Method

A4	LINKED METHOD	A4-3
A4.1	Creating a Linked Method	A4-3
A4.2	Selecting a Linked Method	A4-6
A4.3	Running a Linked Titration Method	A4-7

APPENDIX 4

- Press to link the highlighted method to the current method or to return to the previous screen.

Link Titration Method				
Select the method to be linked.				
----- No Link USER0006 copy of Alkalinity of Wa USER0007 copy of Acidity of Water USER0008 Acidity				
Select	Escape			

- The second titration can be started either "automatically" or "manually". If *Automatically* is selected, the second titration will start immediately after the first titration has completed. No user input is required. If *Manually* is selected, user input is required to start the second method.

View / Modify Method				
Id: USER0009 Modified: Jun 17, 2010 13:30				
Select the option to be modified.				
Calculations: Sample Calc. by Weight Dilution Option: Disabled Titrant Name: 0.02M AgNO3 Titrant Conc.: 2.0000E-2 M (mol/L) Analyte Size: 1.000 g Analyte Entry: Manual Maximum Titrant Volume: 25.000 mL Stirring Speed: 1400 RPM Potential Range: -2000.0 to 2000.0 mV Volume/Flow Rate: 25 mL / 50.0 mL/min Signal Averaging: 1 Reading Significant Figures: XXXXX Linked To: Acidity Start Linked Method: Automatically				
Select	Escape	Print Method	Method 2	Page Down

- To access the method options for the second method, press Method 2.

```

View / Modify Method
Id: USER0008 Modified: Jun 17, 2010 10:33
Select the option to be modified.
Name: Acidity
Method Revision: 1.0
Analog Board: Analog 2
Pump Configuration:
Dosing Type: Dynamic
End Point Mode: Fixed 8.300 pH
Pre-Titration Volume: 0.000 mL
Pre-Titration Stir Time: 0 Sec
Measurement Mode: Signal Stability
Electrode Type: pH
Blank Option: No Blank
Calculations: Sample Calc. by Volume
Dilution Option: Disabled
Titrant Name: 0.1N NaOH
  
```

Select	Escape	Print Method	Method 1	Page Down
--------	--------	-----------------	----------	--------------

- The analyte entry for Method 2 depends on the type of calculation. If both methods are using the same type of calculation (sample calculation by volume or by weight), the option for *Same as Previous* is available for Method 2. This allows the same weight or volume to be used in both calculations. *Same as Previous* is not available if different types of calculations are being used.

```

Analyte Entry
Select the entry mode of analyte.
Fixed Weight or Volume
Manual Weight or Volume
Same as Previous
Verify the correct formula is being used,
I.E. weight or volume analyte type.
  
```

Select	Escape			
--------	--------	--	--	--

APPENDIX 4

- Once the linked method has been set up, press Escape to save the new method and to return to the main screen.


13:31:30 Jun 17, 2010				
Analog 1				
Chloride				
Linked To: Acidity				
0.0 meq/Kg				
No Result				
Pump 1 Selected				
End Point Volume				
Manual	25.0 °C	0.000 mL		mV 125.4
General Options	Select Method	Method Options	Burette	Mode

A4.2 Selecting a linked method

All linked methods are moved to the top of the user method list and are noted with "*" next to the method number. When a method is linked, only the first method can be selected and the second method is indented below it.


Titration Methods																																
Select the method to be activated.																																
<table><tr><td>HI0001EN</td><td>0.1N NaOH Titr. Conc.</td></tr><tr><td>HI0002EN</td><td>0.1N HCl Titr. Conc.</td></tr><tr><td>HI0003EN</td><td>0.1M Na2S2O3 Titr. Conc.</td></tr><tr><td>HI0010EN</td><td>0.1M FAS Titr. Conc.</td></tr><tr><td>HI1004EN</td><td>Alkalinity of Water</td></tr><tr><td>HI1005EN</td><td>Acidity of Water</td></tr><tr><td>HI1007EN</td><td>Chloride in Water</td></tr><tr><td>HI1008EN</td><td>Neutralization w/ H2SO4</td></tr><tr><td>HI1009EN</td><td>Neutralization w/ NaOH</td></tr><tr><td>HI1011EN</td><td>Troubleshooting 1</td></tr><tr><td>HI1012EN</td><td>Troubleshooting 2</td></tr><tr><td>HI1014EN</td><td>Concentration of H3PO4</td></tr><tr><td>* USER0009</td><td>Chloride</td></tr><tr><td>USER0008</td><td>Acidity</td></tr></table>					HI0001EN	0.1N NaOH Titr. Conc.	HI0002EN	0.1N HCl Titr. Conc.	HI0003EN	0.1M Na2S2O3 Titr. Conc.	HI0010EN	0.1M FAS Titr. Conc.	HI1004EN	Alkalinity of Water	HI1005EN	Acidity of Water	HI1007EN	Chloride in Water	HI1008EN	Neutralization w/ H2SO4	HI1009EN	Neutralization w/ NaOH	HI1011EN	Troubleshooting 1	HI1012EN	Troubleshooting 2	HI1014EN	Concentration of H3PO4	* USER0009	Chloride	USER0008	Acidity
HI0001EN	0.1N NaOH Titr. Conc.																															
HI0002EN	0.1N HCl Titr. Conc.																															
HI0003EN	0.1M Na2S2O3 Titr. Conc.																															
HI0010EN	0.1M FAS Titr. Conc.																															
HI1004EN	Alkalinity of Water																															
HI1005EN	Acidity of Water																															
HI1007EN	Chloride in Water																															
HI1008EN	Neutralization w/ H2SO4																															
HI1009EN	Neutralization w/ NaOH																															
HI1011EN	Troubleshooting 1																															
HI1012EN	Troubleshooting 2																															
HI1014EN	Concentration of H3PO4																															
* USER0009	Chloride																															
USER0008	Acidity																															
Select	New Method	Delete	Page Up	Page Down																												

Use the arrow keys to highlight the method to be activated and press .



A4.3 Running a Linked Titration Method

Press to begin the linked titration.



After the completion of the first titration, the selected pump will return to the home position before the second titration is started. The titrator will either start the second titration automatically or wait for user input.

APPENDIX 4

<p>Start Linked Titration</p> <p>Analog 1 Method 1: Chloride</p> <p>6.9869 ppm (mg/L) <u>Titration Completed</u> End Point Volume 0.985 mL</p> <hr/> <p>Analog 2 Method 2: Acidity</p> <p><u>Manual</u> <u>Prepare To Run</u> mV 25.0 °C Pump 1 Selected 484.9</p> <p>Press <Abort> to cancel linked titration.</p>				<p>10:43:39 Jun 17, 2010</p> <p>Analog 2</p> <p>Acidity</p> <p>Linked By: Chloride</p> <p>10.009 mg/L</p> <p><u>Stirrer 1</u> 1400/1400 RPM</p> <p><u>In Progress</u></p> <p><u>Pump 1 Running</u> Burette: 25 mL Volume Delivered</p> <p><u>Manual</u> 25.0 °C 0.200 mL pH 3.853</p>					
Abort				View Curve		Suspend		Stop	

When the second titration has finished, the results for both titrations will be shown on the screen.

<p>10:43:53 Jun 17, 2010</p> <p>Analog 1 Method 1: Chloride</p> <p>6.9869 ppm (mg/L) <u>Titration Completed</u> End Point Volume 0.985 mL</p> <hr/> <p>Analog 2 Method 2: Acidity</p> <p>39.669 mg/L <u>Titration Completed</u> End Point Volume 0.793 mL</p> <p><u>Manual</u> 25.0 °C pH 8.383</p>				
ESCAPE	Method Options	View Report 1	View Report 2	

To start a new titration, press Escape to return to the main screen.

Appendix 5. Contents

A5	ACCESSORIES	A5-3
A5.1	Solutions	A5-3
A5.1.1	pH Calibration Solutions	A5-3
A5.1.2	pH Calibration Solutions in FDA Approved Bottle	A5-3
A5.1.3	pH Technical Calibration Solutions	A5-3
A5.1.4	pH Millesimal Calibration Solutions	A5-3
A5.1.5	Electrode Cleaning Solutions	A5-4
A5.1.6	Electrode Cleaning Solutions in FDA Approved Bottle	A5-4
A5.1.7	Electrode Storage Solutions	A5-4
A5.1.8	Electrode Storage Solutions in FDA Approved Bottle	A5-4
A5.1.9	Refilling Electrolyte Solutions	A5-4
A5.1.10	Refilling Electrolyte Solutions in FDA Approved Bottle	A5-5
A5.1.11	ORP Pretreatment Solutions	A5-5
A5.1.12	Titration Reagents	A5-5
A5.1.13	Ion Selective Electrode Calibration Solutions	A5-5
A5.2	Sensors	A5-6
A5.2.1	pH Electrodes	A5-6
A5.2.2	ORP Electrodes	A5-7
A5.2.3	Half-cell Electrodes	A5-8
A5.2.4	Ion Selective Electrodes	A5-8
A5.2.5	Temperature Sensor	A5-8
A5.3	Titration components	A5-9

A5 ACCESSORIES**A5.1 Solutions****A5.1.1 pH Calibration Solutions**

HI 7001M	—>	pH 1.68 Buffer Solution, 230 mL
HI 7001L	—>	pH 1.68 Buffer Solution, 500 mL
HI 7004M	—>	pH 4.01 Buffer Solution, 230 mL
HI 7004L	—>	pH 4.01 Buffer Solution, 500 mL
HI 7006M	—>	pH 6.86 Buffer Solution, 230 mL
HI 7006L	—>	pH 6.86 Buffer Solution, 500 mL
HI 7007M	—>	pH 7.01 Buffer Solution, 230 mL
HI 7007L	—>	pH 7.01 Buffer Solution, 500 mL
HI 7009M	—>	pH 9.18 Buffer Solution, 230 mL
HI 7009L	—>	pH 9.18 Buffer Solution, 500 mL
HI 7010M	—>	pH 10.01 Buffer Solution, 230 mL
HI 7010L	—>	pH 10.01 Buffer Solution, 500 mL

A5.1.2 pH Calibration Solutions in FDA Approved Bottle

HI 8004L	—>	pH 4.01 Buffer Solution, 500 mL
HI 8006L	—>	pH 6.86 Buffer Solution, 500 mL
HI 8007L	—>	pH 7.01 Buffer Solution, 500 mL
HI 8009L	—>	pH 9.18 Buffer Solution, 500 mL
HI 8010L	—>	pH 10.01 Buffer Solution, 500 mL

A5.1.3 pH Technical Calibration Solutions

HI 5016	—>	pH 1.68 Buffer Solution, 500 mL
HI 5003	—>	pH 3.00 Buffer Solution, 500 mL
HI 5004	—>	pH 4.01 Buffer Solution, 500 mL
HI 5068	—>	pH 6.86 Buffer Solution, 500 mL
HI 5007	—>	pH 7.01 Buffer Solution, 500 mL
HI 5091	—>	pH 9.18 Buffer Solution, 500 mL
HI 5010	—>	pH 10.01 Buffer Solution, 500 mL
HI 5124	—>	pH 12.45 Buffer Solution, 500 mL

A5.1.4 pH Millesimal Calibration Solutions

HI 6016	—>	pH 1.679 Buffer Solution, 500 mL
HI 6016-01	—>	pH 1.679 Buffer Solution, 1 L
HI 6003	—>	pH 3.000 Buffer Solution, 500 mL
HI 6003-01	—>	pH 3.000 Buffer Solution, 1 L
HI 6004	—>	pH 4.010 Buffer Solution, 500 mL

APPENDIX 5

HI 6004-01	—>	pH 4.010 Buffer Solution, 1 L
HI 6068	—>	pH 6.862 Buffer Solution, 500 mL
HI 6068-01	—>	pH 6.862 Buffer Solution, 1 L
HI 6007	—>	pH 7.010 Buffer Solution, 500 mL
HI 6007-01	—>	pH 7.010 Buffer Solution, 1 L
HI 6091	—>	pH 9.177 Buffer Solution, 500 mL
HI 6091-01	—>	pH 9.177 Buffer Solution, 1 L
HI 6010	—>	pH 10.010 Buffer Solution, 500 mL
HI 6010-01	—>	pH 10.010 Buffer Solution, 1 L
HI 6124	—>	pH 12.450 Buffer Solution, 500 mL
HI 6124-01	—>	pH 12.450 Buffer Solution, 1 L

A5.1.5 Electrode Cleaning Solutions

HI 7061M	—>	General Purpose Solution, 230 mL
HI 7061L	—>	General Purpose Solution, 500 mL
HI 7073M	—>	Protein Cleaning Solution, 230 mL
HI 7073L	—>	Protein Cleaning Solution, 500 mL
HI 7074M	—>	Inorganic Cleaning Solution, 230 mL
HI 7074L	—>	Inorganic Cleaning Solution, 500 mL
HI 7077M	—>	Oil & Fat Cleaning Solution, 230 mL
HI 7077L	—>	Oil & Fat Cleaning Solution, 500 mL

A5.1.6 Electrode Cleaning Solutions in FDA Approved Bottle

HI 8061M	—>	General Purpose Solution, 230 mL
HI 8061L	—>	General Purpose Solution, 500 mL
HI 8073M	—>	Protein Cleaning Solution, 230 mL
HI 8073L	—>	Protein Cleaning Solution, 500 mL
HI 8077M	—>	Oil & Fat Cleaning Solution, 230 mL
HI 8077L	—>	Oil & Fat Cleaning Solution, 500 mL

A5.1.7 Electrode Storage Solutions

HI 70300M	—>	Storage Solution, 230 mL
HI 70300L	—>	Storage Solution, 500 mL

A5.1.8 Electrode Storage Solutions in FDA Approved Bottle

HI 80300M	—>	Storage Solution, 230 mL
HI 80300L	—>	Storage Solution, 500 mL

A5.1.9 Refilling Electrolyte Solutions

HI 7071	—>	3.5M KCl + AgCl Electrolyte, 30 mL, for single junction electrodes
HI 7072	—>	1M KNO ₃ Electrolyte, 30 mL
HI 7075	—>	KNO ₃ and KCl Electrolyte, 30 mL

- HI 7076 —> 1M NaCl Electrolyte, 30 mL
HI 7078 —> $(\text{NH}_4)_2\text{SO}_4$ Electrolyte, 30 mL
HI 7082 —> 3.5M KCl Electrolyte, 30 mL, for double junction electrodes

A5.1.10 Refilling Electrolyte Solutions in FDA Approved Bottle

- HI 8071 —> 3.5M KCl + AgCl Electrolyte, 30 mL, for single junction electrodes
HI 8072 —> 1M KNO_3 Electrolyte, 30 mL
HI 8082 —> 3.5M KCl Electrolyte, 30 mL, for double junction electrodes

A5.1.11 ORP Pretreatment Solutions

- HI 7091M —> Reducing Pretreatment Solution, 230 mL
HI 7091L —> Reducing Pretreatment Solution, 500 mL
HI 7092M —> Oxidizing Pretreatment Solution, 230 mL
HI 7092L —> Oxidizing Pretreatment Solution, 500 mL

A5.1.12 Titration Reagents

- HI 70429 —> 0.05 M AgNO_3 Titration Reagent, 1 L
HI 70433 —> 0.01 N Stabilized Iodine Titration Reagent, 1 L
HI 70439 —> 0.1 M $\text{Na}_2\text{S}_2\text{O}_3$ Titration Reagent, 1 L
HI 70440 —> 0.02 N Stabilized Iodine Titration Reagent, 1 L
HI 70441 —> 0.04 N Stabilized Iodine Titration Reagent, 1 L
HI 70448 —> 0.02 M AgNO_3 Titration Reagent, 1 L
HI 70449 —> 0.02 M EDTA Titration Reagent, 1 L
HI 70455 —> 0.01 N NaOH Titration Reagent, 1 L
HI 70456 —> 0.1 N NaOH Titration Reagent, 1 L
HI 70457 —> 1 N NaOH Titration Reagent, 1 L
HI 70458 —> 0.01 M H_2SO_4 Titration Reagent, 1 L
HI 70459 —> 0.05 M H_2SO_4 Titration Reagent, 1 L
HI 70462 —> 0.01 N HCl Titration Reagent, 1 L
HI 70463 —> 0.1 N HCl Titration Reagent, 1 L
HI 70464 —> 1 N HCl Titration Reagent, 1 L

A5.1.13 Ion Selective Electrode Calibration Solutions

- HI 4001-01 —> 0.1 M Ammonia Standard
HI 4001-02 —> 100 ppm Ammonia Standard (as N)
HI 4001-03 —> 1000 ppm Ammonia Standard (as N)
HI 4002-01 —> 0.1 M Bromide Standard
HI 4003-01 —> 0.1 M Cadmium Standard
HI 4004-01 —> 0.1 M Calcium Standard
HI 4005-01 —> 0.1 M Carbon Dioxide Standard
HI 4005-03 —> 1000 ppm Carbon Dioxide Standard (as CaCO_3)
HI 4007-01 —> 0.1 M Chloride Standard

APPENDIX 5

HI 4007-02	—>	100 ppm Chloride Standard
HI 4007-03	—>	1000 ppm Chloride Standard
HI 4008-01	—>	0.1 M Cupric Standard
HI 4010-01	—>	0.1 M Fluoride Standard
HI 4010-02	—>	100 ppm Fluoride Standard
HI 4010-03	—>	1000 ppm Fluoride Standard
HI 4011-01	—>	0.1 M Iodide Standard
HI 4012-01	—>	0.1 M Lead Standard
HI 4012-21	—>	0.1 M Sulfate Standard
HI 4013-01	—>	0.1 M Nitrate Standard
HI 4013-02	—>	100 ppm Nitrate Standard
HI 4013-03	—>	1000 ppm Nitrate Standard
HI 4014-01	—>	0.1 M Potassium Standard
HI 4015-01	—>	0.1 M Silver Standard

A5.2 Sensors

A5.2.1 pH Electrodes

HI 1043B

Glass-body, double junction, refillable, combination pH electrode.

Use: strong acid/alkali

HI 1053B

Glass-body, triple ceramic, conic shape, refillable, combination pH electrode.

Use: emulsions

HI 1083B

Glass-body, micro, Viscolene, nonrefillable, combination pH electrode.

Use: biotechnology, micro titration

HI 1131B

Glass-body, single junction, refillable, combination pH electrode.

Use: general purpose

HI 1330B

Glass-body, semimicro, single junction, refillable, combination pH electrode.

Use: laboratory

HI 1331B

Glass-body, semimicro, single junction, refillable, combination pH electrode.

Use: flasks

HI 1230B

Plastic-body (PEI), double junction, gel-filled, combination pH electrode.

Use: general purpose

HI 2031B

Glass-body, semimicro, conic, refillable, combination pH electrode.

Use: semisolid products

HI 1332B

Plastic-body (PEI), double junction, refillable, combination pH electrode.

Use: general purpose.

FC 100B

Plastic-body (PVDF), double junction, refillable, combination pH electrode.

Use: cheese

FC 200B

Plastic-body (PVDF), single junction, conic, Viscolene, non-refillable, combination pH electrode.

Use: milk, yogurt, dairy products, and semi-solid foods

FC 210B

Glass-body, double junction, conic, Viscolene, combination pH electrode.

Use: milk, yogurt, and cream

FC 220B

Glass-body, single junction, refillable, combination pH electrode.

Use: creams, fruit juices, and sauces

FC 911B

Plastic-body (PVDF), double junction, refillable, combination pH electrode.

Use: creams, fruit juices, and sauces

HI 1413B

Glass-body, single junction, flat tip, Viscolene, combination pH electrode.

Use: surfaces, skin, leather, paper, emulsions

A5.2.2 ORP Electrodes**HI 3131B**

Glass-body, refillable, combination platinum ORP electrode.

Use: titration

HI 3230B

Plastic-body (PEI), gel-filled, combination platinum ORP electrode.

Use: general purpose

HI 4430B

Plastic-body (PEI), gel-filled, combination gold ORP electrode.

Use: general purpose

APPENDIX 5

A5.2.3 Half-cell Electrodes

HI 2110B

Glass-body, single half-cell pH electrode.

Use: general purpose

HI 5311

Glass-body, Ag/AgCl reference half-cell electrode, double junction, refillable with 4mm banana plug with 1m (3.3') cable.

Use: general purpose with wide temperature range

HI 5315

Plastic-body (PEI), double junction, Ag/AgCl reference half-cell electrode, refillable with 4mm plug with 1 m (3.3') cable.

Use: Ion Selective Electrodes

HI 5412

Glass-body, single Calomel reference half-cell electrode, refillable with 4mm plug with 1m (3.3') cable.

Use: general purpose with constant temperature range.

A5.2.4 Ion Selective Electrodes

HI 4101 Ammonia ISE

HI 4002 / HI 4102 Bromide ISE

HI 4003 / HI 4103 Cadmium ISE

HI 4004 / HI 4104 Chloride ISE

HI 4105 Carbon Dioxide ISE

HI 4007 / HI 4107 Chloride ISE

HI 4008 / HI 4108 Cupric ISE

HI 4009 / HI 4109 Cyanide ISE

HI 4010 / HI 4110 Fluoride ISE

HI 4011 / HI 4111 Iodide ISE

HI 4012 / HI 4112 Lead ISE

HI 4013 / HI 4113 Nitrate ISE

HI 4014 / HI 4114 Potassium ISE

HI 4015 / HI 4115 Silver / Sulfide ISE















FC 300B Sodium

A5.2.5 Temperature Sensor


HI 7662-T

Temperature probe with 1 m (3.3') paneled cable.


A5.3 Titrator components

	Pump assembly		Propeller
	Burette (25 mL syringe)		Stirrer support with positioning collar and positioning screw
	Aspiration tube with fitting and protection tube		Stirrer stand
	Dispensing tube with normal dispensing tip, fitting, protection tube and tube guide		Burette blank support
	Tube locks		Pump Assembly and Burette Blank support locking screw with plastic head
	Tool for burette cap removal		Temperature Probe
	Light spectrum protection screen		Power Cable
	Overhead Stirrer		USB Cable


APPENDIX 5



USB Storage Device



Shorting cap



Instruction Manual Binder


MAN902C
10/10

General Titration Applications Brochure

HI 902 Color

AUTOMATIC POTENTIOMETRIC TITRATOR

Revision 2.3



0.1N Sodium Hydroxide Titrant Concentration

Description:

Method for the standardization (titer determination) of 0.1N Sodium Hydroxide (NaOH) titrant solution against Potassium Hydrogen Phthalate (KHP). The results are expressed in **N (eq/L)**.

Reference:

AOAC Official Methods of Analysis, Official Method 936.16

Electrode:

- HI 1131B Combination pH Electrode
- HI 7662-T Temperature Probe

Reagents:

- HI 70456 0.1N Sodium Hydroxide solution (1 L)
- HI 70401 Potassium Hydrogen Phthalate (20 g)
- HI 70436 Distilled Water (1 gal)

Accessories:

- HI 70300L Storage Solution (500 mL)
- HI 7071 Electrode Fill Solution (4*30 mL)
- HI 7004L pH 4.01 Buffer Solution (500 mL)
- HI 7007L pH 7.01 Buffer Solution (500 mL)
- HI 7010L pH 10.01 Buffer Solution (500 mL)
- HI 740036P 100-mL Plastic Beaker (10 pcs)
- Analytical Balance with a minimum resolution of 0.0001g is recommended

NOTE: A pH electrode calibration is not needed for this method.

Procedure:

- Connect the pH electrode and temperature probe to the titrator.
- Press "*Select Method*" from the main screen. Use the arrow keys to highlight 'HI0001EN 0.1N NaOH Titr. Conc.' and press "*Select*".
- Install a 25-mL burette filled with 0.1N sodium hydroxide solution (HI 70456) on pump-one and verify that no air bubbles are present in the burette or tubing. If necessary prime until all air has been removed completely.
- Crush approximately 3 grams of potassium hydrogen phthalate (HI 70401) and dry it for 2 hours at 120°C. Cool to room temperature in a desiccator.
- Place a clean 100-mL plastic beaker on the analytical balance.
- Zero the balance.
- Carefully weigh approximately 0.20 grams of dried potassium hydrogen phthalate into the beaker.

NOTE: Ensure that all of the potassium hydrogen phthalate is on the bottom of the beaker.
- Record the exact weight of the sample once the balance has stabilized with an accuracy of 0.0001 grams.

- Remove the beaker from the balance and add distilled water to the 50-mL mark on the beaker.
- Place the beaker under the stirrer assembly and immerse the pH electrode, temperature probe and stirrer. Ensure that the reference junction of the pH probe is 5-6 mm below the surface. If necessary add extra distilled water.
- Press "*Start*". You will be prompted to enter the analyte size (weight of KHP). Use the numeric keypad to enter the exact weight and press "*Enter*" to start the analysis.

NOTE: Ensure that the KHP dissolves completely during the pre-titration stir time. Erroneous results may occur if the sample does not dissolve completely prior to titration. If necessary the pre-titration stir time can be increased.
- At the end of the titration, after detection of the equivalence point, 'titration completed' will appear with the titrant concentration. The result is expressed in **N (eq/L) of sodium hydroxide**.
- Remove the electrodes and stirrer from the sample and rinse them thoroughly with distilled water.
- Record the result (titer concentration).

NOTE: For improved accuracy, repeat this procedure a minimum of three times and calculate the average value.

For methods utilizing 0.1N NaOH titrant solution, follow the steps below to enter the titer/standardized value.

- Select the method utilizing 0.1N NaOH.
- Select 'Method Options' from the main screen.
- Using the arrow keys, highlight 'Titrant Conc.' and press "Select".
- Use the numeric keypad to enter the standardized (titer) value of the titrant then press 'Accept'.
- Press 'Escape' to exit the Method Options screen and select 'Save Method' option.

Method Parameters:

Name:	0.1N NaOH Titr. Conc.
Method Revision:	2.3
Titration Type:	Standard Titration
Analog Board:	Analog 1
Stirrer Configuration:	Stirrer 1
Pump Configuration:	
Titrant Pump:	Pump 1
Dosing Type:	Dynamic
min Vol:	0.030 mL
max Vol:	0.500 mL
delta E:	4.500 mV
End Point Mode:	pH 1EQ point, 1st Der.
Recognition Options:	
Threshold:	500 mV/mL
Range:	No
Filtered Derivatives:	No
Pre-Titration Volume:	5.000 mL

0.1N Sodium Hydroxide Titrant Concentration

Pre-Titration Stir Time: 60 sec
 Measurement Mode: Signal Stability
 delta E: 0.3 mV
 delta t: 5 sec
 t-min wait: 3 sec
 t-max wait: 30 sec
 Electrode Type: pH
 Blank Option: No Blank
 Calculations: Stdz. Titrant by Weight
 Dilution Option: Disabled
 Titrant Name: 0.1N NaOH
 Analyte Size: 0.200 g
 Analyte Entry: Manual
 Maximum Titrant Volume: 15.000 mL
 Stirring Speed: 1400 rpm
 Potential Range: -2000.0 to 2000.0 mV
 Volume/Flow Rate: 25 mL/50 mL/min
 Signal Averaging: 1 Reading
 Significant Figures: XXXXX

Calculations:

Calculations: Stdz. Titrant by Weight
 Titrant units: N (eq/L)
 Titrant volume dosed: V (L)
 Standard weight: 0.200 g
 (titrant/standard): 1.000 eq/mol
 MW of standard: 204.23 g/mol

$$\frac{\text{eq}}{\text{L}} \text{NaOH} = \frac{0.200 * 1.000}{204.23 * v(\text{L})}$$

Results:

Titration Report

Method Name: 0.1N NaOH Titr. Conc.
 Time & Date: 11:38 April 21, 2010
 Titration ID: Ti_00010

Titration Results

Method Name: 0.1N NaOH Titr. Conc.
 Time & Date: 11:38 April 21, 2010
 Analyte size: 0.2182 g
 End Point Volume: 10.056 mL
 pH Equivalence Point: 8.273
 Results: 0.10625 N (eq/L)
 Initial and Final pH: 4.060 to 8.757
 Titration Duration: 3:21 [mm:ss]
 Titration went to Completion
 Operator name: _____

0.1N Hydrochloric Acid Titrant Concentration

Description:

Method for the standardization (titer determination) of 0.1N Hydrochloric Acid (HCl) titrant solution against standardized 0.1N Sodium Hydroxide (NaOH) solution. The results are expressed in **N (eq/L)**.

Reference:

AOAC Official Methods of Analysis, Official Method 936.15

Electrode:

- HI 1131B Combination pH Electrode
- HI 7662-T Temperature Probe

Reagents:

- HI 70463 0.1N Hydrochloric Acid solution (1 L)
- HI 70456 0.1N Sodium Hydroxide solution (1 L)
- HI 70436 Distilled Water (1 gal)

Accessories:

- HI 70300L Storage Solution (500 mL)
- HI 7071 Electrode Fill Solution (4*30 mL)
- HI 7004L pH 4.01 Buffer Solution (500 mL)
- HI 7007L pH 7.01 Buffer Solution (500 mL)
- HI 7010L pH 10.01 Buffer Solution (500 mL)
- HI 740036P 100-mL Plastic Beakers (10 pcs)
- 10-mL Class-A Volumetric Pipette

NOTE: A pH electrode calibration is not needed for this method.

Procedure:

- Connect the pH electrode and temperature probe to the titrator.
- Press "*Select Method*" from the main screen. Use the arrow keys to highlight 'HI0002EN 0.1N HCl Titr. Conc.' and press "*Select*".
- Install a 25-mL burette filled with 0.1N hydrochloric acid solution (HI 70463) on pump-one and verify that no air bubbles are present in the burette or tubing. If necessary prime until all air has been removed completely.
- Use a class-A volumetric pipette to transfer exactly 10.00 mL of standardized 0.1N sodium hydroxide solution (HI 70456) to a clean 100-mL plastic beaker. Add distilled water to the 50-mL mark on the beaker.
For the determination of the exact concentration of sodium hydroxide, follow **HI0001EN 0.1N Sodium Hydroxide Titrant Concentration**.
- Place the beaker under the stirrer assembly and immerse the pH electrode, temperature probe and stirrer. Ensure that the reference junction of the pH probe is 5-6 mm below the surface. If necessary add extra distilled water.
- Press "*Start*". The titrator will start the analysis.

- At the end of the titration, after detection of the equivalence point, 'titration completed' will appear with the titrant concentration. The result is expressed in **N (eq/L) of hydrochloric acid**.
- Remove the electrodes and stirrer from the sample and rinse them thoroughly with distilled water.
- Record the result (titer concentration).

NOTE: For improved accuracy, repeat this procedure a minimum of three times and calculate the average value.

For methods utilizing 0.1N HCl titrant solution, follow the steps below to enter the titer/standardized value.

- Select the method utilizing 0.1N HCl.
- Select 'Method Options' from the main screen.
- Using the arrow keys, highlight 'Titrant Conc.' and press "Select".
- Use the numeric keypad to enter the standardized (titer) value of the titrant then press 'Accept'.
- Press 'Escape' to exit the Method Options screen and select 'Save Method' option.

Method Parameters:

Name:	0.1N HCl Titr. Conc.
Method Revision:	2.3
Titration Type:	Standard Titration
Analog Board:	Analog 1
Stirrer Configuration:	Stirrer 1
Pump Configuration:	
Titrant Pump:	Pump 1
Dosing Type:	Dynamic
min Vol:	0.030 mL
max Vol:	0.500 mL
delta E:	6.000 mV
End Point Mode:	pH 1EQ point, 1st Der.
Recognition Options:	
Threshold:	500 mV/mL
Range:	No
Filtered Derivatives:	No
Pre-Titration Volume:	5.000 mL
Pre-Titration Stir Time:	0 sec
Measurement Mode:	Signal Stability
delta E:	1.0 mV
delta t: 2.0 sec	
t-min wait: 3 sec	
t-max wait: 15 sec	
Electrode Type:	pH
Blank Option:	No Blank
Calculations:	Stdz. Titrant by Volume
Dilution Option:	Disabled
Titrant Name:	0.1N HCl
Analyte Size:	10.000 mL
Analyte Entry:	Fixed
Maximum Titrant Volume:	15.000 mL
Stirring Speed:	1400 rpm
Potential Range:	-2000.0 to 2000.0 mV
Volume/Flow Rate:	25 mL/50 mL/min
Signal Averaging:	1 Reading

0.1N Hydrochloric Acid Titrant Concentration

Significant Figures: XXXXX

Calculations:

Calculations: Stdz. Titrant by Volume
Titrant units: N (eq/L)
Titrant volume dosed: V (L)
Standard volume: 10.000 mL
Standard conc. 0.100 eq/L

$$\frac{\text{eq}}{\text{L}} \text{HCl} = \frac{10.000 * 0.100}{V(\text{L}) * 1000}$$

Results:

Titration Report

Method Name: 0.1N HCl Titr. Conc.
Time & Date: 14:28 April 21, 2010
Titration ID: Ti_00011

Titration Results

Method Name: 0.1N HCl Titr. Conc.
Time & Date: 14:28 April 21, 2010
Analyte size: 10.000 mL
End Point Volume: 9.979 mL
pH Equivalence Point: 5.059
Results: 0.10215 N (eq/L)
Initial and Final pH: 12.135 to 4.989
Titration Duration: 2:45 [mm:ss]
Titration went to Completion
Operator name: _____

0.1M Sodium Thiosulfate Titrant Concentration

Description:

Method for the standardization (titer determination) of 0.1M Sodium Thiosulfate ($\text{Na}_2\text{S}_2\text{O}_3$) titrant solution, against Potassium Iodate (KIO_3). The results are expressed in **M (mol/L)**.

Reference:

Adaptation of AOAC Official Methods of Analysis, Official Method 942.27

Electrode:

- HI 3131B Combination ORP Electrode

Reagents:

- HI 70439 0.1M Sodium thiosulfate solution (1 L)
- HI 70407 Potassium Iodate (20 g)
- HI 70425 16% Sulfuric Acid (500 mL)
- HI 70468 Potassium Iodide (35 g)
- HI 70436 Distilled Water (1 gal)

Accessories:

- HI 70300L Storage Solution (500 mL)
- HI 7071 Electrode Fill Solution (4*30 mL)
- Analytical Balance with a minimum resolution of 0.0001g is recommended
- 150-mL Glass Beakers
- 100-mL Class-A Volumetric Flask
- 10-mL Class-A Volumetric Pipette

Procedure:

- Connect the ORP probe to the titrator.
- Press “*Select Method*” from the main screen. Use the arrow keys to highlight ‘HI0003EN 0.1M $\text{Na}_2\text{S}_2\text{O}_3$ Titr. Conc.’ and press “*Select*”.
- Install a 25-mL burette with 0.1M sodium thiosulfate solution (HI 70439) on pump-one and verify that no air bubbles are present in the burette or tubing. If necessary prime until all air has been removed completely.
- Crush approximately 2 grams of potassium iodate (HI 70407) and dry it for 2 hours at 120°C. Cool to room temperature in a desiccator.
- Weigh 0.35 g of dried potassium iodate with an accuracy of 0.0001 g. Transfer the salt to a 100-mL volumetric flask. Add approximately 80 mL of distilled water and mix. Dissolve completely before bringing to volume.
- Use a class-A volumetric pipette to transfer exactly 10.00 mL of prepared standard solution to a 150-mL glass beaker and add distilled water to the 100-mL mark on the beaker.
- To the beaker add 5 mL of 16% sulfuric acid (HI 70425) and 1.5 g of potassium iodide (HI 70468).
- Place the beaker under the stirrer assembly and immerse the ORP electrode and stirrer. Ensure that the reference junction of the ORP electrode is

5-6 mm below the surface. If necessary add extra distilled water.

- Press “*Start*”. You will be prompted to enter the analyte size (weight of KIO_3). Use the numeric keypad to enter the exact weight and press “*Enter*” to start the analysis.
- At the end of the titration, after detection of the equivalence point, ‘titration completed’ will appear with the titrant concentration. The result is expressed in **M (mol/L) of sodium thiosulfate**.
- Remove the electrode and stirrer from the sample and rinse them thoroughly with distilled water.
- Record the result (titer concentration).

NOTE: For optimal accuracy repeat this procedure a minimum of three times and calculate the mean value.

For methods utilizing 0.1N $\text{Na}_2\text{S}_2\text{O}_3$ titrant solution, follow the steps below to enter the titer/standardized value.

- Select the method utilizing 0.1N $\text{Na}_2\text{S}_2\text{O}_3$.
- Select ‘Method Options’ from the main screen.
- Using the arrow keys, highlight ‘Titrant Conc.’ and press “*Select*”.
- Use the numeric keypad to enter the standardized (titer) value of the titrant then press ‘*Accept*’.
- Press ‘*Escape*’ to exit the Method Options screen and select ‘*Save Method*’ option.

Method Parameters:

Name:	0.1M $\text{Na}_2\text{S}_2\text{O}_3$ Titr. Conc.
Method Revision:	2.3
Titration Type:	Standard Titration
Analog Board:	Analog 1
Stirrer Configuration:	Stirrer 1
Pump Configuration:	
Titrant Pump:	Pump 1
Dosing Type:	Dynamic
min Vol:	0.030 mL
max Vol:	0.600 mL
delta E:	6.500 mV
End Point Mode:	mV 1EQ point, 1st Der.
Recognition Options:	
Threshold:	50 mV/mL
Range:	No
Filtered Derivatives:	No
Pre-Titration Volume:	5.000 mL
Pre-Titration Stir Time:	0 sec
Measurement Mode:	Signal Stability
delta E:	0.3 mV
delta t:	2.0 sec
t-min wait:	2 sec
t-max wait:	20 sec
Electrode Type:	ORP
Blank Option:	No Blank
Calculations:	Stdz. Titrant by Weight
Dilution Option:	Enabled
Final Dilution Volume:	100.000 mL
Aliquot Volume:	10.000 mL

0.1M Sodium Thiosulfate Titrant Concentration

Analyte size to be diluted: 0.350 g
 Titrant Name: 0.1M Na₂S₂O₃
 Analyte Size: 0.350 g
 Analyte Entry: Manual
 Maximum Titrant Volume: 15.000 mL
 Stirring Speed: 1400 rpm
 Potential Range: -2000.0 to 2000.0 mV
 Volume/Flow Rate: 25 mL/50 mL/min
 Signal Averaging: 1 Reading
 Significant Figures: XXXXX

Calculations:

Calculations: Stdz. Titrant by Weight
 Titrant units: M (mol/L)
 Titrant volume dosed: V (L)
 Standard weight: 0.350 g
 Dilution factor: 0.10
 Final dilution volume: 100.000 mL
 Aliquot volume: 10.000 mL
 (titrant/standard): 6.000 mol/mol
 MW of standard: 214.00 g/mol

$$\frac{\text{mol}}{\text{L}} \text{Na}_2\text{S}_2\text{O}_3 = \frac{0.350 * 0.10 * 6.0}{214.00 * V(\text{L})}$$

Results:

Titration Report

Method Name: 0.1M Na₂S₂O₃ Titr. Conc.
 Time & Date: 08:52 April 22, 2010
 Titration ID: Ti_00003

Titration Results

Method Name: 0.1M Na₂S₂O₃ Titr. Conc.
 Time & Date: 08:52 April 22, 2010
 Analyte size: 0.3523 g
 End Point Volume: 9.871 mL
 mV Equivalence Point: 220.4
 Results: 0.10007 M (mol/L)
 Initial and Final mV: 262.3 to 183.0
 Titration Duration: 3:58 [mm:ss]
 Titration went to Completion
 Operator name: _____

0.1M Ferrous Ammonium Sulfate Titrant Concentration

(Ferrous Ammonium Sulfate – FAS)

Description:

Method for the standardization (titer determination) of 0.1M Ferrous Ammonium Sulfate $[\text{Fe}(\text{NH}_4)_2(\text{SO}_4)_2 \cdot 6\text{H}_2\text{O}]$ titrant solution, against Potassium Dichromate $[\text{K}_2\text{Cr}_2\text{O}_7]$. The results are expressed in **M (mol/L)**.

Reference:

Standard Methods for the Examination of Water and Wastewater 21st edition, Method 5220B.

Electrode:

- HI 3131B Combination ORP Electrode

Reagents:

- HI 70444 25% Sulfuric Acid (500 mL)
- HI 70436 Distilled Water (1 gal)
- Ferrous Ammonium Sulfate Hexahydrate (ACS grade or better)
- Potassium Dichromate (ACS grade or better)

Accessories:

- HI 70300L Storage Solution (500 mL)
- HI 7071 Electrode Fill Solution (4*30 mL)
- HI 740036P100-mL Plastic Beakers (10 pcs.)
- Analytical Balance with a minimum resolution of 0.0001g is recommended
- 100-mL Class-A Volumetric Flasks
- 500-mL Class-A Volumetric Flasks
- 10-mL Class-A Volumetric pipette

Procedure:

- Connect the ORP electrode to the titrator.
- Press “*Select Method*” from the main screen. Use the arrow keys to highlight ‘HI0010EN 0.1M FAS Titr. Conc.’ and press “*Select*”.
- Install a 25-mL burette with 0.1M FAS solution on pump-one and verify that no air bubbles are present in the burette or tubing. If necessary prime until all air has been removed completely.

NOTE: To prepare a 0.1M FAS titrant solution: Add 19.607 g FAS to a 500-mL class-A volumetric flask. Add about 300 mL of distilled water, then add 40 mL of 25% sulfuric acid (HI 70444) to the flask. Invert the solution to mix. Allow the flask to return to room temperature before bringing to volume with distilled water.

- Dry approximately 2 grams of potassium dichromate for 2 hours at 150°C. Cool to room temperature in a desiccator before using.
- Weigh 0.49 g of dried potassium dichromate with an accuracy of 0.0001 g. Transfer the salt to a 100 mL volumetric flask. Add approximately 80 mL of distilled water and mix. Dissolve completely before bringing to volume.

- Use a class-A volumetric pipette to transfer exactly 10.00 mL of prepared standard solution to a clean 100-mL plastic beaker and 25 mL of 25% sulfuric acid solution (HI 70444). Use distilled water to bring the total volume to about 50 mL.
- Place the beaker under the stirrer assembly and immerse the ORP electrode and stirrer. Ensure that the reference junction of the ORP electrode is 5-6 mm below the surface. If necessary add extra distilled water.
- Press “*Start*”. You will be prompted to enter the analyte size (weight of $\text{K}_2\text{Cr}_2\text{O}_7$). Use the numeric keypad to enter the exact weight and press “*Enter*” to start the analysis.
- At the end of the titration, after detection of the equivalence point, the message ‘titration completed’ will appear with the titrant concentration. The result is expressed in **M (mol/L) of ferrous ammonium sulfate**.
- Remove the electrode and stirrer from the sample and rinse them thoroughly with distilled water.
- Record the result (titer concentration).

NOTE: For optimal accuracy repeat this procedure a minimum of three times and calculate the average value.

For methods utilizing 0.1M Ferrous Ammonium Sulfate titrant solution, follow the steps below to enter the titer/standardized value.

- Select the method utilizing 0.1M FAS.
- Select ‘Method Options’ from the main screen.
- Using the arrow keys, highlight ‘Titrant Conc.’ and press “*Select*”.
- Use the numeric keypad to enter the standardized (titer) value of the titrant then press ‘*Accept*’.
- Press ‘*Escape*’ to exit the Method Options screen and select ‘*Save Method*’ option.

NOTE: Standardize the Ferrous Ammonium Sulfate titrant solution daily.

Method Parameters:

Name:	0.1M FAS Titr. Conc.
Method Revision:	2.3
Titration Type:	Standard Titration
Analog Board:	Analog 1
Stirrer Configuration:	Stirrer 1
Pump Configuration:	
Titrant Pump:	Pump 1
Dosing Type:	Dynamic
min Vol:	0.030 mL
max Vol:	0.500 mL
delta E:	4.500 mV
End Point Mode:	mV 1EQ point, 1st Der.
Recognition Options:	
Threshold:	35 mV/mL

0.1M Ferrous Ammonium Sulfate Titrant Concentration (Ferrous Ammonium Sulfate – FAS)

Range: No
 Filtered Derivatives: No
 Pre-Titration Volume: 5.000 mL
 Pre-Titration Stir Time: 0 sec
 Measurement Mode: Signal Stability
 delta E: 0.5 mV
 delta t: 3.0 sec
 t-min wait: 2 sec
 t-max wait: 20 sec
 Electrode Type: ORP
 Blank Option: No Blank
 Calculations: Stdz. Titrant by Weight
 Dilution Option: Enabled
 Final Dilution Volume: 100.000 mL
 Aliquot Volume: 10.000 mL
 Analyte size to be diluted: 0.490 g
 Titrant Name: 0.1M FAS*6H2O
 Analyte Size: 0.490 g
 Analyte Entry: Manual
 Maximum Titrant Volume: 15.000 mL
 Stirring Speed: 1400 rpm
 Potential Range: -2000.0 to 2000.0 mV
 Volume/Flow Rate: 25 mL/50 mL/min
 Signal Averaging: 1 Reading
 Significant Figures: XXXXX

Calculations:

Calculations: Stdz. Titrant by Weight
 Titrant units: M (mol/L)
 Titrant volume dosed: V (L)
 Standard weight: 0.490 g
 Dilution factor: 0.10
 Final dilution volume: 100.000 mL
 Aliquot volume: 10.000 mL
 (titrant/standard): 6.000 mol/mol
 MW of standard: 294.18 g/mol

$$\frac{\text{mol}}{\text{L}} \text{ FAS} = \frac{0.490 * 0.10 * 6.0}{294.18 * V(\text{L})}$$

Results:

Titration Report

Method Name: 0.1M FAS Titr. Conc.
 Time & Date: 13:31 April 22, 2010
 Titration ID: Ti_00010

Titration Results

Method Name: 0.1M FAS Titr. Conc.
 Time & Date: 13:31 April 22, 2010
 Analyte size: 0.491 g
 End Point Volume: 9.879 mL
 mV Equivalence Point: 667.4
 Results: 0.1024 M (mol/L)
 Initial and Final mV: 791.3 to 598.0
 Titration Duration: 3:05 [mm:ss]
 Titration went to Completion
 Operator name: _____

0.02 M Silver Nitrate Titrant Concentration

Description:

Method for the standardization (titer determination) of 0.02M Silver Nitrate (AgNO_3) titrant solution against Sodium Chloride (NaCl). The results are expressed in **mol/L (M)**.

Reference:

AOAC Official Methods of Analysis, Official Method 941.18

Electrode:

- HI 4015 Silver Sulfide Half-cell ISE
- HI 5315 ISE Reference Electrode
- HI 7662-T Temperature Probe
- or -
- HI 4115 Silver Sulfide Combination ISE
- HI 7662-T Temperature Probe

Reagents:

- HI 70448 0.02N Silver Nitrate Solution (1 L)
- HI 70427 1.5M Nitric Acid Solution (500 mL)
- HI 70406 Sodium Chloride (20 g)
- HI 70436 Distilled Water (1 gal)

Other Accessories:

- HI 7072 Electrode Fill Solution (4*30 mL)
- 150-mL Glass Beakers
- 100-mL Class-A Volumetric Flask
- 5-mL Class-A Volumetric Pipette
- Analytical Balance with a minimum resolution of 0.0001g is recommended

Procedure:

- Connect the silver sulfide ISE, temperature probe, and reference electrode (if needed) to the titrator.
- Press "*Select Method*" from the main screen. Use the arrow keys to highlight 'HI0200EN 0.02M AgNO_3 Titr. Conc.' and press "*Select*".
- Install a 25-mL burette filled with 0.02M silver nitrate solution (HI 70448) on pump-one and verify that no air bubbles are present in the burette or tubing. If necessary prime until all air has been removed completely.
- Crush approximately 2 grams of sodium chloride (HI 70406) and dry it for 2 hours at 120°C. Cool to room temperature in a desiccator.
- Weigh 0.20 g of dried sodium chloride with an accuracy of 0.0001 g. Transfer the salt to a 100-mL volumetric flask. Add approximately 80 mL of distilled water and mix. Dissolve completely before bringing to volume.
- Use a class-A volumetric pipette to transfer exactly 5.00 mL of prepared standard solution to a 150-mL glass beaker and add distilled water to the 100-mL mark on the beaker.
- To the beaker add 10 mL of 1.5M nitric acid (HI 70427).

- Place the beaker under the stirrer assembly and immerse the electrodes and stirrer. Ensure that the reference junction of the ORP electrode is 5-6 mm below the surface. If necessary add extra distilled water.
- Press "*Start*". You will be prompted to enter the analyte size (weight of NaCl). Use the numeric keypad to enter the exact weight and press "*Enter*" to start the analysis.
- At the end of the titration, after detection of the equivalence point, 'titration completed' will appear with the titrant concentration. The result is expressed in **M (mol/L) of silver nitrate**.
- Remove the electrode and stirrer from the sample and rinse them thoroughly with distilled water.
- Record the result (titer concentration).

NOTE: For optimal accuracy repeat this procedure a minimum of three times and calculate the mean value.

For methods utilizing 0.02M AgNO_3 titrant solution, follow the steps below to enter the titer/standardized value.

- Select the method utilizing 0.02M AgNO_3 .
- Select 'Method Options' from the main screen.
- Using the arrow keys, highlight 'Titrant Conc.' and press "Select".
- Use the numeric keypad to enter the standardized (titer) value of the titrant then press 'Accept'.
- Press 'Escape' to exit the Method Options screen and select 'Save Method' option.

Method Parameters:

Name:	0.02M AgNO_3 Titr. Conc.
Method Revision:	2.3
Titration Type:	Standard Titration
Analog Board:	Analog 1
Stirrer Configuration:	Stirrer 1
Pump Configuration:	
Titrant Pump:	Pump 1
Dosing Type:	Dynamic
min Vol:	0.030 mL
max Vol:	0.500 mL
delta E:	5.000 mV
End Point Mode:	mV 1EQ point, 1st Der.
Recognition Options:	
Threshold:	100 mV/mL
Range:	No
Filtered Derivatives:	Yes
Pre-Titration Volume:	6.000 mL
Pre-Titration Stir Time:	0 sec
Measurement Mode:	Signal Stability
delta E:	1.0 mV
delta t:	1.5 sec
t-min wait:	2 sec
t-max wait:	20 sec
Electrode Type:	Silver Sulfide
Blank Option:	No Blank
Calculations:	Stdz. Titrant by Weight

0.02 M Silver Nitrate Titrant Concentration

Dilution Option: Enabled
 Final Dilution Volume: 100.000 mL
 Aliquot Volume: 5.000 mL
 Analyte size to be diluted: 0.200 g
 Titrant Name: 0.02M AgNO₃
 Analyte Size: 0.200 g
 Analyte Entry: Manual
 Maximum Titrant Volume: 15.000 mL
 Stirring Speed: 1400 rpm
 Potential Range: -2000.0 to 2000.0 mV
 Volume/Flow Rate: 25 mL/50 mL/min
 Signal Averaging: 1 Reading
 Significant Figures: XXXXX

Calculations:

Calculations: Stdz. Titrant by Weight
 Titrant units: M (mol/L)
 Titrant volume dosed: V (L)
 Standard weight: 0.200 g
 Dilution factor: 0.05
 Final dilution volume: 100.000 mL
 Aliquot volume: 5.000 mL
 (titrant/standard): 1.000 mol/mol
 MW of standard: 58.440 g/mol

$$\frac{\text{mol}}{\text{L}} \text{ AgNO}_3 = \frac{0.200 * 0.05 * 1.0}{58.440 * V(\text{L})}$$

Results:

Titration Report

Method Name: 0.02M AgNO₃ Titr. Conc.
 Time & Date: 08:52 June 30, 2010
 Titration ID: Ti_00003

Titration Results

Method Name: 0.02M AgNO₃ Titr. Conc.
 Time & Date: 08:52 June 30, 2010
 Analyte size: 0.2072 g
 End Point Volume: 8.8720 mL
 mV Equivalence Point: 269.3
 Results: 0.02033 M (mol/L)
 Initial and Final mV: 146.7 to 295.3
 Titration Duration: 2:11 [mm:ss]
 Titration went to Completion
 Operator name: _____

Alkalinity of Water

0-2500 mg/L CaCO₃, pH 4.5 Endpoint

Description:

Method for the determination of Total (Methyl Red) Alkalinity in water by titration of a sample to pH 4.5. The results are expressed in **mg/L (ppm) as Calcium Carbonate**.

For the determination of Phenolphthalein Alkalinity, set the endpoint to pH 8.3.

Reference:

Standard Methods for the Examination of Water and Wastewater 21st edition, Method 2320B.

Electrode:

- HI 1131B Combination pH Electrode
- HI 7662-T Temperature Probe

Reagents:

- HI 70463 0.1N Hydrochloric Acid solution (1 L)
- HI 70436 Distilled Water (1 gal)

Accessories:

- HI 70300L Storage Solution (500 mL)
- HI 7071 Electrode Fill Solution (4*30 mL)
- HI 7004 pH 4.01 buffer solution (500 mL)
- HI 7007 pH 7.01 buffer solution (500 mL)
- HI 7010 pH 10.01 buffer solution (500 mL)
- HI 740036P 100-mL Plastic Beakers (10 pcs)
- 50-mL Class-A Volumetric Pipette

Procedure:

- Connect the pH electrode and temperature probe to the titrator.
- Calibrate the electrode using pH 4.01, pH 7.01 and pH 10.01 buffer. Refer to the instruction manual for calibration procedure.
- Press "*Select Method*" from the main screen. Use the arrow keys to highlight 'HI1004EN Alkalinity of Water' and press "*Select*".
- Install a 25-mL burette with 0.1N hydrochloric acid solution (HI 70463) on pump-one and verify that no air bubbles are present in the burette or tubing. If necessary prime the burette until all air has been removed.
For the determination of the exact concentration of the hydrochloric acid follow, **HI0002EN 0.1N Hydrochloric Acid Titrant Concentration**.
- Use a class-A volumetric pipette to transfer exactly 50.00 mL of sample to a clean 100-mL plastic beaker.
- Place the beaker under the stirrer assembly and immerse the pH electrode, temperature probe and stirrer. Ensure that the reference junction of the pH probe is 5-6 mm below the surface.
- Press "*Start*". The titrator will begin the analysis.
- At the end of titration, when pH 4.50 is reached, 'titration completed' will appear with the alkalinity

concentration. The result is expressed in **mg/L of calcium carbonate**.

- Remove the electrodes and stirrer from the sample and rinse them thoroughly with distilled water.
- Record the result.

Method Parameters:

Name:	Alkalinity of Water
Method Revision:	2.3
Titration Type:	Standard Titration
Analog Board:	Analog 1
Stirrer Configuration:	Stirrer 1
Pump Configuration:	
Titrant Pump:	Pump 1
Dosing Type:	Dynamic
min Vol:	0.050 mL
max Vol:	0.500 mL
delta E:	5.000 mV
End Point Mode:	Fixed 4.500 pH
Pre-Titration Volume:	0.000 mL
Pre-Titration Stir Time:	0 sec
Measurement Mode:	Signal Stability
delta E:	1.0 mV
delta t:	2.0 sec
t-min wait:	2 sec
t-max wait:	20 sec
Electrode Type:	pH
Blank Option:	No Blank
Calculations:	Sample Calc. by Volume
Dilution Option:	Disabled
Titrant Name:	0.1 N HCl
Titrant Conc.:	0.1000 N (eq/L)
Analyte Size:	50.000 mL
Analyte Entry:	Fixed
Maximum Titrant Volume:	25.000 mL
Stirring Speed:	1400 rpm
Potential Range:	-2000.0 to 2000.0 mV
Volume/Flow Rate:	25 mL/50 mL/min
Signal Averaging:	1 Reading
Significant Figures:	XXXXX

Calculations:

Calculations:	Sample Calc. by Volume
Titrant units:	N (eq/L)
Titrant volume dosed:	V (L)
Final result units:	mg/L
Titrant conc.:	0.100 eq/L
(sample/titrant):	1.000 mol/eq
MW of sample:	100.09 g/mol
Sample volume:	50.000 mL

$$\frac{\text{mg}}{\text{L}} \text{CaCO}_3 = \frac{V(\text{L}) * 1000 * 0.10 * 0.5 * 100.09 * 1000}{50.00}$$

Alkalinity of Water
0-2500 mg/L CaCO₃, pH 4.5 Endpoint

Results:

Titration Report

Method Name: Alkalinity of Water
Time & Date: 09:04 April 18, 2010
Titration ID: Ti_00004

Titration Results

Method Name: Alkalinity of Water
Time & Date: 09:04 April 18, 2010
Analyte size: 50.000 mL
End Point Volume: 9.336 mL
pH Fixed End Point: 4.500
Results: 934.44 mg/L
Initial and Final pH: 10.232 to 4.419
Titration Duration: 3:23 [mm:ss]
Titration went to Completion
Operator name: _____

Acidity of Water

0-2500 mg/L CaCO₃, pH 8.3 Endpoint

Description:

Method for the determination of Total (Phenolphthalein) Acidity in water by titration of a sample to pH 8.3. The results are expressed in **mg/L (ppm) as Calcium Carbonate**.

For the determination of Methyl Orange Acidity, set the endpoint to pH 3.7.

Reference:

Standard Methods for the Examination of Water and Wastewater 21st edition, Method 2310B.

Electrode:

- HI 1131B Combination pH Electrode
- HI 7662-T Temperature Probe

Reagents:

- HI 70456 0.1N Sodium Hydroxide solution (1 L)
- HI 70436 Distilled Water (1 gal)

Accessories:

- HI 70300L Storage Solution (500 mL)
- HI 7071 Electrode Fill Solution (4*30 mL)
- HI 7004 pH 4.01 Buffer Solution (500 mL)
- HI 7007 pH 7.01 Buffer Solution (500 mL)
- HI 7010 pH 10.01 Buffer Solution (500 mL)
- HI 740036P100-mL Plastic Beakers (10 pcs)
- 50-mL Class-A Volumetric Pipette

Procedure:

- Connect the pH electrode and temperature probe to the titrator.
- Calibrate the electrode using pH 4.01, pH 7.01 and pH 10.01 buffer. Refer to the instruction manual for calibration procedure.
- Press "*Select Method*" from the main screen. Use the arrow keys to highlight 'HI1005EN Acidity of Water' and press "*Select*".
- Install a 25-mL burette with 0.1N sodium hydroxide acid solution (HI 70456) on pump-one and verify that no air bubbles are present in the burette or tubing. If necessary prime the burette until all air has been removed.
For the determination of the exact concentration of the sodium hydroxide follow, **HI0001EN 0.1N Sodium Hydroxide Titrant Concentration**.
- Use a class-A volumetric pipette to transfer exactly 50.00 mL of sample to a clean 100-mL plastic beaker.
- Place the beaker under the stirrer assembly and immerse the pH electrode, temperature probe and stirrer. Ensure that the reference junction of the pH probe is 5-6 mm below the surface.
- Press "*Start*". The titrator will begin the analysis
- At the end of titration, when pH 8.30 is reached, 'titration completed' will appear with the acidity

concentration. The result is expressed in **mg/L of calcium carbonate**.

- Remove the electrodes and stirrer from the sample and rinse them thoroughly with distilled water.
- Record the result.

Method Parameters:

Name:	Acidity of Water
Method Revision:	2.3
Titration Type:	Standard Titration
Analog Board:	Analog 1
Stirrer Configuration:	Stirrer 1
Pump Configuration:	
Titrant Pump:	Pump 1
Dosing Type:	Dynamic
min Vol:	0.050 mL
max Vol:	0.500 mL
delta E:	5.000 mV
End Point Mode:	Fixed 8.300 pH
Pre-Titration Volume:	0.000 mL
Pre-Titration Stir Time:	0 sec
Measurement Mode:	Signal Stability
delta E:	1.0 mV
delta t:	2.0 sec
t-min wait:	2 sec
t-max wait:	20 sec
Electrode Type:	pH
Blank Option:	No Blank
Calculations:	Sample Calc. by Volume
Dilution Option:	Disabled
Titrant Name:	0.1N NaOH
Titrant Conc.:	0.1000 N (eq/L)
Analyte Size:	50.000 mL
Analyte Entry:	Fixed
Maximum Titrant Volume:	25.000 mL
Stirring Speed:	1400 rpm
Potential Range:	-2000.0 to 2000.0 mV
Volume/Flow Rate:	25 mL/50 mL/min
Signal Averaging:	1 Reading
Significant Figures:	XXXXX

Calculations:

Calculations:	Sample Calc. by Volume
Titrant units:	N (eq/L)
Titrant volume dosed:	V (L)
Final result units:	mg/L
Titrant conc.:	0.100 eq/L
(sample/titrant):	0.500 mol/eq
MW of sample:	100.09 g/mol
Sample volume:	50.000 mL

$$\frac{\text{mg}}{\text{L}} \text{CaCO}_3 = \frac{V(\text{L}) * 1000 * 0.10 * 0.5 * 100.09 * 1000}{50.00}$$

Acidity of Water
0-2500 mg/L CaCO₃, pH 8.3 Endpoint

Results:

Titration Report

Method Name: Acidity of Water
Time & Date: 09:15 April 22, 2010
Titration ID: Ti_00005

Titration Results

Method Name: Acidity of Water
Time & Date: 09:15 April 22, 2010
Analyte size: 50.000 mL
End Point Volume: 5.879 mL
pH Fixed End Point: 8.300
Results: 588.43 mg/L
Initial and Final pH: 2.465 to 8.398
Titration Duration: 3:42 [mm:ss]
Titration went to Completion
Operator name: _____

Chloride in Water

0.00-150.00 mg/L

Description:

Method for the determination of chloride in water. The results are expressed in **mg/L (ppm) as Chloride**.

Reference:

Standard Methods for the Examination of Water and Wastewater 21st edition, Method 4500-Cl⁻.

Electrode:

- HI 4015 Silver Sulfide Half-cell ISE
- HI 5315 ISE Reference Electrode
- HI 7662-T Temperature Probe
- or -
- HI 4115 Silver Sulfide Combination ISE
- HI 7662-T Temperature Probe

Reagents:

- HI 70448 0.02N Silver Nitrate solution (1 L)
- HI 70427 1.5M Nitric Acid Solution (500 mL)
- HI 70436 Distilled Water (1 gal)

Accessories:

- HI 7072 Electrode Fill Solution (4*30 mL)
- 150-mL Glass Beakers
- 100-mL Class-A Volumetric Pipette
- 10-mL Class-A Volumetric Pipette

Procedure:

- Connect the silver sulfide ISE, reference electrode and temperature probe to the titrator.
- Press "*Select Method*" from the main screen. Use the arrow keys to highlight 'HI1007EN Chloride in Water' and press "*Select*".
- Install a 25-mL burette with 0.02M silver nitrate solution (HI 70448) on pump-one and verify that no air bubbles are present in the burette or tubing. If necessary prime until all air has been removed completely.
For the determination of the exact concentration of silver nitrate, follow **HI0200EN 0.02M Silver Nitrate Titrant Determination**.
- Use a class-A glass pipette to transfer exactly 100.00 mL of sample to a clean 150-mL glass beaker.
- Add 10 mL of 1.5 M nitric acid solution (HI 70427) to the beaker.
- Place the beaker under the stirrer assembly and immerse the silver sulfide ISE, reference electrode, temperature probe and stirrer. Ensure the electrodes are 5-6 mm below the surface.
- Press "*Start*". The titrator will begin the analysis.
- At the end of the titration, after detection of the equivalence point, 'titration completed' will appear with the chloride concentration. The result is expressed in **mg/L as chloride**.

- Remove the electrodes and stirrer from the sample and rinse thoroughly with distilled water.

Method Parameters:

Name:	Chloride in Water
Method Revision:	2.3
Titration Type:	Standard Titration
Analog Board:	Analog 1
Stirrer Configuration:	Stirrer 1
Pump Configuration:	
Titrant Pump:	Pump 1
Dosing Type:	Dynamic
min Vol:	0.030 mL
max Vol:	0.500 mL
delta E:	5.000 mV
End Point Mode:	mV 1EQ point, 1st Der.
Recognition Options:	
Threshold:	100 mV/mL
Range:	No
Filtered Derivatives:	No
Pre-Titration Volume:	0.000 mL
Pre-Titration Stir Time:	0 sec
Measurement Mode:	Signal Stability
delta E:	1.0 mV
delta t:	1.5 sec
t-min wait:	2 sec
t-max wait:	20 sec
Electrode Type:	Silver Sulfide ISE
Blank Option:	No Blank
Calculations:	Sample Calc. by Volume
Dilution Option:	Disabled
Titrant Name:	0.02M AgNO3
Titrant Conc.:	2.0000E-2 M (mol/L)
Analyte Size:	100.00 mL
Analyte Entry:	Manual
Maximum Titrant Volume:	25.000 mL
Stirring Speed:	1400 rpm
Potential Range:	-2000.0 to 2000.0 mV
Volume/Flow Rate:	25 mL/50 mL/min
Signal Averaging:	1 Reading
Significant Figures:	XXXXX

Calculations:

Calculations:	Sample Calc. by Volume
Titrant units:	M (mol/L)
Titrant volume dosed:	V (L)
Final result units:	ppm (mg/L)
Titrant conc.:	0.02000 mol/L
(sample/titrant):	1.000 mol/mol
MW of sample:	35.453 g/mol
Sample volume:	100.00 mL
$\frac{\text{mg}}{\text{L}}$	$\frac{V(\text{L}) * 1000 * 0.02 * 1.0 * 35.45 * 1000}{100.00}$

Chloride in Water

0.00-150.00 mg/L

Results:

Titration Report

Method Name: Chloride in Water
Time & Date: 09:20 April 23, 2010
Titration ID: Ti_00029

Titration Results

Method Name: Chloride in Water
Time & Date: 09:20 April 23, 2010
Analyte size: 100.00 mL
End Point Volume: 4.781 mL
mV Equivalence Point: 280.3
Results: 33.897 mg/L
Initial and Final mV: 194.8 to 298.5
Titration Duration: 1:24 [mm:ss]
Titration went to Completion
Operator name: _____

Neutralization with Sulfuric Acid

0.00-200.00 meq/L

Description:

Method for the determination of concentration of strong or weak bases by titration of a sample to the equivalence point. The results are expressed in meq/L.

Electrode:

- HI 1131B Combination pH Electrode
- HI 7662-T Temperature Probe

Reagents:

- HI 70459 0.05M Sulfuric Acid solution (1 L)
- HI 70436 Distilled Water (1 gal.)

Accessories:

- HI 70300L Storage Solution (500 mL)
- HI 7071 Electrode Fill Solution (4*30 mL)
- HI 7004L pH 4.01 Buffer Solution (500 mL)
- HI 7007L pH 7.01 Buffer Solution (500 mL)
- HI 7010L pH 10.01 Buffer Solution (500 mL)
- HI 740036P 100-mL Plastic Beakers (10 pcs.)
- 10-mL Class-A Volumetric Pipette

NOTE: For this method a pH electrode calibration is not necessary.

Procedure:

- Connect the pH electrode and temperature probe to the titrator.
- Press "*Select Method*" from the main screen. Use the arrow keys to highlight 'HI1008EN Neutralization w/ H2SO4' and press "*Select*".
- Install a 25-mL burette filled with 0.05N sulfuric acid solution (HI 70459) on pump-one and verify that no air bubbles are present in the burette or tubing. If necessary prime until all air has been removed completely.

For the determination of the exact concentration of sulfuric acid, follow **HI0103EN 0.05M Sulfuric Acid Titrant Determination**.

- Use a class-A volumetric pipette to transfer exactly 10.00 mL of sample to a clean 100-mL plastic beaker, and use distilled water to bring up the volume to about 50 mL.
- Place the beaker under the stirrer assembly and immerse the pH electrode, temperature probe and stirrer. Ensure that the reference junction of the pH probe is 5-6 mm below the surface. If necessary add extra distilled water.
- Press "*Start*". The titrator will start the analysis.
- At the end of titration, after detection of the equivalence point, 'titration completed' will appear together with the base concentration. The result is expressed in meq/L.

- Remove the electrodes and stirrer from the sample and rinse them thoroughly with distilled water.
- Record the result.

Method Parameters:

Name: Neutralization w/ H2SO4
 Method Revision: 2.3
 Titration Type: Standard Titration
 Analog Board: Analog 1
 Stirrer Configuration: Stirrer 1
 Pump Configuration:
 Titrant Pump: Pump 1
 Dosing Type: Dynamic
 min Vol: 0.050 mL
 max Vol: 0.500 mL
 delta E: 20.000 mV
 End Point Mode: pH 1EQ point, 1st Der.
 Recognition Options:
 Threshold: 50 mV/mL
 Range: No
 Filtered Derivatives: No
 Pre-Titration Volume: 0.000 mL
 Pre-Titration Stir Time: 0 sec
 Measurement Mode: Signal Stability
 delta E: 1.0 mV
 delta t: 2.0 sec
 t-min wait: 2 sec
 t-max wait: 15 sec
 Electrode Type: pH
 Blank Option: No Blank
 Calculations: Sample Calc. by Volume
 Dilution Option: Disabled
 Titrant Name: 0.05M H2SO4
 Titrant Conc.: 5.0000E-2 M (mol/L)
 Analyte Size: 10.000 mL
 Analyte Entry: Fixed
 Maximum Titrant Volume: 20.000 mL
 Stirring Speed: 1400 rpm
 Potential Range: -2000.0 to 2000.0 mV
 Volume/Flow Rate: 25 mL/50 mL/min
 Signal Averaging: 1 Reading
 Significant Figures: XXXXX

Calculations:

Calculations: Sample Calc. by Volume
 Titrant units: M (mol/L)
 Titrant volume dosed: V (L)
 Final result units: meq/L
 Titrant conc.: 5.000E-2 mol/L
 (sample/titrant): 2.000 eq/mol
 Sample volume: 10.000 mL

$$\frac{\text{meq}}{\text{L}} = \frac{V(\text{L}) * 1000 * 0.05 * 2.0 * 1000}{10.00}$$

Neutralization with Sulfuric Acid

0.00-200.00 meq/L

Results:

Titration Report

Method Name: Neutralization w/ H2SO4
Time & Date: 13:04 April 23, 2010
Titration ID: Ti_00008

Titration Results

Method Name: Neutralization w/ H2SO4
Time & Date: 13:04 April 23, 2010
Analyte size: 10.000 mL
End Point Volume: 9.562 mL
pH Equivalence Point: 7.966
Results: 95.620 meq/L
Initial and Final pH: 11.655 to 6.248
Titration Duration: 3:26 [mm:ss]
Titration went to Completion
Operator name: _____

Neutralization with Sodium Hydroxide

0.00-200.00 meq/L

Description:

Method for the determination of concentration of strong or weak acids, by titration of a sample to the first equivalence point. The results are expressed in **meq/L**.

Electrode:

- HI 1131B Combination pH Electrode
- HI 7662-T Temperature Probe

Reagents:

- HI 70456 0.1N Sodium Hydroxide solution (1 L)
- HI 70436 Distilled Water (1 gal)

Accessories:

- HI 70300L Storage Solution (500 mL)
- HI 7071 Electrode Fill Solution (4*30 mL)
- HI 7004L pH 4.01 Buffer Solution (500 mL)
- HI 7007L pH 7.01 Buffer Solution (500 mL)
- HI 7010L pH 10.01 Buffer Solution (500 mL)
- HI 740036P 100 mL Plastic Beakers (10 pcs)
- 10-mL Class-A Volumetric Pipette

NOTE: For this method a pH electrode calibration is not necessary.

Procedure:

- Connect the pH electrode and temperature probe to the titrator.
- Press "*Select Method*" from the main screen. Use the arrow keys to highlight 'HI1009EN Neutralization w/ NaOH' and press "*Select*".
- Install a 25-mL burette with 0.1N sodium hydroxide solution (HI 70456) on pump-one and verify that no air bubbles are present in the burette or tubing. If necessary prime until all air has been removed completely.

For the determination of the exact concentration of sodium hydroxide, follow **HI0001EN 0.1N Sodium Hydroxide Titrant Determination**.

- Use a class-A volumetric pipette to transfer exactly 10.00 mL of sample to a clean 100-mL plastic beaker, and use distilled water to bring up the volume to about 50 mL.
- Place the beaker under the stirrer assembly and immerse the pH electrode, temperature probe and stirrer. Ensure that the reference junction of the pH probe is 5-6 mm below the surface. If necessary add extra distilled water.
- Press "*Start*". The titrator will start the analysis.
- At the end of the titration, after detection of the equivalence point, 'titration completed' will appear with the acid concentration. The result is expressed in **meq/L**.

- Remove the electrodes and stirrer from the sample and rinse them thoroughly with distilled water.
- Record the result.

Method Parameters:

Name:	Neutralization w/ NaOH
Method Revision:	2.3
Titration Type:	Standard Titration
Analog Board:	Analog 1
Stirrer Configuration:	Stirrer 1
Pump Configuration:	
Titrant Pump:	Pump 1
Dosing Type:	Dynamic
min Vol:	0.050 mL
max Vol:	0.500 mL
delta E:	20.000 mV
End Point Mode:	pH 1EQ point, 1st Der.
Recognition Options:	
Threshold:	50 mV/mL
Range:	No
Filtered Derivatives:	No
Pre-Titration Volume:	0.000 mL
Pre-Titration Stir Time:	0 sec
Measurement Mode:	Signal Stability
delta E:	1.0 mV
delta t:	2.0 sec
t-min wait:	2 sec
t-max wait:	15 sec
Electrode Type:	pH
Blank Option:	No Blank
Calculations:	Sample Calc. by Volume
Dilution Option:	Disabled
Titrant Name:	0.1N NaOH
Titrant Conc.:	0.1000 N(eq/L)
Analyte Size:	10.000 mL
Analyte Entry:	Fixed
Maximum Titrant Volume:	20.000 mL
Stirring Speed:	1400 rpm
Potential Range:	-2000.0 to 2000.0 mV
Volume/Flow Rate:	25 mL/50 mL/min
Signal Averaging:	1 Reading
Significant Figures:	XXXXX

Calculations:

Calculations:	Sample Calc. by Volume
Titrant units:	N (eq/L)
Titrant volume dosed:	V (L)
Final result units:	meq/L
Titrant conc.:	0.1000 mol/L
(sample/titrant):	1.000 eq/mol
Sample volume:	10.000 mL

$$\frac{\text{meq}}{\text{L}} = \frac{V(\text{L}) * 1000 * 0.1 * 1.0 * 1000}{10.00}$$

Neutralization with Sodium Hydroxide

0.00-200.00 meq/L

Results:

Titration Report

Method Name: Neutralization w/ NaOH
Time & Date: 13:04 April 23, 2010
Titration ID: Ti_00009

Titration Results

Method Name: Neutralization w/ NaOH
Time & Date: 13:04 April 23, 2010
Analyte size: 10.000 mL
End Point Volume: 15.970 mL
pH Equivalence Point: 8.431
Results: 159.70 meq/L
Initial and Final pH: 2.675 to 10.316
Titration Duration: 3:20 [mm:ss]
Titration went to Completion
Operator name: _____

Troubleshooting 1

Description:

A method for verifying the dosing accuracy of the titrator. This method should be used to troubleshoot a titrator equipped with a 25-mL burette. The titrator dispenses a 20.000-mL pre-titration volume, waits 20 seconds, and dispenses an additional 20.000-mL dose, bringing the total volume to 40.000 mL. This procedure can also be used to check the stability of the temperature and mV channels.

Reference:

ISO/TC 48/SC1N 380E to 383E: "Piston and/or Plunger Operated Volumetric Apparatus".

Accessories:

- HI 762000C 0°C Temperature Key
- HI 762070C 70°C Temperature Key
- HI 70436 Distilled Water (1 gal)
- HI 7662-T Temperature Probe
- Shorting cap
- Analytical Balance with a minimum resolution of 0.0001 g is recommended
- Narrow Neck Beaker

Large Dose Dispensing Procedure:

- Connect the shorting cap to the BNC socket on Analog Board 1.
- Press "Select Method" from the main screen. Use the arrow keys to highlight 'HI1011EN Troubleshooting 1' and press "Select".
- Install the 25-mL burette with HI 70436 room temperature (25°C) distilled water on pump-one and verify that no air bubbles are present in the burette or tubing. If necessary prime until all air has been removed completely.
- Add a small amount of water to a narrow neck beaker to have vapor-saturated air space just above the liquid level and minimize evaporation.
- Place the narrow neck beaker on an analytical balance.
- Zero the balance.
- Place the dosing tip through the neck of the beaker. Take care not to immerse it in the liquid during dispensing and not to touch the beaker's walls.
- Press "Start" to start dosing.
- Write down the mass read from the balance after each dose.
- The following information is needed to verify the accuracy of the dosing system:
 - The temperature of the dispensed water
 - The atmospheric air pressure
 - The density of the weight used to calibrate the balance

Note: This procedure can be repeated on pump 2

Method Parameters:

Name:	Troubleshooting 1
Method Revision:	2.3
Titration Type:	Standard Titration
Analog Board:	Analog 1
Stirrer Configuration:	Stirrer 1
Pump Configuration:	
Titrant Pump:	Pump 1
Dosing Type:	Linear - 20.000 mL
End Point Mode:	Fixed 10.0 mV
Pre-Titration Volume:	20.000 mL
Pre-Titration Stir Time:	0 sec
Measurement Mode:	Timed Increment
t-incr. Wait:	20 sec
Electrode Type:	Shorting Cap
Blank Option:	No Blank
Calculations:	No Formula (mL only)
Titrant Name:	DI water
Maximum Titrant Volume:	40.000 mL
Stirring Speed:	0 rpm
Potential Range:	-2000.0 to 2000.0 mV
Volume/Flow Rate:	25 mL/50 mL/min
Signal Averaging:	1 Reading
Significant Figures:	XXXXX

Calculations:

The measured volume of the dispensed liquid is calculated from the measured mass using the following equation:

$$V = m * \frac{1}{\rho} * \left(1 + \frac{\rho_{air}}{\rho_L} - \frac{\rho_{air}}{\rho_{std}} \right)$$

V	Volume of measured mass of water [mL]
M	Measured mass of water [g]
ρ_L	Density of dispensed water [g/mL]
ρ_{air}	Density of ambient air [g/mL]
ρ_{STD}	Density of calibration standard weight [g/mL]

If the actual values of the above parameters are not accessible the following equation can be used:

$$V = M * F$$

V	Volume of measured mass of water [mL]
F	Transformation factor

The transformation factor takes into account the air buoyancy, the water density and their temperature dependence. Standard values can be used to obtain the transformation factor.

The values from the table below have been calculated by correcting the air and water density with temperature, assuming: dry air at 760 torr (dry air at 760 torr and 20°C has a density $\rho_{air} = 0.0012$ g/mL) and density of calibration steel-standard weights $\rho_{STD} = 8$ g/mL.

Troubleshooting 1

Temperature (°C)	FACTOR
17.0	1.002290
18.0	1.002467
19.0	1.002654
20.0	1.002853
21.0	1.003061
22.0	1.003282
23.0	1.003512
24.0	1.003752
25.0	1.004002
26.0	1.004261
27.0	1.004531
28.0	1.004809
29.0	1.005097
30.0	1.005395

The specifications of the dosing accuracy is $\pm 0.1\%$ of full burette volume ($\pm 0.1\% * 25 \text{ mL} = \pm 0.025 \text{ mL}$). For the accuracy of other burette volumes, see the HI 901 / HI 902 Instruction Manual.

If the results are not correct, check all fittings for leakage, and burette and tubing for air bubbles. Repeat the measurements.

The table below contains the settings to be used for other burette volumes:

Burette Volume	Dosing Type	Pre-Titr. Volume	Max. Titr. Volume
5 mL	4.000 mL	4.000 mL	8.000 mL
10 mL	8.000 mL	8.000 mL	16.000 mL
25 mL	20.000 mL	20.000 mL	40.000 mL

Temperature Channel Fast Check Procedure:

- Connect the shorting cap to the BNC socket on Analog Board 1.
- Connect the HI 762000C 0°C temperature key to the RCA socket (temperature sensor input) on Analog Board 1.
- Once on the main screen select “Mode”, ensure that Analog Board 1 is active then select “mV 1”.
- Press “General Options”, use the arrow keys to highlight ‘Temperature’ and press “Select”.
- Select ‘Temperature Source’ and then ‘Automatic Temperature’.
- Press “Escape” two times to return to the main screen.
- The titrator should display ATC $0.0 \pm 0.4^\circ\text{C}$ with no fluctuation or drift.
- Connect the HI 762070C 70°C temperature key to the RCA socket (temperature sensor input) on Analog Board 1
- The titrator should display ATC $70.0 \pm 0.4^\circ\text{C}$ with no fluctuation or drift.

Note: This procedure can be repeated on analog 2

Temperature & mV Channel Logging Procedure:

- Connect the shorting cap to the BNC socket on Analog Board 1.
- Connect the HI 762000C 0°C temperature key to the RCA socket (temperature sensor input) on Analog Board 1.
- Once on the main screen select “Mode”, ensure that Analog Board 1 is active then select “mV 1”.
- Press “General Options”, use the arrow keys to highlight ‘Temperature’ and press “Select”.
- Select ‘Temperature Source’ then ‘Automatic Temperature’.
- Press “Escape” two times to return to the main screen.
- Press “mV Setup”, use the arrow keys to highlight ‘Logging Interval’. Set the logging interval to 15 seconds and press “Accept”. Press “Escape” to return to the main screen.
- Press the “result” key and use the arrow keys to highlight ‘Setup pH/mV/ISE Report’. Press “Select”.
- Select ‘Potential’ and ‘Temperature and Units’ (the selected fields are marked with an “*”). All other fields should be unselected.
- Press “Save Report” to return to the Data Parameters screen.
- Press “Escape” to return to the main screen. Once on the main screen, select “Mode” then “mV 1” to enter mV mode on Analog Board 1.
- Press “Start Log” to start the automatic log.
- Let the log run for about 10 minutes. Press “Stop Log” to stop the automatic log.
- Press “Result”, highlight ‘Review Available Reports’, and press “View Report”.
- The mV column should display $0.0 \pm 0.1 \text{ mV}$, and the temperature column should display $0.0^\circ\text{C} \pm 0.4^\circ\text{C}$.
- This procedure can be repeated using the HI 762070C 70°C temperature key.

Note: This procedure can be repeated on analog 2

If any problems occur, please contact your nearest Hanna Service Office.

Troubleshooting 2

Description:

A method for verifying the dosing accuracy of the titrator. This method should be used to troubleshoot a titrator equipped with a 25-mL burette. The titrator dispenses a 10.000-mL pre-titration volume, waits 10 seconds, and then dispenses an additional 0.500-mL dose 20 times, waiting 10 seconds between each addition, bringing the total volume to 20.000 mL. This procedure can also be used to check the stirrer functionality.

Reference:

ISO/TC 48/SC1N 380E to 383E: "Piston and/or Plunger Operated Volumetric Apparatus".

Accessories:

- HI 762000C 0°C Temperature Key
- HI 70436 Distilled Water (1 gal)
- HI 7662-T Temperature Probe
- Shorting cap
- Analytical Balance with a minimum resolution of 0.0001 g is recommended
- Narrow Neck Beaker

Small Dose Dispensing Procedure:

- Connect the shorting cap to the BNC socket on Analog Board 1.
- Press "*Select Method*" from the main screen. Use the arrow keys to highlight 'HI1012EN Troubleshooting 2' and press "*Select*".
- Install the 25-mL burette with HI 70436 room temperature (25°C) distilled water on pump-one and verify that no air bubbles are present in the burette or tubing. If necessary prime until all air has been removed completely.
- Add a small amount of water to a narrow neck beaker to have vapor-saturated air space just above the liquid level and minimize evaporation.
- Place the narrow neck beaker on an analytical balance.
- Zero the balance.
- Place the dosing tip through the neck of the beaker. Take care not to immerse it in the liquid during dispensing and not to touch the beaker's walls.
- Press "*Start*" to start dosing.
- Write down the mass read from the balance after each dose.
- The following information is needed to verify the accuracy of the dosing system:
 - The temperature of the dispensed water
 - The atmospheric air pressure
 - The density of the weight used to calibrate the balance

Note: This procedure can be repeated on pump 2

Method Parameters:

Name:	Troubleshooting 2
Method Revision:	2.3
Titration Type:	Standard Titration
Analog Board:	Analog 1
Stirrer Configuration:	Stirrer 1
Pump Configuration:	
Titrant Pump:	Pump 1
Dosing Type:	Linear - 0.500 mL
End Point Mode:	Fixed 10.0 mV
Pre-Titration Volume:	10.000 mL
Pre-Titration Stir Time:	0 sec
Measurement Mode:	Timed Increment
t-incr. Wait:	10 sec
Electrode Type:	Shorting Cap
Blank Option:	No Blank
Calculations:	No Formula (mL only)
Titrant Name:	DI water
Maximum Titrant Volume:	20.000 mL
Stirring Speed:	0 rpm
Potential Range:	-2000.0 to 2000.0 mV
Volume/Flow Rate:	25 mL/50 mL/min
Signal Averaging:	1 Reading
Significant Figures:	XXXXX

Calculations:

The measured volume of the dispensed liquid is calculated from the measured mass using the following equation:

$$V = m * \frac{1}{\rho} * \left(1 + \frac{\rho_{air}}{\rho_L} - \frac{\rho_{air}}{\rho_{std}} \right)$$

V	Volume of measured mass of water [mL]
M	Measured mass of water [g]
ρ_L	Density of dispensed water [g/mL]
ρ_{air}	Density of ambient air [g/mL]
ρ_{STD}	Density of calibration standard weight [g/mL]

If the actual values of the above parameters are not accessible the following equation can be used:

$$V = M * F$$

V	Volume of measured mass of water [mL]
F	Transformation factor

The transformation factor takes into account the air buoyancy, the water density and their temperature dependence. Standard values can be used to obtain the transformation factor.

The values from the table below have been calculated by correcting the air and water density with temperature, assuming: dry air at 760 torr (dry air at 760 torr and 20°C has a density $\rho_{air} = 0.0012$ g/mL) and density of calibration steel-standard weights $\rho_{STD} = 8$ g/mL.

Troubleshooting 2

Temperature (°C)	FACTOR
17.0	1.002290
18.0	1.002467
19.0	1.002654
20.0	1.002853
21.0	1.003061
22.0	1.003282
23.0	1.003512
24.0	1.003752
25.0	1.004002
26.0	1.004261
27.0	1.004531
28.0	1.004809
29.0	1.005097
30.0	1.005395

The specifications of the dosing accuracy is $\pm 0.1\%$ of full burette volume ($\pm 0.1\% * 25 \text{ mL} = \pm 0.025 \text{ mL}$). For the accuracy of other burette volumes, see the HI 901 / HI 902 Instruction Manual.

If the results are not correct, check all fittings for leakage, and burette and tubing for air bubbles. Repeat the measurements.

The table below contains the settings to be used for other burette volumes:

Burette Volume	Dosing Type	Pre-Titr. Volume	Max. Titr. Volume
5 mL	4.000 mL	4.000 mL	8.000 mL
10 mL	8.000 mL	8.000 mL	16.000 mL
25 mL	20.000 mL	20.000 mL	40.000 mL

Stirring Device Fast Check Procedure:

- From the titrators main screen, press “stir” on the keyboard and use the arrow keys to set the stir speed at 100 RPM.
- Slowly increase the stir speed up to 2500 RPM.
- Check that the propeller stirrer starts turning faster and faster, following the commands.

Note: This procedure can be repeated on stirrer 2

If any problems occur, please contact your nearest Hanna Service Office.

Concentration of Phosphoric Acid

0.00000–1.0000⁻² mol/L

Description:

Method for the determination of phosphoric acid (H_3PO_4), by titration of a sample to the point of inflection.

The first inflection point corresponds to the H_3PO_4 content, and the difference between the first and second corresponds to H_2PO_4^- . The results are expressed in **M (mol/L)**.

If only phosphoric acid and no other acids or bases are present in the sample, then $[\text{H}_3\text{PO}_4] = [\text{H}_2\text{PO}_4^-]$. If $[\text{H}_3\text{PO}_4] > [\text{H}_2\text{PO}_4^-]$, this means that other strong acids are also present (e.g., hydrochloric or sulfuric acid). If $[\text{H}_3\text{PO}_4] < [\text{H}_2\text{PO}_4^-]$ this means that other weak acids or bases are present (e.g., citric acid / citrate, or ascorbic acid / ascorbate).

Electrode:

- HI 1131B Combination pH Electrode
- HI 7662-T Temperature Probe

Reagents:

- HI 70456 0.1N Sodium Hydroxide solution (1 L)
- HI 70436 Distilled Water (1 gal)

Accessories:

- HI 70300L Storage Solution (500 mL)
- HI 7071 Electrode Fill Solution (4*30 mL)
- HI 7004L pH 4.01 Buffer Solution (500 mL)
- HI 7007L pH 7.01 Buffer Solution (500 mL)
- HI 7010L pH 10.01 Buffer Solution (500 mL)
- 100-mL Class-A Volumetric Pipette
- 150-mL Glass Beakers

NOTE: For this method it is not necessary to calibrate the pH electrode first.

Procedure:

- Connect the pH electrode and temperature probe to the titrator.
- Press “*Select Method*” from the main screen. Use the arrow keys to highlight ‘HI1014EN Concentration of H3PO4’ and press “*Select*”.
- Install a 25-mL burette with 0.1N sodium hydroxide solution (HI 70456) on pump-one and verify that no air bubbles are present in the burette or tubing. If necessary prime until all air has been removed completely.
- For sparkling samples, remove the CO_2 by passing nitrogen through it, using an ultrasonic bath, or by stirring and applying vacuum.
- Use a class-A volumetric pipette to transfer exactly 100.00 mL of sample to a clean 150-mL glass beaker.
- Place the beaker under the stirrer assembly and immerse the pH electrode, temperature probe and

stirrer. Ensure that the reference junction of the pH probe is 5-6 mm below the surface.

- Press “*Start*”. The titrator will start the analysis.
- At the end of the titration, after detection of the second equivalence point, the message ‘titration completed’ will appear with the phosphoric acid concentration. The first equivalence point corresponds to the H_3PO_4 concentration and the second equivalence point corresponds to the H_2PO_4^- concentration.
- Remove the electrodes and stirrer from the sample and rinse them thoroughly with distilled water.
- Record the result.

Method Parameters:

Name:	Concentration of H3PO4
Method Revision:	2.3
Titration type:	Standard Titration
Analog Board:	Analog 1
Stirrer Configuration:	Stirrer 1
Pump Configuration:	
Titrant Pump:	Pump 1
Dosing Type:	Dynamic
min Vol:	0.030 mL
max Vol:	0.500 mL
delta E:	8.000 mV
End Point Mode:	pH 2EQ points, 1st Der.
Recognition Options:	
Threshold:	50 mV/mL
Range:	No
Filtered Derivatives:	No
Pre-Titration Volume:	0.000 mL
Pre-Titration Stir Time:	10 sec
Measurement Mode:	Signal Stability
delta E:	0.8 mV
delta t:	1.5 sec
t-min wait:	2 sec
t-max wait:	20 sec
Electrode Type:	pH
Blank Option:	No Blank
Calculations:	Sample Calc. by Volume
Dilution Option:	Disabled
Titrant Name:	0.1N NaOH
Titrant Conc.:	0.1000 N (eq/L)
Analyte Size:	100.00 mL
Analyte Entry:	Fixed
Maximum Titrant Volume:	20.000 mL
Stirring Speed:	1400 rpm
Potential Range:	-2000.0 to 2000.0 mV
Volume/Flow Rate:	25 mL/50 mL/min
Signal Averaging:	1 Reading
Significant Figures:	XXXXX

Concentration of Phosphoric Acid0.00000–1.0000⁻² mol/L**Calculations:**

Calculations: Sample Calc. by Volume
 Titrant units: N (eq/L)
 Titrant volume dosed: V (L)
 Final result units: M (mol/L)
 Titrant conc.: 0.1000 eq/L
 (sample/titrant): 1.000 mol/eq
 Sample volume: 100.000 mL

$$\frac{\text{mol}}{\text{L}} = \frac{V(\text{L}) * 100 * 0.1 * 1.0}{100.00}$$

Results:

Titration Report

Method Name: Concentration of H3PO4
 Time & Date: 10:11 April 23, 2010
 Titration ID: Ti_00014

Titration Results

Method Name: Concentration of H3PO4
 Time & Date: 10:11 April 23, 2010
 Analyte size: 100.000 mL

Equivalence point 1:

pH: 4.677
 Volume: 4.397 mL
 Results: 4.3972e⁻³ mol/L

Equivalence point 2:

pH: 8.916
 Volume: 4.429 mL
 Results: 4.4293e⁻³ mol/L

Titration Duration: 6:01 [mm:ss]



Titration went to Completion

Operator name: _____

TITRATION THEORY

HI 901 and HI 902

AUTOMATIC POTENTIOMETRIC TITRATOR



Contents

1	GENERAL REVIEW OF TITRATION THEORY	5
1.1	Introduction to Titrations	5
1.2	Uses of Titrations	5
1.3	Advantages and Disadvantages of Titrations	6
2	TYPES OF TITRATION	7
2.1	Titrations According to The Measurement Method	7
2.1.1	Amperometric Titrations	7
2.1.2	Potentiometric Titrations	7
2.1.3	Spectrophotometric Titrations	8
2.2	Titrations According to The Reaction Type	9
2.2.1	Acid - Base Titrations	9
2.2.2	Argentometric Titrations	10
2.2.3	Complexometric Titrations	10
2.2.4	Ion Selective Titrations	11
2.2.5	Non-aqueous Solvent Acid - Base Titrations	11
2.2.6	Precipitation Titrations	12
2.2.7	Redox Titrations	12
2.2.8	Karl Fischer Titration	13
2.3	Titrations According to The Titration Sequence	13
2.3.1	Back Titrations	13
2.3.2	Multiple Endpoint Titrations	14
3	INTRODUCTION TO TITRATION APPARATUS AND TYPICAL TITRATION PROCEDURE	15
3.1	Manual Titration	15
3.2	Automatic Titration	16
4	TITRATION RESULTS	17
4.1	Accuracy	17
4.2	Repeatability	17

TITRATION THEORY

4.3	Sources of Error	17
4.3.1	Sampling Errors	17
4.3.2	Errors with Titrant and Standard	17
4.3.2.1	Preparation Errors	17
4.3.2.2	Dispensing Errors	18
4.3.3	Chemical Reaction Errors	18
4.3.4	Endpoint Determination Errors	18
5	CALCULATIONS	19
5.1	Sample Calculation	19
5.2	Standardize Titrant	19
5.3	Blank Titration	20
5.4	Multiple Endpoint Titration	20
5.5	Back Titration	21
6	GLOSSARY	22

1 GENERAL REVIEW OF TITRATION THEORY

1.1 Introduction to Titrations

A titration is a quantitative, volumetric procedure used in analytical chemistry to determine the concentration of an analyte (the species being measured) in solution. The concentration of the analyte is determined by slowly adding a titrant (reagent) to the solution. As the titrant is added, a chemical reaction occurs between the titrant and the analyte.

Titration reactions are relatively fast, simple reactions that can be expressed using a chemical equation. The titration reaction continues as the titrant is added until all of the analyte is consumed and the analyte reacts completely and quantitatively with the titrant.

The point at which all of the analyte has been reacted is called the equivalence point, also known as the theoretical or stoichiometric endpoint. This point is accompanied by an abrupt physical change in the solution, which sharply defines the endpoint of the reaction. The physical change associated with the titration endpoint can be produced by the titrant or an indicator and can be detected either visually or by some other physical measurement.

Titration cannot be used to determine the quantity of all analytes. The chemical reaction between the titrant and analyte must fulfill four requirements:

- The reaction must be fast and occur within approximately one second after the titrant is added
- The reaction must go to completion
- The reaction must have well-known stoichiometry (reaction ratios)
- A convenient endpoint or inflection point

Titration is highly precise and can provide many advantages over alternative methods. Titrations are quickly performed and require relatively simple apparatus and instrumentation.

1.2 Uses of Titrations

Titration can be used in many applications, including:

- Acid content of plant effluents, food (e.g.: cheese and wine), plating and etching baths, petroleum products, drugs
- Base content of fertilizer (containing ammonia), bleach, minerals
- Hardness in water
- Metal content of alloys, minerals, ores, clays, waters, plating baths, paints, paper, plant materials, biological fluids, petroleum products
- Moisture content in foodstuffs, petrochemicals, pharmaceutical products, and plastics
- Redox reagent concentrations such as available chlorine in potable water, peroxide, traces of oxidants and reductants in food, reductants in high temperature or high pressure boiler water, vitamin analysis

TITRATION THEORY

1.3 Advantages and Disadvantages of Titrations

Some advantages of titrations as an analytical technique are:

- More precise results than many instrumental methods, such as measurement by electrode, the accuracy of the measurement is up to 0.1%
- Simple methods, reasonable capital costs, and easy training
- Suitability to measure major components of a mixture or product
- Automation can reduce time and labor spent on each analysis

Some disadvantages of titrations are:

- Time it takes to prepare standards and titrants
- Good technique is required to achieve precise results (training and practice required)
- Not suitable for determining trace or minor components of a mixture or product
- Limited dynamic range, it may require additional sample preparations (dilution) and repeat analyses

2 TYPES OF TITRATIONS

2.1 Titrations According to The Measurement Method


2.1.1 Amperometric Titrations

An amperometric titration is performed by placing two electrodes (often a metal ISE and a reference electrode) into the sample solution and holding the potential of the metal electrode at a selected voltage. The current that flows, due to the oxidation or reduction of a reactant or product, is plotted vs. volume of titrant to provide the titration curve and locate the equivalence point. Changes in the current are due to changes in the concentration of a particular species (being oxidized or reduced at the electrode).

Generally the reaction between the analyte and titrant forms a new species. Depending on the titration, the reactants are electroactive and the products are not, or vice-versa. Amperometric titration curves look like two straight lines intersecting at the equivalence point, this is due to the change in the electroactivity of the solution.

Many metal ions can be amperometrically titrated using a precipitation, complexation or redox reaction. Some metal ions and species that can be determined in this manner include silver, barium, halides, potassium, magnesium, palladium, molybdate, sulfate, tungstate, zinc, bismuth, cadmium, fluoride, indium, thallium, iodine, and gold.

Figure 1 shows four amperometric titrations and their endpoints. In graph "A" the analyte is electroactive and gives current but the reacted species does not. In "B" the reactant is not active but the titrant is. In "C" both the analyte and titrant are active and both give current flow. Graph "D" shows the same situation as "B"; however, the current has an opposite sign (the titrant is reduced).



2.1.2 Potentiometric Titrations

Potentiometric titrations are done by measuring the voltage across the solution using an electrode system. An electrode system consists of an indicator electrode and a reference electrode. As titrant is added the variations in the potential of the indicator electrode, with respect to the reference electrode, are monitored to show the progress of the titration.


Potentiometry is the measurement of a potential under conditions of zero current flow. The

TITRATION THEORY

measured potential can then be used to determine the analytical quantity of interest, generally a component concentration of the analyte solution. The potential that develops in the electrochemical cell is the result of the free energy change that would occur if the chemical phenomena were to proceed until the equilibrium condition has been satisfied.

There are many types of titrations where potentiometry can be used, e.g., pH electrodes for acid-base titrations, platinum ORP electrodes in redox titrations, ion selective electrodes, such as chloride or fluoride for a specific ion titration, and silver electrodes for argentometric (silver-based) titrations. An example of potentiometric titrations are shown below. Figure 2 "A" is the pH of a solution vs. the volume of titrant and "B" is the potential from a chloride electrode vs. the volume of AgNO_3 .

Figure 2




2.1.3 Spectrophotometric Titrations

The name comes from the method used to detect the endpoint of the titration, not its chemistry. Highly colored indicators that change color during the course of the titration are available for many titrations. More accurate data on the titration curve can be obtained if the light absorption is monitored instrumentally using a light source, a simple monochromator and a photodetector, rather than visually determining the color or light absorption change. Light absorption by either an indicator or by one of the reactants or products can be used to monitor the titration.

In the first titration curve, Figure 3 "A", the absorption of a metal-indicator complex is being monitored. The absorption is constant while the metal is complexed by the EDTA titrant. The metal indicator complex was stripped, causing a sharp break in the titration curve. The point where all the metal is complexed and stripped from the indicator is the equivalence point. This point is marked by "e.p." on the graph.

In the second titration curve, Figure 3 "B", the metal complex is being measured while being titrated with EDTA. The new complex being formed is not colored and does not absorb light. The extrapolated intersection of the two lines determines the equivalence point.

Figure 3



2.2 Titrations According to The Reaction Type

2.2.1 Acid-Base Titrations

Acid–base titrations are the most common type of titrations. They are based upon a reaction between an acid and a base, a stoichiometric neutralization, or the exchange of protons. Virtually all acid–base titrations are carried out using a strong acid or a strong base as the titrant. The endpoint of a titration carried out with a weak acid or a weak base would be difficult to detect due to a small change in pH at the equivalence point.

Chemical indicators can be used to determine the endpoint. The indicator will change color to signify that the end of the titration has been reached. The color of the indicator is dependent upon the concentration of ions in the solution. An acid–base indicator is composed of a conjugate weak acid–weak base pair, where the two forms exhibit different colors depending on the pH of the solution. For an indicator, the acid ionization constant K_a is usually written as:

$$K_a = \frac{[H_3O^+][In^-]}{[HIn]}$$

HIn is the acid form of the indicator and In^- is the base form. At the center of the change region, the ratio of $[In^-]$ to $[HIn]$ is one, $[H_3O^+] = K_a$ and $pH = pK_a$. The color change region is usually ± 1 pH unit around this point. Table 1 contains a list of some aqueous acid–base chemical indicators, as well as the pH range, the pK_a and the expected color (acid and base form). When choosing the proper indicator you should select one that has a pK_a as close to the endpoint of the titration.

When chemical indicators are not suitable, a potentiometric pH titration can also be used. The pH of the solution is plotted versus the volume of titrant added. Figure 4 shows a traditional strong acid–strong base titration curve. The graph shows the

Table 1

pH Range	Indicator	pK_a	Acid Form	Base Form
0.0 - 1.6	Methyl Violet		Yellow	Blue
1.2 - 2.8	Thymol Blue	1.65	Red	Yellow
3.2 - 4.4	Methyl Orange	3.46	Red	Yellow
3.8 - 5.4	Bromocresol Green	4.90	Yellow	Blue
4.8 - 6.0	Methyl Red	5.00	Red	Yellow
5.2 - 6.8	Chlorophenol Blue	6.25	Yellow	Red
6.0 - 7.6	Bromothymol Blue	7.30	Yellow	Blue
6.6 - 8.0	Phenol Red	8.00	Yellow	Red
7.4 - 9.0	Metacresol Purple	8.30	Yellow	Purple
8.0 - 9.6	Thymol Blue	9.20	Yellow	Blue
8.2 - 10.0	Phenolphthalein	9.50	Clear	Pink
9.4 - 10.6	Thymolphthalein		Clear	Blue
10.1 - 12.0	Alizarin Yellow R		Yellow	Red
11.4 - 12.6	Indigo Carmine		Blue	Yellow

TITRATION THEORY

volume of NaOH added to an acidic solution and the resulting pH of the solution. Note the abrupt change in the pH at the equivalence point.

2.2.2 Argentometric Titrations

Argentometric titrations use silver (nitrate) as the titrant and are generally precipitation titrations, as many silver salts are insoluble. These titrations are commonly used to titrate and determine the concentration of bromide, chloride, cyanide, iodide, and sulfide.

Argentometric titrations can be done with Mohr's indicator (when all of the chloride has reacted, a red silver chromate precipitate is formed) or the titration can be easily followed with a silver ISE (or chloride ISE for chloride titrations) and a reference electrode.

Figure 5 shows the titration of 50 mL of 0.1N NaCl with 0.1N AgNO₃. The potentiometric signal is from a chloride ISE and is plotted as pCl (- log [Cl⁻]).

Figure 4





Figure 5




2.2.3 Complexometric Titrations

A complex is a species where a central metal ion is covalently bonded to one or more electron donating groups called ligands. In a complexometric titration, metal ions are titrated using a titrant that binds strongly to it. Often these titrants contain EDTA or CDTA, polydentate ligands that form very stable coordination compounds with metal ions. The complexation reaction must be fast in order to be useful for direct titration. Some metal ions react too slowly with EDTA for a direct titration.

An indicator electrode that responds to the metal ion can be used to monitor the titration progress. The titration curve will appear similar to a usual potentiometric titration. Complexation indicators change color at the endpoint as all metal ions are "consumed", or complexed, by the titrant.

The titration curve will appear similar to a potentiometric titration when using an indicator electrode that responds to the metal ion (see Figure 6).

Figure 6



2.2.4 Ion Selective Titrations

The most popular ion selective titration is an acid-base titration. The hydrogen ion concentration is specifically measured and monitored during the titration process to locate the equivalence point. Using an ion selective electrode (ISE) as the indicator electrode, the potentiometric signal (in mV) is used to directly follow a specific ion's concentration (or activity).

Examples of ISE titrations include titrating fluoride with an aluminum titrant using a fluoride ISE, chloride with silver nitrate using a chloride ISE, sodium with a sodium ISE, etc. The equivalence point can be determined by plotting the mV value vs. the amount of titrant added.


2.2.5 Non-aqueous Solvent Acid-Base Titrations

Non-aqueous solvents must be used to titrate very weak acids and bases due to the inherent leveling effect water has on all acids and bases dissolved in it. A wide variety of weak acids and bases can be titrated using non-aqueous solvents. Mixtures of acids or bases can often be individually analyzed in a single sequential titration.

Titration of Acids

Weak acids with pK_a 's up to about 11 can be titrated in non-aqueous solvents. These include carboxylic acids, enols, phenols, imides, sulfonic acids, and inorganic acids. Water or lower alcohols are suitable for titrating medium to strong acids (pK_a less than 5). Titrating a weaker acid with a strong base titrant requires a solvent less acidic than water or ethanol/methanol. Solvents such as acetone, acetonitrile, t-butyl

Figure 7



TITRATION THEORY

alcohol, dimethylformamide, isopropanol and pyridine have been found to work well for acid-base titrations of strong, medium and weak acids/bases. Titrants include alcoholic potassium hydroxide and various sodium or potassium alkoxides in a 10:1 mixture of benzene/methanol. The best titrants are quaternary ammonium hydroxides (such as tetrabutylammonium hydroxide) due to good solubility of tetraalkylammonium salts of the titrated acids and the clean potentiometric titration curve obtained (see Figure 7).

Titration of Bases

Weak bases with pK_b 's up to about 11, which do not ionize with water, can be titrated in non-aqueous solvents. These bases include aliphatic and aromatic amines, basic nitrogen heterocycles, alkali metal and amine salts of acids, and many other organic basic compounds. Titrating a weak base with a strong acid titrant requires a basic solvent that is as weak as possible. Water and alcohols allow the titration of medium strength bases such as aliphatic amines ($pK_b = 4$ to 5), but not the titration of weaker bases such as pyridine ($pK_b = 8.8$). Glacial acetic acid works well for weak bases and has been used extensively. Less basic solvents such as acetone, acetonitrile, and nitromethane extend the range of titrable compounds.

The endpoint for non-aqueous titrations are usually determined potentiometrically using a pH glass electrode, a modified calomel or double junction reference electrode with a low-flow rate reference junction. Good potentiometric titration curves are obtained in most solvents, except those with very low dielectric constants such as benzene, chloroform and others, when high electrical resistance of the solvent causes unstable potentials.

2.2.6 Precipitation Titrations

Precipitation titrations allow for faster analysis compared to the old gravimetric analysis, where a precipitate is formed, filtered, dried and weighed to analyze a compound. Typically silver halides, silver thiocyanate and a few mercury, lead, and zinc salts are titrated using this method. The chemical reactions must form an insoluble salt and precipitate out quickly in order to be analyzed by this method. When the reaction is not quick, a back titration can be used. A measured excess of the precipitating reagent (titrant) is added to force the reaction to occur, and then unreacted titrant is then titrated with a standard solution of another reagent.

2.2.7 Redox Titrations


There are a number of oxidation-reduction reactions that can be used to determine unknown concentration by titration. If the reaction goes to completion, is fast and has an analytical signal available to follow it, a titration can be performed. The term "fast" means that each addition of titrant is reacted completely and the sensing electrode is able to detect the change in solution in less than one second.

Redox titrations are potentiometric titrations where the mV signal from a combination ORP (redox) electrode (usually with a platinum indicator electrode) is used to follow the reaction of oxidant/reductant. The electrode potential is determined by the Nernst equation and is controlled by the oxidant/reductant ratio.

Visual indicators such as Ferrion are also available. The oxidized and reduced form of the indicator will have different colors and can be used to determine the end point.

Various reductants can be determined by titrants with oxidants such as potassium permanganate, potassium chromate or iodine. Commonly used reductants that are used as titrants include sodium thiosulfate, and ferrous ammonium sulfate.

Figure 8



As with Acid-Base titrations the potential changes dramatically at the equivalence point.

2.2.8 Karl Fischer Titrations

This method is based on a well-defined chemical reaction between water and the Karl Fischer reagent. The chemistry provides excellent specificity for water determination. The method can be used to determine free and bound water in a sample matrix. The Karl Fischer method is widely considered to produce the most rapid, accurate and reproducible results and has the largest detectable concentration range spanning 1 ppm to 100%.

The determination of water content is one of the most commonly practiced methods in laboratories around the world. Knowledge of water content is critical to understanding chemical and physical properties of materials and ascertaining product quality. Water content determination is conducted on many sample types including pharmaceuticals and cosmetics, foods and natural products, organic and inorganic compounds, chemicals, solvents and gases, petroleum and plastic products as well as paints and adhesives. The KF method is verifiable and can be fully documented. As a result, Karl Fischer titration is the standard method for analysis of water in a multitude of samples as specified by numerous organizations including the Association of Official Analytical Chemists, the United States and European Pharmacopoeia, ASTM, American Petroleum Institute, British Standards and DIN.

2.3 Titrations According to The Titration Sequence

2.3.1 Back Titrations

Back titrations are generally used when a reaction is too slow to be directly accomplished using a "direct" titration, where the reaction goes to completion within a few seconds. In a back titration, a large excess of a reagent is added to the sample solution, helping a slow reaction to go to completion. The unreacted, excess reagent is then titrated. The difference in the total volume of the first reagent added and amount determined from the second titration is the quantity of reagent required to complete the first reaction.


TITRATION THEORY

2.3.2 Multiple Endpoint Titrations

Under certain conditions, some titrations can exhibit more than one equivalence point and be titratable to the individual endpoints to determine the concentration of each individual component. Examples of these types of titrations include acid-base (where different strength acid or bases are in a mixture), redox (where each species has a different reduction potential), complexometric (where different species are separately titratable), and acid-base using polyprotic acids (the pK_a of the different protons varies enough to separate them).

Figure 9 shows three different types of multiple endpoint titrations. "A" shows the titration of a polyprotic acid. The different acid strengths of the first and second proton can be determined. "B" illustrates a mixture of two different metal redox species, where the different redox potentials allow the species to be separated. "C" is the titration of a solution containing strong, weak, and very weak acids.

Figure 9



3 INTRODUCTION TO TITRATION APPARATUS AND TYPICAL TITRATION PROCEDURE



3.1 Manual Titration

Apparatus required for manual titration include:

- Volumetric Burette, for precisely controlled delivery of titrant to the reaction vessel
- An Erlenmeyer, or similar flask, that facilitates constant mixing or swirling required to ensure solution homogeneity
- Volumetric pipettes for the precise addition of samples and indicator solutions
- Titrant solutions of known concentration
- A visual or instrumental indicator for detecting the completion of the reaction

A typical manual titration consists of the following steps:

1. A volumetric pipette is typically used to add a known volume of sample to the flask
2. An indicator solution or instrument probe is added to the flask
3. A burette is used to measure the addition of titrant to the flask and dispense titrant in a controlled manner
4. Titrant is added via the burette until the method indication signals the reaction endpoint
5. The concentration of analyte is calculated based on the concentration and volume of titrant required to reach the endpoint



TITRATION THEORY

3.2 Automatic Titration

Automatic titrators are high-precision analytical instruments that deliver the titrant, monitor the physical change associated with the titration reaction, automatically stop at the endpoint and calculates the concentration of the analyte. Automatic titrators are best for repetitive titrations and high-accuracy analyses.

An automatic titrator must have an accurate liquid dispensing system. In high accuracy systems like the HI 900-series titrators, the liquid dispensing system consists of a stepper-motor driven piston syringe burette capable of accurately and precisely dispensing very small volumes of titrant, a valve system to switch between titrant intake and outlet and a dispensing tip. These three main subsystem components must be as accurate as possible, with very low gear backlash in the burette pump, minimal piston seal flexing, precision ground inner diameter of the glass syringe, a low dead volume valve, minimal evaporation/permeation, and chemically resistant tubing.

Apparatus required for automatic titration include:

- An automatic titrator, equipped with a burette
- A beaker
- An electronic stirring system, either a propeller stirrer or a magnetic stir bar and stir plate
- Volumetric pipettes for the precise addition of samples
- Standard titrant solutions of known concentration
- An electrode system that can be used to determine the endpoint of the titration

A typical automatic titration consists of the following steps:

1. Set up the automatic titrator according to the manufacturer's instructions
2. A volumetric pipette is typically used to add a known volume of sample to the beaker
3. Submerge the propeller stirrer or add the stir bar to the beaker, and turn on
4. Start the titration, the titrator will automatically stop at the endpoint and determine the concentration of the analyte

4 TITRATION RESULTS

4.1 Accuracy

The factors most critical to achieving accurate results with the HI 900 titration systems are the concentration of the sample, size of the sample and having an optimized set of method parameters.

4.2 Repeatability

Repeatability, or the agreement between replicate determinations, is expressed quantitatively as the relative standard deviation (RSD).

4.3 Sources of Error

One of the advantages of volumetric analysis is excellent accuracy and precision. The sources of error can be grouped into sampling, titrant and standards, chemical reactions, endpoint determination and calculations.

4.3.1 Sampling Errors

- Selection of a non-homogeneous or non-representative sample
- Sample changed or was contaminated during collection, storage or transfers
- Poor technique when transferring sample to beaker or flask
- Errors in the balance, calibrate and check balance regularly

4.3.2 Errors with Titrant and Standard

4.3.2.1 Preparation Errors

Incorrect preparation due to:

- Poor technique in weighing the salt or when transferring to volumetric glassware
- Low-purity of salts or water used to make titrant and standard
- Dirty or wet glassware
- Improper storage of titrant or standard which allows water gain, evaporation or deterioration
- Failure to standardize frequently to adjust for change in titrant
- Failure to flush titrator tubing with a volume of titrant before standardizing
- Volume errors from pipettes and volumetric flasks, grade A glassware is required
- Balance errors when weighing out salts, calibrate and check balance regularly

TITRATION THEORY

4.3.2.2 Dispensing Errors

Incorrect dispensing due to:

- Dead valve volume and leaking valve
- Inaccuracy in motor drive and gear lash/ backlash
- Poor burette/ piston seal
- Non-uniform diameter of burette glass cylinder
- Chemical incompatibility with tubing or bubble generation
- Density/ temperature changes in titrant

4.3.3 Chemical Reaction Errors

- Inappropriate solvent or sample resulting in side reactions
- Poor mixing of the titrant and solvent or sample in the titration vessel
- Reaction between titrant and sample is not rapid
- Reaction does not go to completion
- Reaction has side reactions

4.3.4 Endpoint Determination Errors

Most manual titrations use a visual indicator to indicate when the endpoint is reached and the titration should be stopped. Automatic titrators use instrumental methods to determine the end of a titration and the equivalence point. There are two predominant methods used to determine the equivalence point, first derivative and second derivative.

The inflection point of the titration curve (mV vs. Volume) is normally assumed to be the equivalence point. The first derivative is often used to determine the inflection point. The maximum value of the first derivative (dmV vs. dV) corresponds to the theoretical equivalence point. During a titration it is rare to have a data point exactly at the first derivative maximum, the maximum value is determined by interpolating the first derivative data points.

The second derivative ($d^2 \text{ mV vs. } dV^2$) can also be used to determine the equivalence point, and can offer advantages over the first derivative method. Second derivatives have increased sensitivity to smaller inflection points and easier numerical evaluation of the actual equivalence point. The value where the second derivative is equal to zero is the equivalence point. The second derivative requires fewer points located near the equivalence point, where data is often not obtained or not as reliable.

Errors in determining the endpoint can result from:

- Incorrect signals from the sensor
- Sensor drift
- Sensor or instrument has slow response, keep sensors in good condition
- Inappropriate setting on the titrator

5 CALCULATIONS

The main variables used in calculating a result from a titration are the sample volume, the concentration of the titrant, and the volume of titrant required to reach the equivalence point. At the equivalence point, an equal number of equivalents of the analyte and titrant has been added.

5.1 Sample Calculation

By Mass

$$C_{sample} = \frac{V_{titrant} \times C_{titrant} \times Ratio \times FW_{analyte}}{m_{sample}} \times 100$$

C sample	Sample Concentration (g/100g)
V titrant	Volume of titrant (L)
C titrant	Titrant Concentration (eq/L)
Ratio	Equivalence ratio of analyte/ titrant (mol analyte/ eq titrant)
FW analyte	Formula Weight of the Analyte (g/mol)
m sample	Mass of sample (g)

By Volume

$$C_{sample} = \frac{V_{titrant} \times C_{titrant} \times Ratio \times FW_{analyte}}{V_{sample}} \times 100$$

C sample	Sample Concentration (g/100mL)
V titrant	Volume of titrant (L)
C titrant	Titrant Concentration (eq/L)
Ratio	Equivalence ratio of analyte/ titrant (mol analyte/ eq titrant)
FW analyte	Formula Weight of the Analyte (g/mol)
V sample	Volume of Sample (mL)

5.2 Standardize Titrant

Titration standardization is the second most important calculation in titrations. A primary standard is titrated in order to determine the concentration of the titrant. This is essentially a typical titration calculated in "reverse", where the concentration of the solution is known and the titrant is the unknown.

By Mass

$$C_{titrant} = \frac{m_{standard} \times Ratio}{FW_{standard} \times V_{titrant}}$$

C titrant	Titrant Concentration (N)
m standard	Mass of Standard (g)
Ratio	Equivalence ratio of titrant/standard (eq titrant/ mol standard)
FW standard	Formula Weight of the Standard (g/mol)
V titrant	Volume of Titrant (L)

TITRATION THEORY

By Volume

$$C_{\text{titrant}} = \frac{V_{\text{standard}} \times (1 \text{ L} / 1000 \text{ mL}) \times C_{\text{standard}}}{V_{\text{titrant}}}$$

C titrant	Concentration of titrant (N)
V standard	Volume of Standard (mL)
C standard	Concentration of standard (eq/L)
V titrant	Volume of Titrant (L)

5.3 Blank Titration

In a blank titration a pre-titration is performed, often times on the solvent to be used for the sample titration, and the titrant volume required to reach the endpoint is noted. This blank value nullifies error due to titrant required to react with the components of the titration solution matrix. The basic titration equation can be used for a blank titration, with the single modification that the volume of titrant used in the blank titration should be subtracted from the regular titration titrant volume.

$$C_{\text{sample}} = \frac{C_{\text{titrant}} \times (V_{\text{sample}} - V_{\text{blank}}) \times \text{Ratio} \times \text{FW}_{\text{analyte}}}{m_{\text{sample}}} \times 100$$

C Sample	Sample Concentration (g/100g)
C titrant	Titrant Concentration (eq/L)
V sample	Volume of Titrant required for the sample (L)
V blank	Volume of Titrant required for the blank (L)
Ratio	Equivalence ratio of analyte/ titrant (mol analyte/ eq titrant)
FW analyte	Formula Weight of the Analyte (g/mol)
m sample	Mass of sample (g)

5.4 Multiple Endpoint Titration

Some titrations have two or more endpoints, each corresponding to the equivalence point for a specific reaction. Multiple endpoint titrations are similar to a blank titration in that the volume of titrant required to reach the first endpoint is subtracted from the titrant volume used to reach the next sequential endpoint.

$$C_{\text{sample 1}} = \frac{V_{\text{titrant 1}} \times C_{\text{titrant}} \times \text{Ratio} \times \text{FW}_{\text{analyte 1}}}{m_{\text{sample}}} \times 100$$

$$C_{\text{sample 2}} = \frac{(V_{\text{titrant 2}} - V_{\text{titrant 1}}) \times C_{\text{titrant}} \times \text{Ratio} \times \text{FW}_{\text{analyte 2}}}{m_{\text{sample}}} \times 100$$

$$C_{\text{sample 3}} = \frac{(V_{\text{titrant 3}} - V_{\text{titrant 2}}) \times C_{\text{titrant}} \times \text{Ratio} \times \text{FW}_{\text{analyte 3}}}{m_{\text{sample}}} \times 100$$

C sample1	Sample 1 Concentration (g/100g)
C sample2	Sample 2 Concentration (g/100g)
C sample3	Sample 3 Concentration (g/100g)
V titrant 1	Volume of titrant required to reach the first end point (L)
V titrant 2	Volume of titrant required to reach the second end point (L)
V titrant 3	Volume of titrant required to reach the third end point (L)
C titrant	Concentration of titrant (N)
Ratio	Equivalence ratio of analyte/ titrant (mol analyte/ eq titrant)
FW analyte 1	Formula Weight of the Analyte 1 (g/mol)
FW analyte 2	Formula Weight of the Analyte 2 (g/mol)
FW analyte 3	Formula Weight of the Analyte 3 (g/mol)
m sample	Weight of Sample (mL)

5.5 Back Titration

The equation used in back titration calculations is also similar to the equation for a blank titration. Instead of subtracting the initial amount of titrant needed to react with the blank, the amount of second titrant needed to react with the excess titrant added in the first titration is subtracted from the amount of the first titrant added. The difference between the two amounts is the amount of titrant necessary to reach the first equivalence point.

$$C_{sample} = \frac{(C_{titrant\ 1} \times V_{titrant\ 1} - C_{titrant\ 2} \times V_{titrant\ 2}) \times Ratio \times FW_{analyte}}{V_{sample}} \times 100$$

C sample	Sample Concentration (g/100mL)
C titrant 1	Concentration of titrant 1 (N)
V titrant 1	Volume of titrant 1 (L)
C titrant 2	Concentration of titrant 2 (N)
V titrant 2	Volume of titrant 2 (L)
Ratio	Equivalence ratio of analyte/ titrant (mol analyte/ eq titrant)
FW analyte	Formula Weight of the analyte (g/mol)
V sample	Volume of sample (mL)

TITRATION THEORY

6 GLOSSARY

Acid

A chemical species that can donate one or more protons (hydrogen ions).

Acid-Base Titration

Stoichiometric neutralization titrations, based upon the reaction that occurs between an acid and base.

Activity

A physical property corresponding to the concentration of all ions in a solution. Electrodes respond to activity.

Amperometric Titration

Titrations where the current flow between two electrodes (often a metal electrode and a reference electrode) are used to monitor the titration progress.

Analyte

The chemical species being measured in a titration.

Argentometric Titration

Titrations that use silver (nitrate) as the titrant. These titrations are typically precipitation titrations.

Automatic Titrator

An instrument designed to automatically carry out a titration. It will add the appropriate amount of titrant, determine the endpoint and calculate the results.

Back Titration

A type of titration where an excess amount of titrant is added to a sample, forcing a sluggish reaction to go to completion. The excess reagent is then "back" titrated with a second titrant.

Base

A chemical species that can accept one or more protons (hydrogen ions).

Biamperometric Indication

Uses a double platinum pin electrode to measure the current flow through a titration solution.

Bivoltametric Indication

Uses a double platinum pin electrode to measure the voltage required to maintain a constant current flow through a titration solution while constant voltage is applied across the platinum elements of the electrode.

Burette

A graduated cylindrical piece of laboratory glassware that is used to dispense precise amounts of solution.

Complex Ion

A species where a central metal ion is covalently bonded to one or more electron donating groups called ligands.

Complexometric Titrations

Metal ions are titrated using a titrant that binds strongly to it. The titrants often contain Ethylenediaminetetraacetic Acid (EDTA) or Cyclohexylenedinitrilotetraacetic Acid (CDTA).

Endpoint

The point where a titration is stopped because a physical change in the solution has indicated a completed titration. Titration endpoints typically coincide with the equivalence point. A fixed value endpoint (pH or mV) can be used as well. The titration will stop at the desired point regardless if the titration is complete.

Equivalence point

The point where the quantity of titrant is stoichiometrically equal to the quantity of analyte.

Formal

The theoretical number of equivalents per liter of the solution. It is used in solutions where the exact concentration of a species may be affected by the other ions present, therefore the stated concentration may not be exactly correct.

Gravimetric Analysis

A quantitative determination of an analyte based on the mass of the solid.

Indicator Electrode

An electrode that responds to the species of interest. The electrode potential is proportional to the concentration or activity of that ion in the solution being measured.

Indicators

Chemical indicators are typically organic dyes that change form under different physical conditions, causing a color change that can be seen by an analyst. Typically used in manual titrations, chemical indicators have been replaced with electrometric indicators, which are used with automatic titrators.

Inflection Point

The point on a titration curve where the second derivative curve changes signs.

Ion Selective Electrode (ISE)

An electrode that responds to a specific ion. The electrode potential is proportional to the concentration or activity of that ion in the solution being measured.

Karl Fischer Titration

A titration that uses a chemical reaction that is specific for determining water.

Manual Titration

A titration that is carried out by hand. The analyst must add the appropriate amount of titrant, determine the endpoint and calculate the results.

Molar

The concentration of a solute in a solution.

Mole (mol)

A quantity of a chemical species. The molecular weight of a substance in grams is equal to the mass of one mole of the substance. One mole is equal to 6.022×10^{23} atoms or molecules.

Monochromator

A device that allows only a narrow range of wavelengths to pass through it by separating the light into different wavelengths.

Multiple Endpoint Titration

A titration that reacts multiple species in solution sequentially using the same titrant. The concentration of each analyte can be determined from their respective endpoints.

Nernst Equation

The fundamental equation relating cell voltage to the concentration of a solution.

TITRATION THEORY

Neutralization

A chemical reaction where an acid and a base react to form a neutral salt and water.

Non-aqueous

A solution that does not contain water.

Non-aqueous Titration

A titration that is performed in non-aqueous solutions, typically used to titrate very weak acids and bases to eliminate the leveling effect water has on all acids and bases dissolved in it.

Normal

The concentration of a solution which accounts for any stoichiometric difference between the various species in a solution.

Oxidation / Reduction Potential (ORP)

The measurement describing whether a species wants to donate or accept electrons from other species in a redox reaction. If a solution's reduction potential is higher than the species it is reacting with, it will typically gain electrons or be reduced. If the potential is lower than the species it is reacting with, it will typically lose electrons or be oxidized.

Oxidant

The species that is accepting electrons in a redox reaction.

Pipette

Scientific apparatus that is used to deliver precise volumes of liquids.

Polyprotic Acid

Acids that are capable of donating more than one proton per acid molecule.

Potentiometric Titration

A titration in which the endpoint is determined by monitoring the voltage of the solution using an electrode.

Precipitation Titration

A titration in which the analyte reacts with the titrant to form an insoluble compound. The endpoint is typically detected with an ISE sensitive to either the analyte or titrant.

Reagent

The chemical added in a titration that causes the given reaction to occur.

Reduction-Oxidation Reaction (redox)

A chemical reaction in which the atoms involved in the reaction have their oxidation numbers changed. Reduction is the gain of electrons, which decreases the oxidation number. Oxidation is the loss of electrons, which increases the oxidation number.

Reductants

The electron donor in a redox reaction.

Reference Electrode

An electrode that supplies a constant electrode potential. It is used in combination with an "indicator" electrode, allowing for the "indicator" electrode potential to be measured.

Relative Standard Deviation (RSD)

A measure of the amount of relative variation in a set of data. It is calculated by dividing the standard deviation by the mean: $RSD = (\text{Standard Deviation of } X) * 100 / (\text{Mean of } X)$

Repeatability

The variation in sample measurements taken by a single person or instrument under the same conditions.

Spectrophotometric Titration

A titration in which the endpoint is marked by a change in the color and/or color intensity.

Stoichiometry

The quantitative relationship of the reactants and products in a chemical reaction.

Titrant

The chemical added in a titration that causes the given reaction to occur.

Titration

A quantitative, volumetric procedure used in analytical chemistry to determine the concentration of an analyte in solution. The concentration of the analyte is determined by slowly adding a titrant to the solution. As the titrant is added, a chemical reaction between the titrant and the analyte occurs.

Titration Curve

A graph containing the physical data obtained for a titration. The data plotted is often an independent variable (volume of titrant) vs. a dependent variable (pH of the solution). From the titration curve, the equivalence point or endpoint can be determined.

TItR THEORY 901/902

04/10